



Event-specific Method for the Quantification of Maize Bt176 Using Real-time PCR

Protocol

13 July 2011

**Joint Research Centre
Institute for Health and Consumer Protection
Molecular Biology and Genomics Unit**

Method development:


Syngenta Seeds S.A.S.

Method validation:

European Union Reference Laboratory for GM Food and Feed (EURL-GMFF)
Molecular Biology and Genomics Unit


Drafted by

C. Savini (scientific officer)




Report Review

1) L. Bonfini



2) M. Querci



Scientific and technical approval

M. Mazzara (scientific officer)




Compliance to EURL Quality System

S. Cordeil (quality manager)



Authorisation to publish

G. Van den Eede (head of MBG Unit)



Address of contact laboratory:

European Commission, Joint Research Centre
Institute for Health and Consumer Protection (IHCP)
Molecular Biology and Genomics Unit
European Union Reference Laboratory for GM Food and Feed
Via E. Fermi 2749, 21027 Ispra (VA) - Italy

Content

1. GENERAL INFORMATION AND SUMMARY OF THE METHODOLOGY	4
2. VALIDATION AND PERFORMANCE CHARACTERISTICS	4
2.1 GENERAL	4
2.2 LIMIT OF DETECTION (LOD)	5
2.3 LIMIT OF QUANTIFICATION (LOQ)	5
2.4 MOLECULAR SPECIFICITY	5
3. PROCEDURE	6
3.1 GENERAL INSTRUCTIONS AND PRECAUTIONS	6
3.2 REAL-TIME PCR FOR QUANTITATIVE ANALYSIS OF Bt176 MAIZE	7
3.2.1 General	7
3.2.2 Calibration	7
3.2.3 Real-time PCR set-up	7
3.3 DATA ANALYSIS	9
3.4 CALCULATION OF RESULTS	9
4. MATERIALS	9
4.1 EQUIPMENT	9
4.2 REAGENTS AND SOLUTIONS	10
4.3 PRIMERS AND PROBES	11

1. General information and summary of the methodology

This protocol describes an event-specific real-time quantitative TaqMan[®] PCR procedure for the determination of the relative content of event Bt176 DNA to total maize DNA in a sample.

The PCR assay was optimised for use in real-time PCR instruments for plastic reaction vessels.

Template DNA extracted by means of suitable methods should be tested for quality and quantity prior to use in the PCR assay. Tests for the presence of PCR inhibitors (e.g. monitor run of diluted series, use of DNA spikes) are recommended.

For the specific detection of event Bt176 DNA, an 82-bp fragment of the integration region of the construct inserted into the plant genome (located at the 3' flanking DNA region) is amplified using specific primers. PCR products are measured at each cycle (real-time) by means of a target-specific oligonucleotide probe labelled with two fluorescent dyes: FAM (6-carboxyfluorescein) as a reporter dye at its 5' end and TAMRA (carboxytetramethylrhodamine) as a quencher dye at its 3' end.

For the relative quantification of event Bt176 DNA, a maize-specific reference system (*adh1*) amplifies a 135-bp fragment of the maize endogenous *alcohol dehydrogenase 1* gene (*adh1*), using specific primers and an *adh1* specific probe labelled with VIC and TAMRA: VIC as a reporter dye at its 5' end and TAMRA as a quencher dye at its 3' end.

The measured fluorescence signal passes a threshold value after a certain number of cycles. This threshold cycle is called the "Ct" value. For quantification of the amount of event Bt176 DNA in a test sample, the normalised Δ Ct values of the calibration samples are used to calculate by linear regression a reference curve Δ Ct-formula. The normalised Δ Ct values of the unknown samples are measured and, by means of the regression formula, the relative amount of Bt176 event DNA is estimated.

2. Validation and performance characteristics

2.1 General

The method has been optimised for suitable DNA extracted from maize leaves or grains containing mixtures of genetically modified and conventional maize.

The repeatability and trueness of the method were tested through an in-house validation study using DNA samples at different GM% contents.

Two replicated samples for each of the five GM-levels were analysed per run in parallel with both the GM-specific system (Bt176) and the target taxon-specific assay (*adh1*). Each replicate was tested in three adjacent wells and the difference of the average Ct values for the Bt176 and *adh1* system (Δ Ct) was plotted against the Δ Ct standard curve for the determination of the relative amount of event Bt176 of each replicate. In total, the quantification of the five GM levels was performed as an average of sixteen replicates per GM level.

2.2 Limit of detection (LOD)

According to the method developer, the relative LOD of the method is $\leq 0.01\%$ in 250 ng of total maize DNA. The relative LOD was not assessed in the in-house validation study.

2.3 Limit of quantification (LOQ)

According to the method developer, the relative LOQ of the method is $\leq 0.10\%$ in 250 ng of total maize DNA. The lowest relative GM content of the target sequence included in the in-house validation study was 0.1%.

2.4 Molecular specificity

According to the method developer, the method exploits a unique DNA sequence in the region of recombination at the 3' insert-to-plant DNA; the sequence is specific to Bt176 event and thus imparts event-specificity to the method.

The specificity of the Bt176 assay (forward/reverse primers and probe) was experimentally tested in real-time PCR by the method developer against DNA extracted from samples containing GM maize Bt11, MON810 and Bt176 maize.

According to the applicant, none of the above mentioned GM lines tested, except the positive control Bt176, produced amplification signals below Ct=37 in replicated samples when various amounts (ng) of total DNA per reaction were used.

The EURL-GMFF verified the specificity of the Bt176 assay by means of a multi-target analytical system for GMO detection in which primers and probes are already present in each well in lyophilized format at concentrations of 900 nM primers/250 nM probes¹.

¹ Querci M., Foti N., Bogni A., Kluga L., Broll H., Van den Eede G. Real-Time PCR-Based Ready-to-Use Multi-Target Analytical System for GMO Detection. Food Anal. Methods 2009; 2: 325–336

The Bt176 assay was tested against 100 ng DNA from 100% GM-lines containing the following events: maize 59122, MIR604, Bt10, Bt11, MON88017, LY038, MON863, TC1507, GA21, Bt176, MON810, NK603, MON89034, event 3272 and T25; soybean A2704-12, 356043, Roundup Ready 40-3-2 and MON89788; cotton 281-24-236, 3006-210-23, LLCotton25, MON531, MON1445, MON15985 and MON88913; rice LLRICE601, LLRICE62 and Bt63; oilseed rape T45, Ms8, Rf3, GT73, Rf1, Rf2, Ms1 and Topas 19/2; potato EH92-527-1; sugarbeet H7-1. Only the positive control sample Bt176 reacted with the Bt176 assay.

3. Procedure

3.1 General instructions and precautions

- The procedures require sterile conditions working experience.
- Laboratory organisation, e.g. "flow direction" during PCR-setup, should follow international guidelines, e.g. ISO 24276:2006.
- PCR-reagents should be stored and handled in a separate room where no nucleic acids (with exception of PCR primers or probes) or DNA degrading or modifying enzymes have been handled previously. All handling of PCR reagents and controls requires dedicated equipment, especially pipettes.
- Equipment used should be sterilised prior to use and any residue of DNA should be removed. All material used (e.g. vials, containers, pipette tips, etc.) should be suitable for PCR and molecular biology applications; it should be DNase-free, DNA-free, sterile and unable to adsorb protein or DNA.
- Filter pipette tips (protect against aerosol) should be used.
- Powder-free gloves should be used and changed frequently.
- Laboratory benches and equipment should be cleaned periodically with 10% sodium hypochlorite solution (bleach).
- Pipettes should be checked regularly for precision and calibrated, if necessary.
- All handling steps - unless specified otherwise - should be carried out at 0 - 4 °C.
- In order to avoid repeated freeze/thaw cycles aliquots should be prepared.

3.2 Real-time PCR for quantitative analysis of Bt176 maize

3.2.1 General

The PCR set-up for the taxon specific target sequence (*adh1*) and for the GMO (Bt176) target sequence is carried out in separate vials. Multiplex PCR (using differential fluorescent labels for the probes) has not been tested or validated.

The use of maximum 250 ng of template DNA per reaction well is recommended.

The method is developed for a total volume of 25 μ L per reaction mixture with the reagents as listed in Table 1 and Table 2.

3.2.2 Calibration

The calibration curve consists of five samples containing fixed percentages of Bt176 DNA in a total amount of 250 ng maize DNA. The GM content of the standard samples ranges from 4.5% to 0.1%.

A calibration curve is produced by plotting the Δ Ct values of calibration samples against the logarithm of the respective GM% contents; the slope (a) and the intercept (b) of the calibration curve ($y = ax + b$) are then used to calculate the mean GM% content of the blind samples based on their normalised Δ Ct values.

3.2.3 Real-time PCR set-up

1. Thaw, mix gently and centrifuge the required amount of components needed for the run. Keep thawed reagents on ice.
2. To prepare the amplification reaction mixtures add the following components (Tables 1 and 2) in two reaction tubes (one for the Bt176 assay and one for the *adh1* assay) on ice in the order mentioned below (except DNA).

Table 1. Amplification reaction mixture in the final volume/concentration per reaction well for *adh1*.

Component	Final concentration	μ L/reaction
Supplemented 2x Sigma JumpStart ReadyMix	1x	12.5
50x Endogenous Assay Stock	1x	0.5
Nuclease free water	#	7.0
Template DNA (max 250 ng)	#	5
Total reaction volume:		25

Table 2. Amplification reaction mixture in the final volume/concentration per reaction well for Bt176.

Component	Final concentration	$\mu\text{L}/\text{reaction}$
Supplemented 2x Sigma JumpStart ReadyMix	1x	12.5
50x Bt176Assay Stock	1x	0.5
Nuclease free water	#	7.0
Template DNA (max 250 ng)	#	5
Total reaction volume:		25

- Mix gently and centrifuge briefly.
- Prepare two reaction tubes (one for the Bt176 and one for the *adh1* master mixes) for each DNA sample to be tested (standard curve samples, unknown samples and control samples).
- Add to each reaction tube the correct amount of master mix (e.g. $20 \times 3 = 60 \mu\text{L}$ master mix for three PCR repetitions). Add to each tube the correct amount of DNA (e.g. $5 \times 3 = 15 \mu\text{L}$ DNA for three PCR repetitions). Vortex each tubes for approx. 10 sec. This step is mandatory to reduce to a minimum the variability among the repetitions of each sample.
- Spin down the tubes in a microcentrifuge. Aliquot $25 \mu\text{L}$ in each well. Seal the reaction plate with optical cover or optical caps. Centrifuge the plate at low speed (e.g. approximately $250 \times g$ for 1 minute at $4 \text{ }^\circ\text{C}$ to room temperature) to spin down the reaction mixture.
- Place the plate into the instrument.
- Run the PCR with cycling conditions described in Table 3

Table 3. Cycling program for maize Bt176/*adh1*

Step	Stage	T $^\circ\text{C}$	Time (sec)	Acquisition	Cycles
1	uracil-N-glycosylase (UNG)	$50 \text{ }^\circ\text{C}$	120	No	1
2	Initial denaturation	$95 \text{ }^\circ\text{C}$	600	No	1
3	Denaturation	$95 \text{ }^\circ\text{C}$	15	No	40
	Amplification Annealing & Extension	$60 \text{ }^\circ\text{C}$	60	Yes	

3.3 Data analysis

After the real-time PCR, analyse the run following the procedure below:

- a) Set the threshold: display the amplification curves of one assay (e.g. Bt176) in logarithmic mode. Locate the threshold line in the area where the amplification profiles are parallel (exponential phase of PCR) and where there is no "fork effect" between repetitions of the same sample. Press the "update (or apply)" button to ensure changes affect Ct values. Switch to the linear view mode by clicking on the Y axis of the amplification plot, and check that the threshold previously set falls within the geometric phase of the curves.
- b) Set the baseline: determine the cycle number at which the threshold line crosses the first amplification curve and set the baseline three cycles before that value (e.g. earliest Ct = 25, set the baseline crossing at $Ct = 25 - 3 = 22$).
- c) Save the settings.
- d) Repeat the procedure described in a) and b) on the amplification plots of the other assay (e.g. *adh1*).
- e) Save the settings and export all the data into an Excel file for further calculations.

3.4 Calculation of results

After having defined a threshold value within the logarithmic phase of amplification as described above, the instrument's software calculates the Ct values for each reaction.

The Reference ΔCt curve is generated by plotting the ΔCt values measured for the calibration points against the logarithm of the GM% content, and by fitting a linear regression line into these data.

Thereafter, the regression formula is used to estimate the relative amount (%) of Bt176 event in the unknown samples of DNA.

4. Materials

4.1 Equipment

- Real-time PCR instrument for plastic reaction vessels (glass capillaries are not recommended for the described buffer composition)

- Plastic reaction vessels suitable for real-time PCR instrument (enabling undisturbed fluorescence detection)
- Software for run analysis (mostly integrated in the software of the real-time PCR instrument)
- Microcentrifuge
- Micropipettes
- Vortex
- Rack for reaction tubes
- 1.5/2.0 mL reaction tubes

4.2 Reagents and solutions

- Sigma JumpStart Taq ReadyMix (2x), Sigma Aldrich Ltd Cat No P-2893
- Sulforhodamine 101, Sigma Cat No S-7635
- 1 M MgCl₂, Sigma Aldrich Ltd Cat No M-1028

10000x Sulforhodamine 101 stock:

Re-suspend 227.5 mg of Sulforhodamine 101 in 250 mL nuclease free water to make a 1.5 mM stock solution.

Vortex well and store at -20 °C.

Supplemented 2x Sigma JumpStart ReadyMix:

For 50 mL: to Sigma Jumpstart Taq ReadyMix (2X), add:

- 550 µL of 1 M MgCl₂
- 20 µL 10000x Sulforhodamine 101.

Vortex well and store at 4 °C for up to 1 year.

50x *adh1* Stock:

(1x concentration = 300 nM Forward primer, 300 nM Reverse primer, 200 nM Probe)

For 1 mL of 50x concentration, in amber (or other light-protective) reaction tube, mix:

- 15 µL of Forward Primer (1000 pmol/µL)
- 15 µL of Reverse Primer (1000 pmol/µL)
- 100 µL of Probe (100 pmol/µL)
- 870 µL nuclease-free water.

Vortex well and store at 4 °C for up to 1 year.

50x Bt176 Assay Stock:

(1x concentration = 100 nM Forward primer, 600 nM Reverse primer, 200 nM Probe)

For 1 mL of 50x concentration, in amber (or other light-protective) reaction tube, mix:

- 5 µL of Forward Primer (1000 pmol/µL)
- 30 µL of Reverse Primer (1000 pmol/µL)
- 100 µL of Probe (100 pmol/µL)
- 865 µL nuclease-free water.

Vortex well and store at 4 °C for up to 1 year.

4.3 Primers and Probes

Name	Oligonucleotide DNA Sequence (5' to 3')
<i>Bt176</i> target sequence	
Bt176 F primer	5' – GGC CGT GAA CGA GCT GTT -3'
Bt176 R primer	5' – GGG AAG AAG CCT ACA TGT TTT CTA A -3'
Bt176 Probe	FAM 5'- AGC AAC CAG ATC GGC CGA CAC C -3' TAMRA
<i>adh1</i> target sequence	
Zm <i>adh1</i> F primer	5' – CGT CGT TTC CCA TCT CTT CCT CC-3'
Zm <i>adh1</i> R primer	5' – CCA CTC CGA GAC CCT CAG TC -3'
Zm <i>adh1</i> Probe	VIC 5' – AAT CAG GGC TCA TTT TCT CGC TCC TCA-3' TAMRA