



Event-specific Method for the Quantification of oilseed rape MS11 by Real-time PCR

Validated Method

Method development:

Bayer CropScience AG

1. General information and summary of the methodology

This protocol describes an event-specific real-time quantitative TaqMan[®] PCR (polymerase chain reaction) procedure for the determination of the relative content of oilseed rape event MS11 DNA to total oilseed rape DNA in a sample.

Template DNA extracted by means of suitable methods should be tested for quality and quantity prior to use in PCR assays. Tests for the presence of PCR inhibitors (e.g. monitor run of diluted series, use of DNA spikes) are also recommended to ensure suitability of the extracted DNA.

For the detection of GM event MS11, a 124 bp fragment of the region spanning the 3' insert-to-plant junction in oilseed rape MS11 is amplified using specific primers. PCR products are measured during each cycle (real-time) by means of a target-specific oligonucleotide probe labelled with FAM (6-carboxyfluorescein) as reporter dye at its 5' end and MBG (minor groove binder) as non-fluorescent quencher dye at its 3' end.

For the relative quantification of GM event MS11, an oilseed rape taxon-specific system amplifies a 101 bp fragment of an oilseed rape *cruciferin A* (*CruA*) endogenous gene (Accession number, GeneBank: X14555), using *CruA* gene-specific primers and a *CruA* gene-specific probe labelled with JOE as reporter dye at its 5' end and BHQ-1 (Black Hole Quencher[®] 1) as non-fluorescent quencher dye at its 3' end

The measured fluorescence signal passes a threshold value after a certain number of cycles. This threshold cycle is called the "Cq" value. For quantification of the amount of MS11 DNA in a test sample, Cq values for the MS11 and the *CruA* systems are determined for the sample. Standard curves are then used to estimate the relative amount of MS11 DNA to total oilseed rape DNA.

2. Validation and performance characteristics

2.1 General

The method was optimised for suitable DNA extracted from genetically modified and conventional oilseed rape seeds and grain. Precision and trueness of the method were tested through an international collaborative ring trial using DNA samples at different GM contents.

2.2 Collaborative trial

The method was validated in an international collaborative study by the European Union Reference Laboratory for GM Food and Feed (EURL GMFF). The study was undertaken with twelve participating laboratories in November 2017.

A detailed validation report can be found at <http://gmo-crl.jrc.ec.europa.eu/StatusOfDossiers.aspx>.

2.3 Limit of detection (LOD)

According to the method developer, the relative LOD of the method is at least 0.023 % (related to mass fraction of GM material) in 200 ng of total suitable oilseed rape DNA. The relative LOD was not assessed in the collaborative study.

2.4 Limit of quantification (LOQ)

According to the method developer, the relative LOQ of the method is at least 0.08 % (related to mass fraction of GM material) in 200 ng of total suitable oilseed rape DNA. The lowest relative GM content of the target sequence included in the collaborative trial was 0.1 % (mass fraction of GM-material).

2.5 Molecular specificity

The method exploits a unique DNA sequence in the region spanning the 3' insert-to-plant junction in oilseed rape MS11 and is therefore event-specific for the event MS11. This was confirmed in the validation study.

3. Procedure

3.1 General instructions and precautions

- The procedures require experience of working under sterile conditions.
- Laboratory organisation, e.g. "forward flow direction" during PCR-setup, should follow international guidelines, e.g. ISO 24276:2006.
- PCR reagents should be stored and handled in a separate room where no nucleic acids (with exception of PCR primers or probes) or DNA degrading or modifying enzymes have been handled previously. All handling of PCR reagents and controls requires dedicated equipment, especially pipettes.
- All the equipment should be sterilised prior to use and any residue of DNA should have been removed. All material used (e.g. vials, containers, pipette tips, etc.) must be suitable for PCR and molecular biology applications. They must be DNase-free, DNA-free, sterile and unable to adsorb protein or DNA.
- Filter pipette tips protected against aerosol should be used.
- Powder-free gloves should be used and changed regularly
- Laboratory benches and equipment should be cleaned periodically, with 10% sodium hypochlorite solution (bleach).

- Pipettes should be checked regularly for precision and calibrated, if necessary.
- All handling steps, unless specified otherwise, should be carried out at room temperature.
- In order to avoid repeated freeze/thaw cycles aliquots should be prepared.

3.2 Real-time PCR for quantitative analysis of oilseed rape event MS11

3.2.1 General

The real-time PCR set-up for the taxon (*CruA*) and the GMO (event MS11) target sequences are carried out in separate vials. Multiplex qPCR (using differential fluorescent labels for the probes) has not been tested or validated by the EURL GMFF.

The method is developed for a total volume of 25 µL per reaction mixture for the GM (event MS11) and the taxon (*CruA*) assay with the reagents as listed in Table 2 and Table 3.

3.2.2 Calibration

The calibration curves have to be established on at least five samples. The first point of the calibration curve (S1) should be established for a sample containing 5 % oilseed rape MS11 DNA in a total of 300 ng of oilseed rape DNA (corresponding to 260870 oilseed rape haploid genome copies with one haploid genome assumed to correspond to 1.15 pg of oilseed rape genomic DNA)⁽¹⁾. Standards S2 to S5 are to be prepared by serial dilutions (dilution factor 4.0 for samples S2-S4 and dilution factor 3.0 for standard S5) according to Table 1 below.

Table 1. Copy number values of the standard curve samples

Sample code	S1	S2	S3	S4	S5
Total amount of oilseed rape DNA in reaction (ng)*	300	75	18.8	4.7	1.6
oilseed rape haploid genome copies	260870	65217	16304	4076	1359
MS11 copies	13043	3261	815	204	68

* Total nanograms are rounded to the integral value

A calibration curve is to be produced by plotting the C_q values against the logarithm of the target copy number for the calibration points. This can be done by means of spreadsheet software, e.g. Microsoft Excel, or directly by options available with the software.

The copy number measured for each unknown sample DNA is obtained by interpolation from the standard curves.

3.2.3 Real-time PCR set-up

1. Thaw, mix and centrifuge the components needed for the run. Keep thawed reagents on ice.
2. In two tubes on ice, add the components in the order mentioned below (except DNA) to prepare the reaction mixes for the MS11 oilseed rape specific system (Table 2) and the *CruA* reference gene system (Table 3). Please note that additional volume is included in the total to cover pipetting variability due to the viscosity of the solution.

Table 2. Amplification reaction mixture in the final volume/concentration per reaction well for the MS11 assay.

Component	Final concentration	µL/reaction
TaqMan [®] Universal PCR Master Mix (2x)	1x	12.5
SHA086 (10 µM)	400 nM	1.00
MDB371 (10 µM)	400 nM	1.00
TM280* (10 µM)	200 nM	0.50
Nuclease free water	-	5.0
DNA	-	5.0
Total reaction volume:		25 µL

*TaqMan[®] probe labelled with 6-FAM at its 5'-end and MGB at its 3'-end

Table 3. Amplification reaction mixture in the final volume/concentration per reaction well for the oilseed rape *CruA* assay.

Component	Final concentration	µL/reaction
TaqMan [®] Universal PCR Master Mix (2x)	1x	12.5
MDB510 (10 µM)	200 nM	0.50
MDB511 (10 µM)	200 nM	0.50
TM458* (10 µM)	200 nM	0.50
Nuclease free water	-	6.0
DNA	-	5.0
Total reaction volume:		25 µL

*TaqMan[®] probe is labelled with JOE at its 5'-end and BHQ1 at its 3'-end

3. Mix well and centrifuge briefly.
4. Prepare two 0.5 mL reaction tubes (one for the oilseed rape MS11 and one for the *CruA* system) for each DNA sample to be tested (standard curve samples, unknown samples and control samples).

5. Add into each reaction tube the amount of reaction mix for 3.5 PCR repetitions (70 μ L for the MS11 oilseed rape system and 70 μ L for the *CruA* system). Add to each tube the correct amount of DNA for 3.5 PCR repetitions (17.5 μ L DNA). The volume for the additional 0.5 repetition will ensure adequate volume when loading the samples. Vortex each tube for approx. 10 seconds. This step is mandatory to reduce to a minimum the variability among the repetitions of each sample.
6. Spin down the tubes. Aliquot 25 μ L for MS11 system and for the *CruA* reference system in each well.
7. Place an optical cover on the reaction plate and briefly centrifuge the plate.
8. Place the reaction plate in the real-time PCR apparatus (possibly apply a compression pad, depending on the model), according to the manufacturer's instructions and your Standard Operating Procedures and start the run.
9. Select FAM as reporter dye for the MS11 and JOE for the *CruA* reference system. Define MGB or non-fluorescent as quencher dye for MS11 specific system and BHQ or non-fluorescent for *CruA* reference system. Select ROX as the passive reference dye. Enter the correct reaction volume (25 μ L).
10. Run the PCR with the cycling program described in Table 4. Users who plan to use the second derivative maximum analysis method (an option e.g. on Roche LC480 instruments) are advised to program 45 cycles instead of 40, in order to be able to quantify down to Cq 40.

Table 4. Cycling program for MS11/*CruA* assays.

Step	Stage	T (°C)	Time (s)	Acquisition	Cycles
1	UNG*	50	120	No	1
2	Initial denaturation	95	600	No	1
3	Denaturation	95	15	No	40**
	Amplification Annealing & Extension	60	60	Yes	

*UNG: Uracil-N-glycosylase

** see comment above for users of second derivative maximum analysis method

3.3 Data analysis

After the real-time PCR, analyse the run following the procedure below:

- a) Analyse data: automatic baseline and threshold settings have given the best results at the EURL GMFF.
- b) Save the settings and export all the data for further calculations.

3.4 Calculation of results

After having defined a threshold value within the logarithmic phase of amplification as described above, the instrument's software calculates the C_q values for each reaction.

The standard curves are generated both for the *CruA* and the MS11 specific assays by plotting the C_q values measured for the calibration points against the logarithm of the DNA copy numbers and by fitting a linear regression line into these data.

Thereafter, the standard curves are used to estimate the DNA copy number in the unknown samples.

To obtain the percentage value of event MS11 DNA in the unknown sample, the MS11 copy number is divided by the copy number of the oilseed rape endogenous gene *CruA* and multiplied by 100 ($GM\% = MS11/CruA \times 100$).

4. Equipment and Materials

4.1 Equipment

- Real-time PCR instrument for plastic reaction vessels (glass capillaries are not recommended for the described buffer composition) and appropriate analysis software
- 96-well reaction plates
- Optical caps/adhesion covers
- Microcentrifuge
- Micropipettes
- Standard bench top centrifuge with rotor or standard microfuge fit for 0.5 mL reaction tubes, centrifuge for 96-Well reaction plates
- Vortex
- Racks for reaction tubes, also cooled
- 0.5, 1.5 mL and 5 or 15 mL DNase free reaction tubes

4.2 Reagents

- TaqMan[®] Universal PCR Master Mix. Applied Biosystems Part No 4318157.

4.3 Primers and Probes

Table 5. Primers and probes for the MS11 and *CruA* methods

	Name	DNA Sequence (5' to 3')	Length (nt)
<i>MS11</i>			
Forward primer	SHA086	CAA gAT ggg AAT TAA CAT CTA CAA ATT g	28
Reverse primer	MDB371	gAA ATC CAT gTA AAg CAg CAg gg	23
Probe	TM280	6- FAM-CgA CCA TgT ACA TCC TAC CA -MGB	20
<i>CruA</i>			
Forward primer	MDB510	ggC CAg ggT TTC CgT gAT	18
Reverse primer	MDB511	CCg TCg TTg TAg AAC CAT Tgg	21
Probe	TM458	JOE- AgT CCT TAT gTg CTC CAC TTT CTg gTg CA -BHQ1	29

FAM: 6-carboxyfluorescein; MGB: minor groove binder; JOE: 4,5-dichloro-dimethoxy-fluorescein; BHQ1: black hole quencher.

5. References

1. Plant DNA C-values Database. Royal Botanic Gardens, Kew, <http://data.kew.org/cvalues/>