

# JRC VALIDATED METHODS, REFERENCE METHODS AND MEASUREMENTS REPORT

## **Event-specific Method for the Quantification of Oilseed Rape MS11 Using Real-time PCR Validation Report**

Jacchia S., Sacco M. G.,  
Savini C., Mazzara M., Emons H.

European Union Reference Laboratory for  
Genetically Modified Food and Feed

2019



This publication is a Validated Methods, Reference Methods and Measurements report by the Joint Research Centre (JRC), the European Commission's science and knowledge service. It aims to provide evidence-based scientific support to the European policymaking process. The scientific output expressed does not imply a policy position of the European Commission. Neither the European Commission nor any person acting on behalf of the Commission is responsible for the use that might be made of this publication.

**Contact information**

Name: EURL GMFF

Email: JRC-EURL-GMFF@ec.europa.eu

**JRC Science Hub**

<https://ec.europa.eu/jrc>

JRC 115617

Ispra: European Commission, 2019

© European Union, 2019

Reuse is authorised provided the source is acknowledged. The reuse policy of European Commission documents is regulated by Decision 2011/833/EU (OJ L 330, 14.12.2011, p. 39).

For any use or reproduction of photos or other material that is not under the EU copyright, permission must be sought directly from the copyright holders.

How to cite this report: European Union Reference Laboratory for GM Food and Feed, Joint Research Centre. "Event-specific Method for the Quantification of Oilseed Rape MS11 Using Real-time PCR", 2019. <http://gmo-crl.jrc.ec.europa.eu/StatusOfDossiers.aspx>.

All images © European Union 2019



# Event-specific Method for the Quantification of Oilseed Rape MS11 Using Real-time PCR

## Validation Report

24 January 2019

European Union Reference Laboratory for GM Food and Feed

### Executive Summary

In line with its mandate<sup>a</sup> the European Union Reference Laboratory for GM Food and Feed (EURL GMFF), in collaboration with the European Network of GMO Laboratories (ENGL), validated an event-specific real-time polymerase chain reaction (qPCR) method for detecting and quantifying oilseed rape event MS11 (unique identifier BCS-BNØ12-7). The validation study was conducted according to the EURL GMFF validation procedure [<http://gmo-crl.jrc.ec.europa.eu/guidancedocs.htm>] and the relevant internationally accepted guidelines<sup>(1-5)</sup>.

In accordance with current EU legislation<sup>b</sup>, Bayer CropScience AG provided the detection method and the positive and negative control samples (genomic DNA from leaves of MS11 oilseed rape as positive control DNA, and genomic DNA from leaves of conventional oilseed rape as negative control DNA). The EURL GMFF verified the method performance data provided by the applicant, where necessary experimentally, prepared the validation samples (calibration samples and blind samples at different GM percentage [copies GM/total oilseed rape haploid genome copies]), organised an international collaborative study and analysed the results.

The EURL GMFF in-house verification and the collaborative study confirmed that the method meets the method performance requirements as established by the EURL GMFF and the ENGL, in line with the provisions of Annex III-3.C.2 to Regulation (EU) No 503/2013, and it fulfils the analytical requirements of Regulation (EU) No 619/2011<sup>c</sup>. This validation report is published at <http://gmo-crl.jrc.ec.europa.eu/StatusOfDossiers.aspx>.

<sup>a</sup> Regulation (EC) No 1829/2003 of 22 September 2003 "on genetically modified food and feed".

<sup>b</sup> Regulation (EC) No 503/2013 of 3 April 2013 "on applications for authorisation of genetically modified food and feed in accordance with Regulation (EC) No 1829/2003 of the European Parliament and of the Council and amending Commission Regulations (EC) No 641/2004 and (EC) No 1981/2006".

<sup>c</sup> Regulation (EU) No 619/2011 of 24 June 2011 laying down the methods of sampling and analysis for the official control of feed as regards presence of genetically modified material for which an authorisation procedure is pending or the authorisation of which has expired.

# Content

<b>1. INTRODUCTION.....</b>	<b>4</b>
<b>2. STEP 1 (DOSSIER ACCEPTANCE) AND STEP 2 (SCIENTIFIC DOSSIER ASSESSMENT AND BIOINFORMATICS ANALYSIS).....</b>	<b>4</b>
<b>3. STEP 3 (EXPERIMENTAL TESTING OF THE SAMPLES AND METHOD).....</b>	<b>7</b>
3.1 DNA EXTRACTION.....	7
3.2 METHOD PROTOCOL FOR THE PCR ANALYSIS .....	7
3.3 EURL GMFF EXPERIMENTAL TESTING .....	8
3.3.1 <i>Determination of the zygosity ratio in the positive control sample .....</i>	<i>8</i>
3.3.2 <i>In-house verification of the method performance against ENGL method acceptance criteria .....</i>	<i>9</i>
3.4 INTERNATIONAL COLLABORATIVE STUDY (STEP 4) .....	10
3.4.1 <i>List of participating laboratories .....</i>	<i>10</i>
3.4.2 <i>Real-time PCR equipment used in the study.....</i>	<i>11</i>
3.4.3 <i>Materials used in the international collaborative study .....</i>	<i>12</i>
3.4.4 <i>Design of the collaborative study .....</i>	<i>13</i>
3.4.5 <i>Deviations reported from the protocol .....</i>	<i>13</i>
<b>4. RESULTS .....</b>	<b>14</b>
4.1 EURL GMFF EXPERIMENTAL TESTING .....	14
4.1.1 <i>Zygosity ratio in the positive control sample .....</i>	<i>14</i>
4.1.2 <i>In-house verification of method performance against ENGL method acceptance criteria .....</i>	<i>14</i>
4.2 RESULTS OF THE INTERNATIONAL COLLABORATIVE STUDY.....	17
4.2.1 <i>PCR efficiency and linearity.....</i>	<i>17</i>
4.2.2 <i>GMO quantification.....</i>	<i>18</i>
4.2.3 <i>Method performance requirements.....</i>	<i>19</i>
<b>5. COMPLIANCE OF THE METHOD FOR DETECTION AND QUANTIFICATION OF EVENT MS11 WITH THE REQUIREMENTS OF REGULATION (EU) NO 619/201121</b>	
<b>6. CONCLUSION .....</b>	<b>22</b>
<b>7. REFERENCES.....</b>	<b>22</b>
<b>ANNEX 1: EVENT-SPECIFIC METHOD FOR THE QUANTIFICATION OF OILSEED RAPE MS11 BY REAL-TIME PCR.....</b>	<b>23</b>

## Quality assurance

The EURL GMFF is ISO 17025:2005 accredited [certificate number: Belac 268 TEST (Flexible Scope for DNA extraction, DNA identification and real Time PCR)] and ISO 17043:2010 accredited (certificate number: Belac 268 PT, proficiency test provider).

The original version of the document containing evidence of internal checks and authorisation for publication is archived within the EURL GMFF quality system.

### **Address of contact laboratory:**

European Commission  
Directorate General Joint Research Centre  
Directorate F – Health, Consumers and Reference Materials  
European Union Reference Laboratory for GM Food and Feed  
Food & Feed Compliance (F.5)  
Via E. Fermi, 2749. TP201  
I-21027 Ispra (VA), Italy

Functional mailbox: [JRC-EURL-GMFF@ec.europa.eu](mailto:JRC-EURL-GMFF@ec.europa.eu)

## 1. Introduction

In line with Regulation (EC) No 1829/2003, Bayer CropScience AG provided the EURL GMFF with an event-specific method for detection and quantification of oilseed rape event MS11 (unique identifier BCS-BNØ12-7) together with genomic DNA as positive and negative control samples.

The dossier was found complete (step 1 of the EURL GMFF validation procedure) and the scientific dossier assessment (step 2) concluded that the reported method performance characteristics, assessed against the ENGL method acceptance criteria<sup>d</sup>, allowed moving the method forward to step 3 of the procedure (experimental testing), where the EURL GMFF verified the purity of the control samples provided and conducted an in-house testing of samples and method.

The positive and negative control DNA, originally submitted in accordance with Art 5(3)(j) and Article 17(3)(j) of Regulation (EC) No 1829/2003, were not compliant with the requirements and were replaced, after a request from the EURL GMFF, with suitable controls samples, which were found of good quality.

Step 3 was completed with the conclusion that the method could be submitted to the collaborative study (step 4). This study confirmed that the method is well suited for quantifying genomic DNA of GM oilseed rape MS11, appropriately extracted from food or feed, down to a GM content level of 0.1 % m/m.

The preparation of the report (step 5) was aligned with the timeline communicated by EFSA for its risk assessment.

## 2. Step 1 (dossier acceptance) and step 2 (scientific dossier assessment and bioinformatics analysis)

Documentation and data supplied by the applicant were evaluated by the EURL GMFF for completeness (step 1) and compliance with the ENGL method acceptance criteria (step 2).

The specificity of the event-specific assay was verified by the applicant and confirmed by the EURL GMFF by means of bioinformatics analysis, on the basis of the sequence data provided by the applicant.

---

<sup>d</sup> EURL/ENGL guidance doc "Definition of Minimum Performance Requirements for Analytical Methods of GMO Testing" (<http://gmo-crl.jrc.ec.europa.eu/guidancedocs.htm>)

### Specificity assessment by the applicant

The specificity of the event-specific assay was assessed by the applicant in duplicate real-time PCR reactions, according to the method described, using 200 ng genomic DNA extracted from: oilseed rape Ms1, Ms8, Rf1, Rf2, Rf3, Topas19/2, T45, Oxy-235, RT73, MON 88302, 73496; soybean A2704-12, A5547-127, FG72, 10 % GTS40-3-2, MON 89788, MON 87701, MON 87708, 10 % 356043, 10 % 305423, CV127, MON 87769, MON 87705, 10 % DAS 68416-4, DAS-81419-2, 44406; maize T25, 5 % Bt176, 4.9 % Bt11, 9.9 % MON 810, 4.3 % GA21, 4.9 % NK603, 9.8 % MON 863, 9.9 % 1507, 9.8 % 3272, MIR604, MIR162, 9.9 % 59122, 10 % 98140, MON 88017, MON 89034, 10 % 40278; cotton LLCotton25, T304-40, GHB614, GHB119, MON 1445, MON 531, MON 15985, MON 88913, 281-24-236 x 3006-210-23; rice LLRICE62; non-GM oilseed rape, soybean, maize, cotton and rice. According to the method developer the MS11 method did not react with any sample except the positive control. In addition, an in-silico specificity test was performed against public sequence databases such as NCBI and the patent database: according to the applicant no significant cross-reactivity was detected.

A previously validated oilseed rape-specific PCR method ([http://gmo-crl.jrc.ec.europa.eu/summaries/Ms8\\_validated\\_Method\\_Corrected%20version%201.pdf](http://gmo-crl.jrc.ec.europa.eu/summaries/Ms8_validated_Method_Corrected%20version%201.pdf)), which amplifies a 101 base pair (bp) fragment of the *cruciferin A (CruA)* of *Brassica napus*, was used as a taxon-specific method.

The specificity of the taxon-specific method was assessed by the applicant in duplicate real-time PCR reactions, according to the method described (Tables 1, 2, 3 and 4), using 200 ng genomic DNA extracted from: non-GM oilseed rape, soybean, maize, cotton and rice. According to the method developer the *CruA* method did not react with any sample except the positive control.

### Bioinformatics specificity assessment by the EURL GMFF

The detection method spans the 3' insert-to-plant junction in oilseed rape MS11. The forward primer "SHA086" binds to the insert. The reverse primer "MDB371" binding site was found in the oilseed rape (*Brassica napus*) genomic border adjacent to the insertion. The probe "TM280" binds to the junction between the insert and the 3' genomic region of *Brassica napus*.

The amplicon size is expected to be 124 bp, consistent to what is reported by the applicant. The sequence of the amplicon was analysed by BLAST (NCBI) against local copies of the "nt" and "patents" databases, and no significant similarity was found with any other published sequence except six entries in the Patnt db from Bayer-Aventis on 'Male-sterile brassica plants and methods for producing same'. In addition, the primers were tested against the sequences of the other GMO events present in the Central Core Sequence Information System (CCSIS) of the JRC, as well as the whole genomes of more than 80 plants (including *Brassica rapa*, *Glycine max*, *Oryza sativa*, *Solanum lycopersicum* and *Zea mays*) using the e-PCR prediction tool (NCBI), and no potential amplicon was identified.

### Verification of the ENGL acceptance parameters

The parameters of the calibration curves (slope, R<sup>2</sup> coefficient) were determined by the applicant by quantifying three test samples at different GM levels (see Table 1).

Table 1. Summary of the slope and R<sup>2</sup> values obtained by the applicant

<b>MS11</b>		<b><i>CruA</i></b>	
<b>Slope</b>	<b>R<sup>2</sup></b>	<b>Slope</b>	<b>R<sup>2</sup></b>
-3.3	0.9964	-3.5	0.9967
-3.2	0.9917	-3.5	0.9960
-3.3	0.9963	-3.5	0.9986
-3.2	0.9976	-3.5	0.9985
-3.3	0.9990	-3.5	0.9956
-3.3	0.9868	-3.6	0.9969

According to the ENGL method acceptance criteria, the average value of the slope of the standard curve shall range from -3.1 to -3.6 and the R<sup>2</sup> coefficient shall be  $\geq 0.98$ .

Table 1 indicates that the slope and R<sup>2</sup> coefficient of the standard curves for the GM system (MS11) and the oilseed rape-specific *cruciferin A* (*CruA*) system, as established by the applicant, were within the ENGL acceptance criteria.

Also precision and trueness of the method were established by the applicant and 18 values for each GM level (expressed as mass fraction of GM-material) were provided. Table 2 reports precision and trueness values for the three GM levels as provided by the applicant. Both parameters were within the ENGL acceptance criteria (trueness  $\pm 25$  %, RSD<sub>r</sub>  $\leq 25$  % across the entire dynamic range).

Table 2. Mean %, precision and trueness values provided by the applicant estimated for single measurements

<b>Expected GM %</b>	<b>Test results</b>		
	<b>0.08</b>	<b>0.9</b>	<b>4.5</b>
Measured mean GM %	0.070	0.80	4.54
Precision (RSD <sub>r</sub> %)	20.60	14.29	13.93
Trueness (bias %)	-13.09	-10.85	0.85

\* Numbers are not rounded but are presented as reported by the applicant

### 3. Step 3 (experimental testing of the samples and method)

#### 3.1 DNA extraction

Genomic DNA was isolated from ground oilseed rape seeds, using a "Dellaporta-derived" protocol previously submitted for detection of oilseed rape event Rf3. This DNA extraction method was assessed earlier by the EURL GMFF in order to confirm its performance characteristics. The protocol for DNA extraction and a report on testing were published in 2007 at [http://gmo-crl.jrc.ec.europa.eu/summaries/OSR\\_DNAExtr\\_sampl\\_correctedversion1\\_CRL\\_VL\\_07\\_04.pdf](http://gmo-crl.jrc.ec.europa.eu/summaries/OSR_DNAExtr_sampl_correctedversion1_CRL_VL_07_04.pdf).

In agreement with the ENGL position, endorsing the modularity principle (see also Annex III to Reg. (EU) No 503/2013), and given the similarity in the matrix, the EURL GMFF considers the above mentioned DNA extraction protocol applicable in the context of the validation of the method for oilseed rape event MS11.

Whenever DNA is extracted from more complex and difficult matrices, a thorough control of the quality of the DNA is recommended in order to ensure that it has the required quality for subsequent PCR analysis.

#### 3.2 Method protocol for the PCR analysis

The PCR method provided by the applicant is an event-specific, quantitative, real-time TaqMan<sup>®</sup> PCR method for the determination of the relative content of GM event MS11 DNA to total oilseed rape DNA. The procedure is a simplex system, in which an oilseed rape specific assay targeting the endogenous gene *cruciferin A (CruA)*, and the GM target assay for MS11 are performed in separate wells. The validated method protocol can be found in Annex 1 and is published at <http://gmo-crl.jrc.ec.europa.eu/StatusOfDossiers.aspx>.

For the detection of GM event MS11, a 124 bp fragment of the region spanning the 3' insert-to-plant junction in oilseed rape MS11 is amplified using specific primers. PCR products are measured during each cycle (real-time) by means of a target-specific oligonucleotide probe labelled with FAM (6-carboxyfluorescein) as reporter dye at its 5' end and MBG (minor groove binder) as non-fluorescent quencher dye at its 3' end.

For the relative quantification of GM event MS11, an oilseed rape taxon-specific system amplifies a 101 bp fragment of an oilseed rape *cruciferin A (CruA)* endogenous gene, using *CruA* gene-specific primers and a *CruA* gene-specific probe labelled with JOE as reporter dye at its 5' end and BHQ-1 (Black Hole Quencher<sup>®</sup> 1) as non-fluorescent quencher dye at its 3' end.

Standard curves are generated for both the MS11 and the *CruA* systems by plotting the C<sub>q</sub> values measured for the calibration points against the logarithm of the DNA copy numbers and by fitting a

regression line into these data. Thereafter, the standard curves are used to estimate the copy numbers in the test sample DNA by interpolation from the standard curves.

For the relative quantification of event MS11 DNA in a test sample, the MS11 copy number is divided by the copy number of the oilseed rape haploid genome and multiplied by 100 to obtain the percentage value (GM % = MS11/ oilseed rape haploid genome x 100).

The absolute copy numbers of the calibration curve samples are calculated by dividing the sample DNA mass (nanograms) by the published average 1C value for the oilseed rape genome (1.15 pg)<sup>(6)</sup>. The copy number values used in the quantification, the GMO contents of the calibration samples, and the total DNA quantity used in the PCR reactions are listed in Table 3.

*Note: Numerical values presented in the following tables were rounded keeping two digits for values  $\leq 1$ , one digit for values between 1 and 10 and no digit for values  $\geq 10$ , unless otherwise stated. The calculations in the MS Excel files however were done over not rounded data. This approach might create small inconsistencies in the numerical values reported in the tables but it allows a higher precision in the final results.*

Table 3. Copy number values of the standard curve samples

<b>Sample code</b>	<b>S1</b>	<b>S2</b>	<b>S3</b>	<b>S4</b>	<b>S5</b>
Total amount of oilseed rape DNA in the reaction (ng) *	300	75	19	4.7	1.6
Target taxon haploid genome copies	260870	65217	16304	4076	1359
Target MS11 copies	13043	3261	815	204	68

### **3.3 EURL GMFF experimental testing**

#### **3.3.1 Determination of the zygosity ratio in the positive control sample**

Annex II of Regulation (EU) No 619/2011 requires that "when results are primarily expressed as GM-DNA copy numbers in relation to target taxon-specific DNA copy numbers calculated in terms of haploid genomes, they shall be translated into mass fraction in accordance with the information provided in each validation report of the EURL GMFF." This requires knowledge of the zygosity of the event. In order to satisfy this requirement, the EURL GMFF assessed the zygosity (GM-target to reference target ratio) in the positive control sample submitted by the applicant.

The copy number of the MS11 and of the *CruA* targets in the positive control sample were determined by digital PCR (dPCR) performed on the BioMark HD System using the 12.765 digital arrays (Fluidigm).

Reaction mixes were prepared in order to test the zygosity in five replicates to a final volume of 9  $\mu\text{L}$  and contained 1X TaqMan<sup>®</sup> Universal PCR Master Mix (Applied Biosystems, Cat. number 4318157), 1X GE sample loading reagent (Fluidigm PN 85000746), primers and probes at concentrations indicated in the corresponding validated method (SHA086 and MDB371 primers at 400 nM each, TM280 probe at 200 nM; MDB510 and MDB511 primers at 200 nM each, TM458 probe at 200 nM), and 1  $\mu\text{L}$  of DNA at a concentration of 1.5 ng/ $\mu\text{L}$ ; the DNA concentration was chosen in order to avoid panel saturation (optimal between 200<positive partitions<700).

Loading of the digital chip was performed according to the manufacturer's instructions by using the IFC controller (Fluidigm). Approximately 4.6  $\mu\text{L}$  of the reaction mixes were loaded into each well and distributed into the 765 partitions constituting one panel. The experiment was repeated three times for a total of fifteen data sets for the GM target and fifteen for the reference target. 'No template controls' were included. Amplification conditions were as reported in Annex 1 or in the Validated Method document at <http://gmo-crl.jrc.ec.europa.eu/StatusOfDossiers.aspx>. Data analysis and copy number calculations were performed using the BioMark digital PCR Analysis software. The range of Cq retention was from 15 to 45.

Calculations of means and variances were carried out according to the procedure outlined for random variables in the Annex 4 of the ENGL guidance document 'Verification of analytical methods for GMO testing when implementing inter-laboratory validated methods'<sup>e</sup>.

### **3.3.2 In-house verification of the method performance against ENGL method acceptance criteria**

The method performance characteristics were verified by quantifying on a copy number basis five blind test samples distributed over a range of GM levels (0.05 % - 4.5 %). The experiments were performed on an ABI 7500, an ABI 7900HT and a Roche LC480 real-time platform under repeatability conditions and followed the protocol provided by the applicant. Test samples at GM levels 0.4 %, 0.9 %, 2.0 % and 4.5 % were tested in two real-time PCR runs with two replicates for each GM-level on each plate (total of four replicates per GM-level). The test sample with GM level 0.05 % was tested in 15 replicates in an additional run for each platform. Average values of the slope and of the R<sup>2</sup> coefficient of the standard curves and method trueness and precision over the dynamic range were evaluated against the ENGL method acceptance criteria. On Roche LC480 platform the method was run at 45 cycles as described in the validated method published at <http://gmo-crl.jrc.ec.europa.eu/StatusOfDossiers.aspx> and in Annex 1 below, and analysed with the second derivative maximum method.

---

<sup>e</sup> Verification of analytical methods for GMO testing when implementing inter-laboratory validated methods. European Network of GMO Laboratories (ENGL), 2011.  
<http://gmo-crl.jrc.ec.europa.eu/doc/ENGL%20MV%20WG%20Report%20July%202011.pdf>

### **3.4 International collaborative study (step 4)**

The international collaborative trial involved twelve randomly selected laboratories, all being "national reference laboratories, assisting the EURL GMFF for testing and validation of methods for detection", as listed in annex to Regulation (EC) No 1981/2006 who had expressed their interest in participation. The study was carried out in accordance with the following internationally accepted guidelines:

- The IUPAC "Protocol for the design, conduct and interpretation of method-performance studies" (Horwitz, 1995) <sup>(1)</sup>
- ISO 5725 "Accuracy (trueness and precision) of measurement methods and results", Part 1 and Part 2 (ISO, 1994); ISO 5725-1:1994/Cor 1 (ISO 1998) and ISO 5725-2:1994/Cor 1 (ISO, 2002) <sup>(2-5)</sup>

The objective of the international collaborative study was to verify in experienced laboratories the trueness and precision of the PCR analytical method provided by the applicant and verified in-house by the EURL GMFF.

#### **3.4.1 List of participating laboratories**

The twelve laboratories participating in MS11 international collaborative study were randomly selected from 27 national reference laboratories (NRL) that offered to participate.

Clear guidance was given to the selected laboratories for strictly following the validation protocol that was provided to them. The participating laboratories are listed in Table 4.

Table 4. Laboratories participating in the validation of the detection method for oilseed rape event MS11

Laboratory	Country
Institute for Diagnosis and Animal Health, Molecular Biology and GMOs Unit - National Reference Laboratory for GMOs in food and feed	RO
National Institute of Biology	SI
State Sanitary and Epidemiological Station, Regional Laboratory of Genetically Modified Food	PL
Science and Advice for Scottish Agriculture	UK
Crop Research Institute - Reference Laboratory for GMO Detection and DNA fingerprinting	CZ
BioGEVES - Groupement d'Intérêt Public – Groupe d'Etude et de contrôle des Variétés et des Semences	FR
Istituto Zooprofilattico Sperimentale delle Regioni Lazio e Toscana	IT
Federal Office of Consumer Protection and Food Safety- Berlin	DE
National Veterinary Research Institute in Pulawy, Department of Hygiene of Animal Feedingstuff	PL
National Food and Veterinary Risk Assessment Institute Molecular Biology and GMO Department	LT
State Office for Agriculture, Foodsafety and Fisheries - Mecklenburg Western Pomerania- Rostock	DE
Walloon Agricultural Research Centre - Department Valorization des productions (D4) - Unit 16 - Authentication and traceability	BE

### 3.4.2 Real-time PCR equipment used in the study

Laboratories involved in the collaborative study used the following real-time PCR equipment: five laboratories used ABI 7900HT, three used ABI 7500, one used ABI 7500 fast, one Roche LC 480, one ABI QuantStudio 6 and one used Stratagene Mx3005P.

This variability of equipment, with its known potential influence on PCR results, reflects the real-life situation in the control laboratories and provides additional assurance that the method is robust and usable under real conditions.

### 3.4.3 Materials used in the international collaborative study

For the validation of the quantitative event-specific method, test samples were provided by the EURL GMFF to the participating laboratories. They were derived from:

- i) genomic DNA extracted by the applicant from hemizygous oilseed rape leaves harbouring the MS11 event, and
- ii) genomic DNA extracted by the applicant from conventional oilseed rape leaves genetically similar to those harbouring the MS11 event.

The control samples were prepared by the EURL GMFF from the genomic DNA provided by the applicant in accordance to Regulation (EC) No 1829/2003, Art 2.11<sup>f</sup>.

The control samples were used by the EURL GMFF to prepare standards (of known GMO content) and test samples (of undisclosed GM content = blind samples) by mixing MS11 oilseed rape DNA and non-GM oilseed rape DNA.

The calibration sample S1 was prepared by mixing the appropriate amount of MS11 DNA with control non-GM oilseed rape DNA to obtain a 5 % (in copy number ratio related to haploid genome copies) GM sample. Calibration samples S2-S4 were prepared by 4.0-fold serial dilutions from the S1 sample and sample S5 by a 3.0-fold serial dilution from the S4 sample.

The twelve NRLs participating in the validation study received the following materials:

- ✓ Five calibration samples with known concentrations of GM-event (175 µL of DNA solution each) labelled from S1 to S5 (Table 3).
- ✓ Twenty blinded test DNA samples (87.5 µL of DNA solution, each at 40 ng/µL) labelled from U1 to U20, representing five GM levels, each in four replicates (Table 5)

Table 5. MS11 blinded samples GM % contents

MS11 GM %
GM copy number/oilseed rape haploid genome copy number x 100
4.5
2.0
0.90
0.40
0.05

<sup>f</sup> Control sample defined as the GMO or its genetic material (positive sample) and the parental organism or its genetic material that has been used for the purpose of the genetic modification (negative sample).

- ✓ Reaction reagents:
  - TaqMan<sup>®</sup> Universal PCR Master Mix (2x), one vial: 8 mL
  - distilled sterile water, one vial: 4 mL
  
- ✓ Primers and probes (1 tube each) as follows:
  - CruA* taxon-specific assay
    - MDB510 primer (10 µM): 160 µL
    - MDB511 primer (10 µM): 160 µL
    - TM458 probe (10 µM): 135 µL
  
  - MS11 assay
    - SHA086 primer (10 µM): 320 µL
    - MDB371 primer (10 µM): 320 µL
    - TM280 probe (10 µM): 160 µL

#### 3.4.4 Design of the collaborative study

Participating laboratories received a detailed validation protocol that included the exact design of the PCR plates, ensuring that on each PCR plate the samples were analysed for the MS11 event-specific system and for the *CruA* taxon-specific system. In total, two plates were run by each participating laboratory.

The laboratories prepared the PCR master-mixes for the oilseed rape event MS11 and the *CruA* assay in accordance with the description provided in the validation protocol. Calibration and test samples were loaded on the PCR plates as per pre-determined plate layout.

The amplification reaction followed the cycling program specified in the protocol. Participants determined the GM % in the test samples according to the instructions and also reported the raw data to the EURL GMFF on an Excel sheet that was designed, validated and distributed by the EURL GMFF. All data are stored by the EURL GMFF on a dedicated and protected server.

The EURL GMFF analysed the data against the parameters and the limits set by the ENGL, i.e. trueness, precision, amplification efficiency and linearity.

#### 3.4.5 Deviations reported from the protocol

Ten laboratories reported no deviations from the validation protocol. One laboratory reported an exchange in position between two blind samples in one plate; one laboratory reported mixing reaction mix components at room temperature instead of on ice, and using 1.5 mL tubes instead of 0.5 mL tubes when adding the samples to the reaction mix.

## 4. Results

### 4.1 EURL GMFF experimental testing

#### 4.1.1 Zygoty ratio in the positive control sample

The results of the digital PCR analysis conducted by the EURL GMFF on the MS11 and *CruA* targets to determine the zygoty ratio in the positive control samples are shown in Table 6.

Table 6. Zygoty ratio of the MS11 and *CruA* targets in the positive control sample.

Mean ratio (MS11/ <i>CruA</i> )	0.22
Standard deviation	0.02
RSD <sub>r</sub> (%)	8.6
Standard error of the mean	0.005
Upper 95 % CI of the mean	0.23
Lower 95 % CI of the mean	0.21

The mean ratio (MS11/*CruA*) equals 0.22. Recently, thanks to more advanced bioinformatics tools, it became evident that the *CruA* target is present in two copies per haploid genome. Despite the 95 % confidence interval (CI) does not include 0.25 for an alpha = 0.05, the mean ratio can be considered equivalent to the expected ratio of 0.025 for an oilseed rape heterozygous control sample and assuming a double-copy endogenous gene target. The zygoty of the positive control sample provides the evidence for a hemizygous status of the MS11 event in *B. napus*. The samples of the standard curve and test samples have been prepared as ratio of GM copies to haploid genome copy numbers.

Therefore, for the control sample used in the present study, results expressed as GM-DNA copy numbers in relation to target taxon specific DNA copy numbers can be translated into mass fraction as follows:

$$0.05 \text{ GM \% in DNA copy number ratio related to haploid genome copy numbers} = 0.1 \text{ GM \% in mass fraction}$$

#### 4.1.2 In-house verification of method performance against ENGL method acceptance criteria

Test samples with GM levels from 0.40 % to 4.5 % were tested in three real-time PCR runs with two replicates for each GM-level on each plate (total of four replicates per GM-level). The sample at 0.05 % GM level was tested for precision in quantification in 15 replicates in separate runs.

Tests were conducted on ABI 7500, ABI 7900 and Roche LC480 for robustness.

The standard curve parameters and the results of efficiency, linearity, trueness and precision obtained in the three real-time PCR runs with the test samples are shown in Tables 7A, 7B, 8, 9 and 10.

According to the ENGL method acceptance criteria, the average value of the slope of the standard curve shall range from -3.1 to -3.6 and the R<sup>2</sup> coefficient shall be ≥ 0.98. Table 7A and 7B document that the slopes of the standard curves and the R<sup>2</sup> coefficients were within the limits established by the ENGL. The EURL GMFF in-house results confirm the data provided by the applicant.

Table 7A. Standard curve parameters of the real-time PCR tests, carried out on ABI 7900HT, ABI 7500, and Roche LC480 to quantify GM levels in the range 0.40 % to 4.5 % in four replicates each. Slope and R<sup>2</sup> coefficient values were rounded to two digits.

	<b>MS11 system</b>			<b>CruA system</b>		
	Slope	PCR efficiency*	R <sup>2</sup>	Slope	PCR efficiency*	R <sup>2</sup>
Run A	-3.23	104	0.99	-3.42	96	0.99
Run B	-3.36	99	0.99	-3.34	99	0.99
Run C	-3.37	98	1.00	-3.40	97	1.00
Run D	-3.36	99	1.00	-3.43	96	1.00
Run E	-3.32	100	0.99	-3.46	95	1.00
Run F	-3.30	101	0.99	-3.42	96	1.00

\* PCR efficiency (%) is calculated using the formula  $\text{Efficiency} = (10^{(-1/\text{slope})} - 1) \times 100$

Runs A-B were carried out on ABI 7900HT; Runs C-D were carried out on ABI 7500; Runs E and F were carried out on Roche LC480.

Table 7B. Standard curve parameters of the real-time PCR tests, carried out on ABI 7500, ABI 7900, and Roche LC480 to quantify the GM level 0.05 % in 15 replicates. Slope and R<sup>2</sup> coefficient values were rounded to two digits.

	<b>MS11 system</b>			<b>CruA system</b>		
	Slope	PCR efficiency*	R <sup>2</sup>	Slope	PCR efficiency*	R <sup>2</sup>
Run G	-3.41	96	0.99	-3.42	96	1.00
Run H	-3.37	98	1.00	-3.45	95	1.00
Run I	-3.24	104	0.99	-3.44	95	1.00

\* PCR efficiency (%) is calculated using the formula  $\text{Efficiency} = (10^{(-1/\text{slope})} - 1) \times 100$

Run G was carried out on ABI 7900HT; Run H was carried out on ABI 7500; Run I was carried out on Roche LC480.

According to the ENGL method acceptance criteria the method trueness (measured as bias in % of the target GM level) should be within ± 25 % of the accepted reference value over the entire

dynamic range and the precision, expressed as  $RSD_r$  % (relative standard deviation of repeatability), should be  $\leq 25$  %, also over the entire dynamic range.

Tables 8, 9 and 10 show that trueness and precision of quantification were within the limits established by the ENGL for the PCR machines used.

Table 8. Values of trueness and precision as established by the EURL GMFF in its in-house verification using an ABI 7900HT. GM % in copy/copy haploid genomes.

<b>Target GM-levels %</b>	<b>Measured GM-level %</b>	<b>Bias % of the target GM-level</b>	<b>Precision (<math>RSD_r</math> %)</b>
4.5	4.8	7.1	13
2.0	2.0	0.65	6.8
0.90	0.89	-1.3	3.6
0.40	0.34	-16	9.8
0.05	0.05	1.9	18

Table 9. Values of trueness and precision as established by the EURL GMFF in its in-house verification using an ABI 7500. GM % in copy/copy haploid genomes.

<b>Target GM-levels %</b>	<b>Measured GM-level %</b>	<b>Bias % of the target GM-level</b>	<b>Precision (<math>RSD_r</math> %)</b>
4.5	4.3	-5.2	7.3
2.0	1.9	-4.9	2.5
0.90	0.78	-14	3.6
0.40	0.36	-11	5.3
0.05	0.05	-4.6	12

Table 10. Values of trueness and precision as established by the EURL GMFF in its in-house verification using a Roche LC480. GM % in copy/copy haploid genomes.

<b>Target GM-levels %</b>	<b>Measured GM-level %</b>	<b>Bias % of the target GM-level</b>	<b>Precision (<math>RSD_r</math> %)</b>
4.5	4.7	4.2	10
2.0	1.9	-3.4	4.5
0.90	0.82	-8.6	7.3
0.40	0.35	-13	6.9
0.05	0.04	-20	17

## 4.2 Results of the international collaborative study

### 4.2.1 PCR efficiency and linearity

The PCR efficiency (%) and  $R^2$  values (expressing the linearity of the regression) for the standard curve, reported by participating laboratories are displayed in Table 11. The PCR efficiency (%) was calculated from the standard curve slopes using the formula:

$$\text{Efficiency} = (10 * (-1/\text{slope})) - 1) \times 100$$

Table 11 indicates that the efficiency of amplification for the MS11 system ranges from 92 % to 111 % and the linearity from 0.97 to 1.00; the amplification efficiency for the oilseed rape-specific system ranges from 91 % to 102 % and the linearity from 0.99 to 1.00. The mean PCR efficiency was 100 % for MS11 assay and 95 % for the *CruA* one. The average  $R^2$  of the methods was 0.99 and 1.00 for the MS11 and *CruA* assays, respectively. Both PCR efficiency and linearity values were within the ENGL acceptance criteria.

Table 11. Values of slope, PCR efficiency and  $R^2$  obtained during the international collaborative trial. Slope and  $R^2$  coefficient values were rounded to two digits.

Lab	Plate	MS11			CruA		
		Slope	PCR Efficiency (%)	$R^2$	Slope	PCR Efficiency (%)	$R^2$
1	A	-3.27	102	0.99	-3.47	94	1.00
	B	-3.34	99	1.00	-3.46	94	1.00
2	A	-3.53	92	0.99	-3.45	95	1.00
	B	-3.45	95	0.99	-3.46	94	1.00
3	A	-3.33	99	1.00	-3.42	96	1.00
	B	-3.36	98	1.00	-3.45	95	1.00
4	A	-3.42	96	1.00	-3.55	91	1.00
	B	-3.43	96	1.00	-3.47	94	1.00
5	A	-3.34	99	0.99	-3.49	93	1.00
	B	-3.33	100	0.99	-3.45	95	1.00
6	A	-3.24	104	0.97	-3.39	97	0.99
	B	-3.09	111	0.97	-3.42	96	0.99
7	A	-3.32	100	1.00	-3.54	91	1.00
	B	-3.38	98	1.00	-3.42	96	1.00
8	A	-3.25	103	0.99	-3.54	92	1.00
	B	-3.38	97	1.00	-3.52	92	1.00
9	A	-3.30	101	1.00	-3.40	97	1.00
	B	-3.32	100	1.00	-3.39	97	1.00
10	A	-3.40	97	1.00	-3.46	94	1.00
	B	-3.43	96	0.99	-3.48	94	1.00

<b>11</b>	A	-3.22	104	0.98	-3.31	100	0.99
	B	-3.17	107	0.97	-3.28	102	0.99
<b>12</b>	A	-3.34	99	0.99	-3.40	97	1.00
	B	-3.34	99	1.00	-3.39	97	1.00
	<b>Mean</b>	-3.33	100	0.99	-3.44	95	1.00

These results confirm the appropriate performance characteristics of the methods tested in terms of efficiency and linearity.

#### 4.2.2 GMO quantification

Table 12 reports the values of quantification for the four replicates of each GM level as reported by each of the twelve participating laboratories.

Table 12. GM % values determined by laboratories for test samples


LAB	GMO content (%) *																			
	0.05				0.40				0.90				2.0				4.5			
	REP 1	REP 2	REP 3	REP 4	REP 1	REP 2	REP 3	REP 4	REP 1	REP 2	REP 3	REP 4	REP 1	REP 2	REP 3	REP 4	REP 1	REP 2	REP 3	REP 4
<b>1</b>	0.04	0.04	0.04	0.04	0.34	0.38	0.32	0.42	0.86	0.91	0.91	0.92	2.0	1.7	2.2	2.0	4.5	4.2	5.3	4.8
<b>2</b>	0.05	0.05	0.06	0.05	0.38	0.38	0.36	0.45	0.88	0.93	0.94	0.95	2.0	1.9	2.0	2.0	4.3	4.4	4.3	4.4
<b>3</b>	0.05	0.04	0.05	0.04	0.31	0.38	0.29	0.29	0.80	0.84	0.75	0.77	1.7	1.8	1.8	1.8	4.3	3.8	4.7	5.0
<b>4</b>	0.06	0.06	0.05	0.06	0.29	0.37	0.42	0.40	0.92	0.93	0.90	1.01	2.0	2.1	2.2	2.1	4.7	4.7	4.8	5.0
<b>5</b>	0.04	0.05	0.05	0.05	0.40	0.39	0.36	0.42	0.92	0.91	0.94	0.97	1.9	2.1	2.1	2.2	4.4	4.5	5.0	4.6
<b>6</b>	0.05	0.06	0.07	0.05	0.50	0.51	0.58	0.48	0.87	1.26	0.81	1.12	2.7	2.7	2.3	2.3	6.1	6.4	6.0	7.2
<b>7</b>	0.06	0.04	0.06	0.06	0.37	0.46	0.41	0.40	0.85	0.93	0.88	0.86	2.1	2.0	2.0	2.0	4.3	4.7	4.9	5.0
<b>8</b>	0.05	0.06	0.04	0.06	0.36	0.41	0.40	0.36	0.96	0.93	0.83	0.87	2.2	2.4	2.1	2.0	5.1	5.2	5.2	4.7
<b>9</b>	0.05	0.05	0.05	0.04	0.36	0.38	0.35	0.36	0.86	0.87	0.85	0.83	1.9	1.8	1.9	1.9	4.4	4.7	4.9	4.5
<b>10</b>	0.05	0.05	0.05	0.05	0.38	0.38	0.35	0.34	0.90	0.79	0.89	0.81	1.9	1.9	2.1	2.0	4.1	4.4	4.5	4.4
<b>11</b>	0.04	0.03	0.05	0.03	0.25	0.34	0.40	0.28	0.92	0.92	0.80	0.94	1.8	1.9	2.6	2.0	5.9	6.0	6.7	4.4
<b>12</b>	0.04	0.04	0.03	0.04	0.34	0.22	0.27	0.37	0.70	0.72	0.89	0.84	1.6	1.6	1.8	2.2	3.6	4.3	4.9	4.5

\* GMO % = (GMO copy number/oilseed rape haploid genome copy number) x 100

n.a. not available

A graphical representation of the data reported in Table 12 is provided in Figure 1 that shows the relative deviation from the true value for each GM level tested for the participating laboratory. The coloured bars represent the deviation of the GM level measured in % of the true GM level; the green bar on the right represents the mean relative deviation over all data before eliminating outliers for each GM level.

Figure 1. Relative deviation (%) from the true value of GM level \*



\* For PL10 at level 2.0 % a very small relative deviation from the true value was observed and therefore the corresponding histogram does not show up in Figure 1. PL = participating laboratory.

Overall, most laboratories' mean relative deviations from the true values were within a maximum of  $\pm 25\%$ . At GM level 0.05 % and 0.40 % eleven laboratories were within the limit; at GM level 2.0 % and 0.9 % twelve laboratories; at GM level 4.5 % ten laboratories were within the limit. One laboratory underestimated GM level 0.05 % by more than 25 %, one laboratory overestimated GM level 0.40 % by more than 25 %, and two laboratories overestimated GM level 4.5 % by more than 25 %. Results from one laboratory (PL06) were identified as outlier with Cochran's test for GM level 0.9 % and with Double Grubbs' test for GM level 4.5 %; results from another laboratory (PL11) were identified as outlier with Cochran's test for GM level 4.5 %. No clear trend for over- or underestimation was observed.

The mean bias generated by all laboratories ranged between -6.4 % and 2.2 %. The method is within the accepted limits established by the ENGL ( $\pm 25\%$  over the entire dynamic range).

#### 4.2.3 Method performance requirements

Among the performance requirements established by ENGL and adopted by the EURL GMFF (<http://gmo-crl.jrc.ec.europa.eu/guidancedocs.htm>), repeatability and reproducibility are to be assessed through an international collaborative trial. Table 13 illustrates the estimation of repeatability and reproducibility at the various GM levels tested during the study (see Table 4 for a list of the participant laboratories).

According to the ENGL method performance requirements the relative reproducibility standard deviation ( $RSD_R$ ), that describes the inter-laboratory variation, should be below 35 % at the target

concentration and over the majority of the dynamic range, while it should be below 50 % at the lower end of the dynamic range.

As it can be observed in Table 13, the method satisfies this requirement at all GM levels tested. In fact, the highest value of  $RSD_R$  % is 19 % at the 0.05 % GM level, thus within the acceptance criterion.

Table 13. Summary of validation results for the MS11 method, expressed as GM copy numbers in relation to target taxon haloid genome copy numbers.

	Test Sample Expected GMO %				
	0.05	0.40	0.90	2.0	4.5
Laboratories having returned valid results	12	12	12	12	12
Samples per laboratory	4	4	4	4	4
Number of outliers	0	0	1	0	2
Reason for exclusion	-	-	C	-	C/DG
Mean value	0.05	0.37	0.88	2.0	4.6
Relative repeatability standard deviation, $RSD_r$ (%)	12	12	5.6	8.3	7.4
Repeatability standard deviation	0.01	0.04	0.05	0.17	0.34
Relative reproducibility standard deviation, $RSD_R$ (%)	19	18	7.7	12	8.2
Reproducibility standard deviation	0.01	0.07	0.07	0.24	0.38
Bias** (absolute value)	0.00	-0.03	-0.02	0.03	0.10
Bias (%)	-3.4	-6.4	-2.6	1.6	2.2

\* C= Cochran's test; DG = Double Grubbs' test; identification and removal of outliers through Cochran and Grubbs tests, according to ISO 5725-2.

\*\* Bias is estimated according to ISO 5725 data analysis protocol.

Table 13 also documents the relative repeatability standard deviation ( $RSD_r$ ) estimated for each GM level. In order to accept methods for collaborative study, the EURL GMFF and ENGL require that the  $RSD_r$  value indicated by the applicant and confirmed by the EURL GMFF through in-house experiments, is below 25 % (see ENGL document "Definition of Minimum Performance Requirements for Analytical Methods of GMO Testing" <http://gmo-crl.jrc.ec.europa.eu/guidancedocs.htm>). As it can be observed from the values reported, the repeatability standard deviation is below 25 % at all GM levels, with the highest value of 12 % at the 0.05 % and the 0.40 % GM level.

The trueness of the method is estimated using the measures of the method bias for each GM level. According to ENGL method performance requirements, trueness should be  $\pm 25$  % across the entire dynamic range. The method satisfies this requirement across the dynamic range tested, with the highest absolute value of bias (%) of -6.4 % at the 0.40 % GM level.

## 5. Compliance of the method for detection and quantification of event MS11 with the requirements of Regulation (EU) No 619/2011

To verify the compliance of the method under validation with the requirements of Regulation (EU) No 619/2011, the following steps were carried out and their outcome is summarised in Table 14:

- at step 2 of the validation process (scientific assessment of the dossier), the EURL GMFF acknowledged that the  $RSD_r$  value at the 0.08 % level shown by the applicant's dossier (expressed as mass fraction of GM-material) was 20.60 %, based on 18 replicates (Table 2), hence below the maximum value of 25 % required by the ENGL. The EURL GMFF therefore concluded that it could accept the applicant's data on method performance;
- at step 3 of the validation process (experimental testing of samples and methods), the EURL GMFF determined the  $RSD_r$  % value at the level of 0.1 % in mass fraction of GM-material (corresponding to 0.05 % expressed in terms of copy number ratio to haploid genome copy numbers). The experiments were carried out under repeatability conditions on fifteen replicates. The  $RSD_r$  resulted to range between 12 % and 18 % (Table 8, 9 and 10) depending on the qPCR platform applied, hence also below 25 %;
- the collaborative study (step 4 of the validation process) established that over the twelve participating laboratories at the level of 0.1 % related to mass fraction of GM material the  $RSD_r$  of the method was 12 %, therefore also below 25 % and in line with the previous data.

The outcome of the different steps is summarised Table 14.

Table 14. Precision of the event-specific method for quantitative detection of MS11 at or around 0.1 % level related to mass fractions of GM material.

Source	$RSD_r$ %	GM %
Applicant's method optimisation	20.60 %	0.08 %
EURL GMFF tests	12 - 18 %	0.1 %
Collaborative study	12 %	0.1 %

Based on the results of the EURL GMFF in-house verification and of the international collaborative study, it is concluded that the method  $RSD_r$  % is lower than 25 % at the level of 0.1 % related to mass fraction of GM material, hence the method meets the requirement laid down in Regulation (EU) No 619/2011.

## 6. Conclusion

The method provided by the applicant has been validated in accordance to the EURL GMFF validation process, respecting all requirements of the relevant EU legislation and international standards for method validation.

This validation study confirmed that the method is applicable to the control samples provided by the applicant (see paragraph 3.4.3), in accordance with the requirements of Annex I-3.C.2 to Commission Regulation (EU) No 503/2013 and (EU) No 619/2011 and meets all method performance requirements established by the ENGL and the EURL GMFF. The method is therefore valid to be used for regulatory purposes, including the quantification of low level presence of 0.1 % (m/m) of the GM event. It can be assumed that it is applicable to any appropriately extracted oilseed rape genomic DNA.

In any case the user of the method is advised to verify the quality of the extracted genomic DNA in order to ensure that it is suitable for the subsequent PCR analysis. This is particularly relevant for more complex matrices of samples from food and feed products.

The validated method is described in detail as "Validated Method" at <http://gmo-crl.jrc.ec.europa.eu/StatusOfDossiers.aspx> and in Annex 1.

## 7. References

1. Horwitz W. Protocol for the design, conduct and interpretation of method- performance studies, *Pure and Appl. Chem.* 1995; 67: 331-343.
2. International Standard (ISO) 5725-1, 1994. Accuracy (trueness and precision) of measurement methods and results. Part 1: General principles and definitions. International Organization for Standardization, Genève, Switzerland.
3. ISO 5725-1:1994/Cor 1:1998.
4. International Standard (ISO) 5725-2, 1994. Accuracy (trueness and precision) of measurement methods and results. Part 2: Basic method for the determination of repeatability and reproducibility of a standard measurement method. International Organization for Standardization, Genève, Switzerland.
5. ISO 5725-2:1994/Cor 1:2002.
6. Plant DNA C-values Database, <http://data.kew.org/cvalues/>

# **Annex 1: Event-specific Method for the Quantification of oilseed rape MS11 by Real- time PCR**

## **Validated Method**

### **Method development:**

Bayer CropScience AG

## 1. General information and summary of the methodology

This protocol describes an event-specific real-time quantitative TaqMan<sup>®</sup> PCR (polymerase chain reaction) procedure for the determination of the relative content of oilseed rape event MS11 DNA to total oilseed rape DNA in a sample.

Template DNA extracted by means of suitable methods should be tested for quality and quantity prior to use in PCR assays. Tests for the presence of PCR inhibitors (e.g. monitor run of diluted series, use of DNA spikes) are also recommended to ensure suitability of the extracted DNA.

For the detection of GM event MS11, a 124 bp fragment of the region spanning the 3' insert-to-plant junction in oilseed rape MS11 is amplified using specific primers. PCR products are measured during each cycle (real-time) by means of a target-specific oligonucleotide probe labelled with FAM (6-carboxyfluorescein) as reporter dye at its 5' end and MBG (minor groove binder) as non-fluorescent quencher dye at its 3' end.

For the relative quantification of GM event MS11, an oilseed rape taxon-specific system amplifies a 101 bp fragment of an oilseed rape *cruciferin A* (*CruA*) endogenous gene (Accession number, GeneBank: X14555), using *CruA* gene-specific primers and a *CruA* gene-specific probe labelled with JOE as reporter dye at its 5' end and BHQ-1 (Black Hole Quencher<sup>®</sup> 1) as non-fluorescent quencher dye at its 3' end

The measured fluorescence signal passes a threshold value after a certain number of cycles. This threshold cycle is called the "Cq" value. For quantification of the amount of MS11 DNA in a test sample, Cq values for the MS11 and the *CruA* systems are determined for the sample. Standard curves are then used to estimate the relative amount of MS11 DNA to total oilseed rape DNA.

## 2. Validation and performance characteristics

### 2.1 General

The method was optimised for suitable DNA extracted from genetically modified and conventional oilseed rape seeds and grain. Precision and trueness of the method were tested through an international collaborative ring trial using DNA samples at different GM contents.

### 2.2 Collaborative trial

The method was validated in an international collaborative study by the European Union Reference Laboratory for GM Food and Feed (EURL GMFF). The study was undertaken with twelve participating laboratories in November 2017.

A detailed validation report can be found at <http://gmo-crl.jrc.ec.europa.eu/StatusOfDossiers.aspx>.

### **2.3 Limit of detection (LOD)**

According to the method developer, the relative LOD of the method is at least 0.023 % (related to mass fraction of GM material) in 200 ng of total suitable oilseed rape DNA. The relative LOD was not assessed in the collaborative study.

### **2.4 Limit of quantification (LOQ)**

According to the method developer, the relative LOQ of the method is at least 0.08 % (related to mass fraction of GM material) in 200 ng of total suitable oilseed rape DNA. The lowest relative GM content of the target sequence included in the collaborative trial was 0.1 % (mass fraction of GM-material).

### **2.5 Molecular specificity**

The method exploits a unique DNA sequence in the region spanning the 3' insert-to-plant junction in oilseed rape MS11 and is therefore event-specific for the event MS11. This was confirmed in the validation study.

## **3. Procedure**

### **3.1 General instructions and precautions**

- The procedures require experience of working under sterile conditions.
- Laboratory organisation, e.g. "forward flow direction" during PCR-setup, should follow international guidelines, e.g. ISO 24276:2006.
- PCR reagents should be stored and handled in a separate room where no nucleic acids (with exception of PCR primers or probes) or DNA degrading or modifying enzymes have been handled previously. All handling of PCR reagents and controls requires dedicated equipment, especially pipettes.
- All the equipment should be sterilised prior to use and any residue of DNA should have been removed. All material used (e.g. vials, containers, pipette tips, etc.) must be suitable for PCR and molecular biology applications. They must be DNase-free, DNA-free, sterile and unable to adsorb protein or DNA.
- Filter pipette tips protected against aerosol should be used.
- Powder-free gloves should be used and changed regularly

- Laboratory benches and equipment should be cleaned periodically, with 10% sodium hypochlorite solution (bleach).
- Pipettes should be checked regularly for precision and calibrated, if necessary.
- All handling steps, unless specified otherwise, should be carried out at room temperature.
- In order to avoid repeated freeze/thaw cycles aliquots should be prepared.

## 3.2 Real-time PCR for quantitative analysis of oilseed rape event MS11

### 3.2.1 General

The real-time PCR set-up for the taxon (*CruA*) and the GMO (event MS11) target sequences are carried out in separate vials. Multiplex qPCR (using differential fluorescent labels for the probes) has not been tested or validated by the EURL GMFF.

The method is developed for a total volume of 25  $\mu$ L per reaction mixture for the GM (event MS11) and the taxon (*CruA*) assay with the reagents as listed in Table 2 and Table 3.

### 3.2.2 Calibration

The calibration curves have to be established on at least five samples. The first point of the calibration curve (S1) should be established for a sample containing 5 % oilseed rape MS11 DNA in a total of 300 ng of oilseed rape DNA (corresponding to 260870 oilseed rape haploid genome copies with one haploid genome assumed to correspond to 1.15 pg of oilseed rape genomic DNA)<sup>(1)</sup>. Standards S2 to S5 are to be prepared by serial dilutions (dilution factor 4.0 for samples S2-S4 and dilution factor 3.0 for standard S5) according to Table 1 below.

Table 1. Copy number values of the standard curve samples

Sample code	S1	S2	S3	S4	S5
Total amount of oilseed rape DNA in reaction (ng)*	300	75	18.8	4.7	1.6
oilseed rape haploid genome copies	260870	65217	16304	4076	1359
MS11 copies	13043	3261	815	204	68

\* Total nanograms are rounded to the integral value

A calibration curve is to be produced by plotting the C<sub>q</sub> values against the logarithm of the target copy number for the calibration points. This can be done by means of spreadsheet software, e.g. Microsoft Excel, or directly by options available with the software.

The copy number measured for each unknown sample DNA is obtained by interpolation from the standard curves.

### 3.2.3 Real-time PCR set-up

1. Thaw, mix and centrifuge the components needed for the run. Keep thawed reagents on ice.
2. In two tubes on ice, add the components in the order mentioned below (except DNA) to prepare the reaction mixes for the MS11 oilseed rape specific system (Table 2) and the *CruA* reference gene system (Table 3). Please note that additional volume is included in the total to cover pipetting variability due to the viscosity of the solution.

Table 2. Amplification reaction mixture in the final volume/concentration per reaction well for the MS11 assay.

Component	Final concentration	µL/reaction
TaqMan <sup>®</sup> Universal PCR Master Mix (2x)	1x	12.5
SHA086 (10 µM)	400 nM	1.00
MDB371 (10 µM)	400 nM	1.00
TM280* (10 µM)	200 nM	0.50
Nuclease free water	-	5.0
DNA	-	5.0
Total reaction volume:		25 µL

\*TaqMan<sup>®</sup> probe labelled with 6-FAM at its 5'-end and MGB at its 3'-end

Table 3. Amplification reaction mixture in the final volume/concentration per reaction well for the oilseed rape *CruA* assay.

Component	Final concentration	µL/reaction
TaqMan <sup>®</sup> Universal PCR Master Mix (2x)	1x	12.5
MDB510 (10 µM)	200 nM	0.50
MDB511 (10 µM)	200 nM	0.50
TM458* (10 µM)	200 nM	0.50
Nuclease free water	-	6.0
DNA	-	5.0
Total reaction volume:		25 µL

\*TaqMan<sup>®</sup> probe is labelled with JOE at its 5'-end and BHQ1 at its 3'-end

3. Mix well and centrifuge briefly.
4. Prepare two 0.5 mL reaction tubes (one for the oilseed rape MS11 and one for the *CruA* system) for each DNA sample to be tested (standard curve samples, unknown samples and control samples).
5. Add into each reaction tube the amount of reaction mix for 3.5 PCR repetitions (70 µL for the MS11 oilseed rape system and 70 µL for the *CruA* system). Add to each tube the correct amount of DNA for 3.5 PCR repetitions (17.5 µL DNA). The volume for the additional 0.5 repetition will ensure adequate volume when loading the samples. Vortex each tube for approx. 10 seconds. This step is mandatory to reduce to a minimum the variability among the repetitions of each sample.
6. Spin down the tubes. Aliquot 25 µL for MS11 system and for the *CruA* reference system in each well.
7. Place an optical cover on the reaction plate and briefly centrifuge the plate.
8. Place the reaction plate in the real-time PCR apparatus (possibly apply a compression pad, depending on the model), according to the manufacturer's instructions and your Standard Operating Procedures and start the run.
9. Select FAM as reporter dye for the MS11 and JOE for the *CruA* reference system. Define MGB or non-fluorescent as quencher dye for MS11 specific system and BHQ or non-fluorescent for *CruA* reference system. Select ROX as the passive reference dye. Enter the correct reaction volume (25 µL).
10. Run the PCR with the cycling program described in Table 4. Users who plan to use the second derivative maximum analysis method (an option e.g. on Roche LC480 instruments) are advised to program 45 cycles instead of 40, in order to be able to quantify down to Cq 40.

Table 4. Cycling program for MS11/*CruA* assays.

Step	Stage	T (°C)	Time (s)	Acquisition	Cycles
1	UNG*	50	120	No	1
2	Initial denaturation	95	600	No	1
3	Denaturation	95	15	No	40**
	Amplification Annealing & Extension	60	60	Yes	

\*UNG: Uracil-N-glycosylase

\*\* see comment above for users of second derivative maximum analysis method

### 3.3 Data analysis

After the real-time PCR, analyse the run following the procedure below:

- a) Analyse data: automatic baseline and threshold settings have given the best results at the EURL GMFF.
- b) Save the settings and export all the data for further calculations.

### 3.4 Calculation of results

After having defined a threshold value within the logarithmic phase of amplification as described above, the instrument's software calculates the C<sub>q</sub> values for each reaction.

The standard curves are generated both for the *CruA* and the MS11 specific assays by plotting the C<sub>q</sub> values measured for the calibration points against the logarithm of the DNA copy numbers and by fitting a linear regression line into these data.

Thereafter, the standard curves are used to estimate the DNA copy number in the unknown samples.

To obtain the percentage value of event MS11 DNA in the unknown sample, the MS11 copy number is divided by the copy number of the oilseed rape endogenous gene *CruA* and multiplied by 100 (GM% = MS11/*CruA* x 100).

## 4. Equipment and Materials

### 4.1 Equipment

- Real-time PCR instrument for plastic reaction vessels (glass capillaries are not recommended for the described buffer composition) and appropriate analysis software
- 96-well reaction plates
- Optical caps/adhesion covers
- Microcentrifuge
- Micropipettes
- Standard bench top centrifuge with rotor or standard microfuge fit for 0.5 mL reaction tubes, centrifuge for 96-Well reaction plates
- Vortex
- Racks for reaction tubes, also cooled
- 0.5, 1.5 mL and 5 or 15 mL DNase free reaction tubes

## 4.2 Reagents

- TaqMan<sup>®</sup> Universal PCR Master Mix. Applied Biosystems Part No 4318157.

## 4.3 Primers and Probes

Table 5. Primers and probes for the MS11 and *CruA* methods

	Name	DNA Sequence (5' to 3')	Length (nt)
<i>MS11</i>			
Forward primer	SHA086	CAA gAT ggg AAT TAA CAT CTA CAA ATT g	28
Reverse primer	MDB371	gAA ATC CAT gTA AAg CAg CAg gg	23
Probe	TM280	6- FAM-CgA CCA TgT ACA TCC TAC CA -MGB	20
<i>CruA</i>			
Forward primer	MDB510	ggC CAg ggT TTC CgT gAT	18
Reverse primer	MDB511	CCg TCg TTg TAg AAC CAT Tgg	21
Probe	TM458	JOE- AgT CCT TAT gTg CTC CAC TTT CTg gTg CA -BHQ1	29

FAM: 6-carboxyfluorescein; MGB: minor groove binder; JOE: 4,5-dichloro-dimethoxy-fluorescein; BHQ1: black hole quencher.

## 5. References

1. Plant DNA C-values Database. Royal Botanic Gardens, Kew, <http://data.kew.org/cvalues/>

## JRC Mission

As the science and knowledge service of the European Commission, the Joint Research Centre's mission is to support EU policies with independent evidence throughout the whole policy cycle.



**EU Science Hub**

[ec.europa.eu/jrc](https://ec.europa.eu/jrc)



@EU\_ScienceHub



EU Science Hub - Joint Research Centre



Joint Research Centre



EU Science Hub