



Event-specific Method for the Quantification of Oilseed Rape Line RT73 Using Real-time PCR

Validation Report

Biotechnology & GMOs Unit
Institute for Health and Consumer Protection
DG Joint Research Centre

07 February 2007

Executive Summary

The JRC as Community Reference Laboratory for GM Food and Feed (CRL-GMFF) (see Regulation EC 1829/2003), in collaboration with the European Network of GMO Laboratories (ENGL), has carried out an international collaborative study to assess the performance of a quantitative event-specific method to detect and quantify the RT73 transformation event in oilseed rape DNA (unique identifier MON-ØØØ73-7). The collaborative trial was conducted according to internationally accepted guidelines ^(1, 2).

In accordance to Regulation (EC) No 1829/2003 of 22 September 2003 on genetically modified food and feed and to Regulation (EC) No 641/2004 of 6 April 2004 on detailed rules for the implementation of Regulation (EC) No 1829/2003, Monsanto provided the detection method and the samples (whole oilseed rape seed containing the transformation event and whole conventional oilseed rape seeds), whereas the JRC extracted the DNA and prepared the validation samples (calibration samples and blind samples at unknown GM percentage [DNA/DNA]). The collaborative trial involved fourteen laboratories from twelve European Countries.

The results of the international collaborative trial met the ENGL performance requirements and the scientific understanding about satisfactory method performance. Therefore, the CRL-GMFF considers the method validated as fit for the purpose of regulatory compliance.

The results of the collaborative study are made publicly available under <http://gmo-crl.jrc.it/>. The method will also be submitted to ISO 21570 for consideration as an international standard.

Drafted by:
C. Savini



Report Verification Team:
1) S. Langrell



2) W. Moens



Scientific and technical approval:
M. Mazzara



Compliance with CRL Quality System:
S. Cordeil



Authorisation to publish:
G. Van den Eede



Address of contact laboratory:

European Commission, Directorate General-Joint Research Centre
Institute for Health and Consumer Protection (IHCP)
Biotechnology and GMOs Unit – Community Reference Laboratory
Via Fermi 1, 21020 Ispra (VA) - Italy

Report on Steps 1-3 of the Validation Process

Monsanto submitted the detection method and control samples for oilseed rape event RT73 (unique identifier MON-ØØØ73-7) under Article 8 and 20 of Regulation (EC) No 1829/2003 of the European Parliament and of the Council "on genetically modified food and feed".

The Community Reference Laboratory for GM Food and Feed (CRL-GMFF), following reception of the documentation and material, including control samples, (step 1 of the validation process) carried out the scientific assessment of documentation and data (step 2) in accordance with Commission Regulation (EC) No 641/2004 "on detailed rules for the implementation of Regulation (EC) No 1829/2003 of the European Parliament and of the Council as regards the application for the authorisation of new genetically modified food and feed, the notification of existing products and adventitious or technically unavoidable presence of genetically modified material which has benefited from a favourable risk evaluation" and according to its operational procedures ("Description of the CRL-GMFF Validation Process", <http://gmo-crl.jrc.it/guidancedocs.htm>).

The scientific assessment focused on the method performance characteristics assessed against the method acceptance criteria set out by the European Network of GMO Laboratories and listed in the "Definition of Minimum Performance Requirements for Analytical Methods of GMO Testing" (<http://gmo-crl.jrc.it/guidancedocs.htm>) (see Annex 1 for a summary of method acceptance criteria and method performance requirements). During step 2, six scientific assessments were performed and requests of complementary information addressed to the applicant. Upon reception of complementary information, the scientific evaluation of the detection method for event RT73 was positively concluded in September 2006.

In August-October 2006, the CRL-GMFF verified experimentally the method characteristics (step 3, experimental testing of the samples and methods) by quantifying five blind GM-levels within the range 0.1%-8.0% on a copy number basis. The experiments were performed in repeatability conditions and demonstrated that the PCR efficiency, linearity, accuracy and precision of the quantifications were within the limits established by the ENGL.

A Technical Report summarising the results of tests carried out by the CRL-GMFF (step 3) is available on request.

Content

1. INTRODUCTION	5
2. LIST OF PARTICIPATING LABORATORIES	6
3. MATERIALS	7
4. EXPERIMENTAL DESIGN	8
5. METHOD	8
DESCRIPTION OF OPERATIONAL STEPS FOLLOWED	8
6. DEVIATIONS REPORTED	9
7. SUMMARY OF RESULTS	9
PCR EFFICIENCY AND LINEARITY	9
GMO QUANTIFICATION	10
8. METHOD PERFORMANCE REQUIREMENTS	12
9. CONCLUSIONS	13
10. QUALITY ASSURANCE	14
11. REFERENCES	14
12. ANNEX 1: METHOD ACCEPTANCE CRITERIA AND METHOD PERFORMANCE REQUIREMENTS AS SET BY THE EUROPEAN NETWORK OF GMO LABORATORIES (ENGL)	15

1. Introduction

Monsanto submitted the detection method and control samples for oilseed rape event RT73 (unique identifier MON-ØØØ73-7) in accordance to Articles 8 and 20 of Regulation (EC) No 1829/2003 of the European Parliament and of the Council "on genetically modified food and feed".

The Directorate General-Joint Research Centre (JRC, Biotechnology and GMOs Unit of the Institute for Health and Consumer Protection) as Community Reference Laboratory for GM Food and Feed (see Regulation EC 1829/2003) organised the international collaborative ring trial for the event-specific method for the detection and quantification of RT73 oilseed rape. The study involved fourteen laboratories, each members of the European Network of GMO Laboratories (ENGL).

Upon reception of method, samples and related data (step 1), the JRC carried out the assessment of the documentation (step 2) and the in-house evaluation of the method (step 3), according to the requirements of Regulation (EC) 641/2004 and following its operational procedures.

The internal in-house experimental evaluation of the method was carried out in August-October 2006.

Following the evaluation of the data and the results of the in-house laboratory tests, the international ring trial was organised (step 4) and took place in October-November 2006.

A method for DNA extraction from oilseed rape seeds, submitted by the applicant, was evaluated by the CRL-GMFF; laboratory testing of the method was carried out in order to confirm its performance characteristics. The protocol for DNA extraction and a report on method testing is available at <http://gmo-crl.jrc.it/>.

The operational procedure of the collaborative study included the following module:

- ✓ Quantitative real-time PCR (Polymerase Chain Reaction). The methodology is an event-specific real-time quantitative TaqMan[®] PCR procedure for the determination of the relative content of event RT73 DNA to total oilseed rape DNA. The procedure is a simplex system, in which an oilseed rape (OSR) *CruA* (*Cruciferin A*) endogenous assay (reference gene) and the target assay (RT73) are performed in separate wells.

The international collaborative ring-trial was carried out in accordance with the following internationally accepted guidelines:

- ✓ ISO 5725 (1994).
- ✓ The IUPAC "Protocol for the design, conduct and interpretation of method-performance studies" (Horwitz, 1995).

2. List of participating laboratories

As part of the international collaborative ring trial the method was tested in fourteen ENGL laboratories to determine its performance. Clear guidance was given to the laboratories with regards to the standard operational procedures to follow for the common execution of the protocol. The participating laboratories are listed in alphabetical order in Table 1.

Table 1. ENGL laboratories participating in the validation study of oilseed rape RT73.

Laboratory	Country
AGES- Austrian Agency for Health and Food Safety	Austria
Central Science Laboratory (CSL)	UK
Centro Nacional de Alimentación - Agencia Española de Seguridad Alimentaria	Spain
Chemisches und Veterinäruntersuchungsamt Freiburg	Germany
Finnish Customs Laboratory - Tullilaboratorio	Finland
GIP-GEVES	France
Institute of Public Health	Belgium
Istituto Superiore di Sanità (ISS)	Italy
Laboratoire National de la Protection des Végétaux d'Orléans	France
LGC	UK
National Food Administration	Sweden
National Institute for Food Safety and Nutrition	Hungary
National Institute of Biology	Slovenia
Research Institute of Crop Production	Czech Rep.

3. Materials

For the validation of the quantitative event-specific method, genomic DNA was extracted from samples consisting of:

- i) seeds of oilseed rape harbouring the RT73 event in homozygous state (Line "SW RaideRR", Lot Number RDR-0103-11058-S) and
- ii) seeds of conventional oilseed rape (Line "Defender", Lot Number CON-0103-11061-S) in accordance to the provisions of Regulation (EC) No1829/2003, Art 2.11 ["control sample defined as the GMO or its genetic material (positive sample) and the parental organism or its genetic material that has been used for the purpose of the genetic modification (negative sample)].

Samples containing mixtures of 100% RT73 and non-GM oilseed rape genomic DNA at different GMO concentrations were prepared by the CRL-GMFF, using the DNA extracted from control samples provided, in a constant amount of total oilseed rape DNA.

Participants received the following materials:

- ✓ Five calibration samples (160 µl of DNA solution each) for the preparation of the standard curve, denominated from S1 to S5.
- ✓ Twenty unknown DNA samples (80 µl of DNA solution each), denominated from U1-2604 to U20-2604.
- ✓ Amplification reagent control for use on each PCR plate.
- ✓ Reaction reagents, primers and probes for the *CruA* reference gene and for the RT73 specific systems as follows:
 - Universal PCR Master Mix 2X, 3 vials: 5 ml each
 - Sterile distilled water: 8.5 ml
- ✓ Primers and probes (1 tube each) as follows:
 - CruA reference system*
 - BncruA primer 1(10 µM): 160 µl
 - BncruA primer 2(10 µM): 160 µl
 - BncruA probe (10 µM): 160 µl
 - RT73 oilseed rape system*
 - RT73 primer 1 (10 µM): 240 µl
 - RT73 primer 2 (10 µM): 240 µl
 - RT73 probe (5 µM): 160 µl

Table 2 shows the GM contents of the unknown samples (denominated from U1-U20) distributed to ring trial participants.

Table 2. RT73 GM contents

RT73 GM % (GM copy number/OSR genome copy number *100)
0.1
0.4
0.9
4.0
8.0

4. Experimental design

Twenty unknown samples (U1-U20), representing five GM levels, were used in the validation study. On each PCR plate, samples were analysed in parallel with both the RT73 and *CruA* specific system. In total, two plates were run per participating laboratory, with two replicates for each GM level analysed on each run. In total, four replicates for each GM level were analysed. PCR analysis was performed in triplicate for all samples. Participating laboratories carried out the determination of the GM% according to the instructions provided in the protocol and using the electronic tool provided (Excel spreadsheet).

5. Method

Description of operational steps followed

For specific detection of event RT73 genomic DNA, a 108-bp fragment of the recombination region of parts of the construct inserted into the plant genome is amplified using two specific primers. PCR products are measured during each cycle (real-time) by means of a target-specific oligonucleotide probe labelled with two fluorescent dyes: FAM is used as the reporter dye at its 5' end and TAMRA as a quencher dye at its 3' end.

For relative quantification of event RT73 DNA, an OSR-specific reference system amplifies a 101-bp fragment of *CruA* (*Cruciferin A*) oilseed rape endogenous gene (GenBank X14555), using a pair of *CruA* gene-specific primers and a *CruA* gene-specific probe labelled with VIC and TAMRA.

Standard curves are generated both for the RT73 and the *CruA* specific systems by plotting the Ct values measured for the calibration points against the logarithm of the DNA copy numbers and by fitting a regression line into these data. Thereafter, the standard curves are used to estimate the copy numbers in the unknown sample DNA by interpolation from the standard curves.

For relative quantification of event RT73 DNA in a test sample, the RT73 copy number is divided by the copy number of the oilseed rape reference gene (*CruA*) and multiplied by 100 to obtain the percentage value ($GM\% = RT73 / CruA * 100$).

Calibration sample called S1 was prepared by mixing the appropriate amount of RT73 DNA in control non-GM oilseed rape DNA to obtain 10% RT73 in a total of 200 ng oilseed rape DNA. Samples called S2 and S3 were prepared by 1:4 serial dilutions of the S1 sample and samples S4 S5 were prepared by 1:3 serial dilutions of the S3 sample.

The absolute copy numbers in the calibration curve samples are determined by dividing the sample DNA weight (nanograms) by the published average 1C value for oilseed rape genome (1.15)⁽³⁾. The copy number values used in the quantification are provided in Table 3.

Table 3. Copy number values of the standard curve samples.

Sample code	S1	S2	S3	S4	S5
Total amount of DNA in reaction (ng/5 µl)	200	50	12.5	4.2	1.4
Oilseed rape genome copies	173910	43480	10870	3623	1208
RT73 oilseed rape genome copies	17391	4348	1087	362.3	120.8

6. Deviations reported

Twelve laboratories reported no deviations.

One laboratory ran the samples in twenty microliters, as only a 384-well plate machine was available; final concentrations of the supplied primers/probe, buffer and the DNA amount were maintained according to the original protocol.

One laboratory prepared the reaction master-mix in a 15 ml tube and did not centrifuge it.

7. Summary of results

PCR efficiency and linearity

The values of the slopes [from which the PCR efficiency is calculated using the formula $((10^{(-1/\text{slope})} - 1) * 100)$] of the standard curves and of the R^2 (expressing the linearity of the regression) reported by participating laboratories for both PCR systems and runs (reference and GM specific systems, plates A and B) are summarised in Table 4.

Table 4. Values of standard curves' slopes, PCR efficiency and linearity (R^2)

Lab	Plate	RT73			CruA		
		Slope	PCR Efficiency (%)	R^2	Slope	PCR Efficiency (%)	R^2
1	A	-3.37	98.14	0.99	-3.39	97.23	0.99
	B	-3.22	95.36	0.99	-3.42	96.12	0.99
2	A	-3.30	99.12	1.00	-3.92	79.83	0.99
	B	-3.46	94.45	0.98	-3.21	95.21	0.98
3	A	-3.49	93.38	0.98	-3.18	93.49	0.99
	B	-3.35	98.74	0.98	-3.55	91.37	0.99
4	A	-3.47	94.02	0.99	-3.49	93.44	1.00
	B	-3.36	98.63	1.00	-3.45	94.82	0.99
5	A	-3.61	89.08	0.99	-3.44	95.46	0.99
	B	-3.40	96.99	0.99	-3.33	99.59	0.99
6	A	-3.43	95.56	1.00	-3.35	98.98	1.00
	B	-3.39	97.12	1.00	-3.32	99.91	1.00
7	A	-3.77	84.16	0.99	-3.44	95.40	0.99
	B	-3.51	92.63	0.98	-3.49	93.59	0.99
8	A	-3.38	97.58	0.99	-3.14	91.58	0.99
	B	-3.43	95.63	0.99	-3.30	99.07	0.99
9	A	-3.57	90.74	1.00	-3.43	95.72	0.99
	B	-3.45	95.08	1.00	-3.47	94.35	0.99
10	A	-3.53	91.91	1.00	-3.59	89.80	1.00
	B	-3.60	89.52	0.99	-3.51	92.67	0.99
11	A	-3.45	94.87	0.99	-3.30	99.28	1.00
	B	-3.39	97.23	0.99	-3.30	98.97	1.00
12	A	-3.38	97.58	0.99	-3.51	92.58	0.99
	B	-3.44	95.35	1.00	-3.52	92.31	0.99
13	A	-3.16	92.86	0.98	-3.36	98.45	1.00
	B	-3.21	95.07	0.98	-3.26	97.56	1.00
14	A	-3.23	96.13	0.98	-3.29	98.51	0.99
	B	-3.43	95.81	0.99	-3.38	97.52	1.00
Mean		-3.42	94.7	0.99	-3.40	95.1	0.99

The mean PCR efficiency was around 95% for both systems and the linearity of the method was on average 0.99.

Data reported confirm the appropriate performance characteristics of the method tested.

GMO quantification

Table 5 shows the mean values of the four replicates for each GM level as provided by all laboratories. Each mean value is the average of three PCR repetitions.

Table 5. GM% mean values determined by laboratories for all unknown samples (U1-U20)

LAB	Sample GMO content (GM% = GM copy number/OSR genome copy number *100)																			
	0.10				0.40				0.90				4.00				8.00			
	REP 1	REP 2	REP 3	REP 4	REP 1	REP 2	REP 3	REP 4	REP 1	REP 2	REP 3	REP 4	REP 1	REP 2	REP 3	REP 4	REP 1	REP 2	REP 3	REP 4
1	0.06	0.08	0.09	0.04	0.33	0.35	0.32	0.39	0.87	1.03	1.08	0.78	4.24	4.25	4.43	3.04	9.66	9.69	8.15	7.00
2	0.05	0.06	0.04	0.03	0.18	0.29	0.25	0.26	0.64	1.04	0.45	0.55	1.88	5.25	3.42	1.99	6.30	3.27	4.29	4.90
3	0.08	0.13	0.11	0.08	0.48	0.46	0.51	0.53	1.14	1.10	0.84	0.82	4.70	4.07	5.73	5.62	8.36	9.50	11.82	9.07
4	0.07	0.12	0.08	0.07	0.29	0.43	0.38	0.28	0.66	0.83	0.91	0.92	3.44	4.62	4.96	3.95	9.99	9.00	6.84	6.82
5	0.06	0.06	0.11	0.06	0.38	0.69	0.28	0.34	0.72	0.88	1.01	1.00	4.82	5.59	5.86	4.07	7.38	7.22	6.63	10.05
6	0.11	0.06	0.09	0.08	0.39	0.32	0.34	0.41	0.88	0.89	1.01	0.88	5.32	4.62	3.74	5.10	7.89	9.20	8.70	8.47
7	0.05	0.10	0.09	0.06	0.32	0.41	0.54	0.31	0.74	1.00	0.79	1.02	3.69	4.41	4.11	4.18	9.59	7.26	8.02	8.93
8	0.07	0.07	0.07	0.04	0.24	0.37	0.29	0.24	0.67	0.58	0.46	0.83	3.74	3.12	3.46	3.17	6.52	5.58	5.85	5.93
9	0.09	0.10	0.06	0.08	0.28	0.33	0.32	0.28	0.91	1.05	0.69	0.85	3.88	4.12	5.84	3.47	8.91	9.45	7.03	8.22
10	0.10	0.10	0.11	0.08	0.35	0.41	0.38	0.43	0.86	0.78	0.79	1.17	3.77	3.39	4.29	4.42	6.41	7.26	9.18	7.28
11	0.07	0.07	0.07	0.08	0.23	0.45	0.36	0.42	0.90	0.92	0.85	0.84	4.22	4.01	3.90	4.49	8.42	9.05	5.98	8.61
12	0.08	0.10	0.08	0.07	0.32	0.40	0.38	0.24	0.74	0.95	0.80	0.95	4.69	3.77	4.19	4.52	8.73	8.58	8.42	9.90
13	0.05	0.06	0.04	0.07	0.25	0.26	0.25	0.27	0.80	0.64	0.65	0.60	3.26	3.76	3.69	3.49	8.89	9.76	8.60	8.08
14	0.09	0.09	0.07	0.05	0.39	0.40	0.36	0.47	0.94	1.00	0.78	0.87	3.97	4.24	4.35	4.31	9.57	8.36	10.01	10.74


In Figure 1 the relative deviation from the true value for each GM level tested is shown for each laboratory. The coloured bars represent the relative GM quantification obtained by the participating laboratories; green bars represent the overall mean.

As observed in Figure 1, the majority of laboratories underestimated the true value at 0.1%, 0.4% and 0.9%, whereas GM levels of 4% and 8% were overestimated by the majority of participants.

In three laboratories the underestimation was over 25% both at 0.1% and at 0.4%. However, the corresponding values were not indicated as outliers according to Cochran and Grubbs tests (ISO 5725-2).

Overall, the average relative deviation (represented by the green bar) was clearly acceptable at all GM levels tested, indicating a satisfactory accuracy of the method.

Figure 1. Relative deviation (%) from the true value of RT73 for all laboratories



8. Method performance requirements

Among the performance criteria established by ENGL and adopted by the CRL-GMFF (<http://gmo-crl.jrc.it/guidancedocs.htm>, see also Annex 1), repeatability and reproducibility were assessed through an international collaborative trial, carried out with the support of ENGL laboratories (see Table 1). Table 6 illustrates the estimation of repeatability and reproducibility at various GM levels, according to the range of GM percentages tested during the collaborative trial.

The *relative reproducibility standard deviation* (RSD_R), that describes the inter-laboratory variation, should be below 33% at the target concentration and over the majority of the dynamic range, while it should be below 50% at the lower end of the dynamic range.

As can be observed in Table 6, the method fully satisfies this requirement at all GM levels tested. In fact, the highest value of RSD_R (%) is 28% at the 0.1% level, well within the acceptance criterion.

Table 6. Oilseed rape RT73: summary of validation results.

Unknown sample GM%	Expected value (GMO %)				
	0.1	0.4	0.9	4.0	8.0
Laboratories having returned results	14	14	14	14	14
Samples per laboratory	4	4	4	4	4
Number of outliers	0	1	0	1	1
Reason for exclusion	-	1 C	-	1 C	1 DG
Mean value	0.08	0.35	0.85	4.23	8.36
Relative repeatability standard deviation, RSD_r (%)	23	17	17	14	14
Repeatability standard deviation	0.02	0.06	0.14	0.58	1.13
Relative reproducibility standard deviation, RSD_R (%)	28	24	19	17	16
Reproducibility standard deviation	0.02	0.09	0.16	0.71	1.37
Bias (absolute value)	-0.02	-0.05	-0.05	0.23	0.36
Bias (%)	-25	-13	-6.0	5.8	4.5

C = Cochran's test; DG= Grubbs' test; identification and removal of outliers through Cochran and Grubbs tests, according to ISO 5725-2.

Bias is estimated according to ISO 5725 data analysis protocol.

Table 6 further documents the *relative repeatability standard deviation* (RSD_r), as estimated for each GM level. In order to accept methods for collaborative ring trial evaluation, the CRL requires that RSD_r values be below 25%, as indicated by ENGL (Definition of Minimum Performance Requirements for Analytical Methods of GMO Testing" [<http://gmo-crl.jrc.it/guidancedocs.htm>]).

As can be observed from the values reported in Table 6, the method satisfies this requirement throughout the whole dynamic range tested.

The *trueness* of the method is estimated using the measures of the method bias for each GM level. According to ENGL method performance requirements, trueness should be $\pm 25\%$ across the entire dynamic range. In this case the method satisfies this requirement; in fact, the highest value of bias (%) is -25% at the 0.1% level, and the bias (%) for the remaining levels is well within the acceptance criterion.

9. Conclusions

The overall method performance has been evaluated with respect to the method acceptance criteria and method performance requirements recommended by the ENGL (as detailed under <http://gmo-crl.jrc.it/guidancedocs.htm>). The method acceptance criteria were reported by the applicant and used to evaluate the method prior to the international collaborative ring trial (please see Annex 1 for a summary of method acceptance criteria and method performance requirements).

The results obtained during the collaborative study indicate that the analytical module of the method submitted by the applicant complies with ENGL performance criteria. The method is therefore applicable to the control samples provided (see paragraph 3 "Materials"), in accordance with the requirements of Annex I-2.C.2 to Commission Regulation (EC) No 641/2004.

10. Quality assurance

The CRL-GMFF carries out all operations according to ISO 9001:2000 (certificate number: CH-32232) and ISO 17025:2005 (certificate number: DAC-PL-0459-06-00) [DNA extraction, qualitative and quantitative PCR in the area of Biology (DNA extraction and PCR method validation for the detection and identification of GMOs in food and feed materials)]

11. References

1. Horwitz, W. (1995) Protocol for the design, conduct and interpretation of method performance studies, *Pure and Appl. Chem*, 67, 331-343.
2. International Standard (ISO) 5725. 1994. Accuracy (trueness and precision) of measurement methods and results. International Organization for Standardization, Genève, Switzerland.
3. Royal Botanic Garden, Kew. Plant DNA C-values Database (release 4.0, October 2005). MD Bennett and IJ Leitch (<http://www.rbgekew.org.uk/cval/homepage.htm>)

12. Annex 1: method acceptance criteria and method performance requirements as set by the European Network of GMO Laboratories (ENGL)

Method Acceptance Criteria should be fulfilled at the moment of submission of a method (Phase 1: acceptance for the collaborative study).

Method Performance Requirements should be fulfilled in a collaborative study in order to consider the method as fit for its purpose (Phase 2: evaluation of the collaborative study results).

Method Acceptance Criteria

Applicability

Definition: The description of analytes, matrices, and concentrations to which a method can be applied.

Acceptance Criterion: The applicability statement should provide information on the scope of the method and include data for the indices listed below for the product/s for which the application is submitted. The description should also include warnings to known interferences by other analytes, or inapplicability to certain matrices and situations.

Practicability

Definition: The ease of operations, the feasibility and efficiency of implementation, the associated unitary costs (e.g. Euro/sample) of the method.

Acceptance Criterion: The practicability statement should provide indication on the required equipment for the application of the method with regards to the analysis *per se* and the sample preparation. An indication of costs, timing, practical difficulties and any other factor that could be of importance for the operators should be indicated.

Specificity

Definition: Property of a method to respond exclusively to the characteristic or analyte of interest.

Acceptance Criterion: The method should be event-specific and be functional only with the GMO or GM based product for which it was developed. This should be demonstrated by empirical results from testing the method with non-target transgenic events and non-transgenic material. This testing should include closely related events and cases where the limit of the detection is tested.

Dynamic Range

Definition: The range of concentrations over which the method performs in a linear manner with an acceptable level of accuracy and precision.

Acceptance Criterion: The dynamic range of the method should include the 1/10 and at least 5 times the target concentration. Target concentration is intended as the threshold relevant for legislative requirements. The acceptable level of accuracy and precision are described below. The range of the standard curve(s) should allow testing of blind samples throughout the entire dynamic range, including the lower (10%) and upper (500%) end.

Accuracy

Definition: The closeness of agreement between a test result and the accepted reference value.

Acceptance Criterion: The accuracy should be within $\pm 25\%$ of the accepted reference value over the whole dynamic range.

Amplification Efficiency

Definition: The rate of amplification that leads to a theoretical slope of -3.32 with an efficiency of 100% in each cycle. The efficiency of the reaction can be calculated by the following equation: $\text{Efficiency} = [10^{(1/\text{slope})} - 1]$

Acceptance Criterion: The average value of the slope of the standard curve should be in the range of $(-3.1 \geq \text{slope} \geq -3.6)$

R² Coefficient

Definition: The R^2 coefficient is the correlation coefficient of a standard curve obtained by linear regression analysis.

Acceptance Criterion: The average value of R^2 should be ≥ 0.98 .

Repeatability Standard Deviation (RSD_r)

Definition: The standard deviation of test results obtained under repeatability conditions. Repeatability conditions are conditions where test results are obtained with the same method, on identical test items, in the same laboratory, by the same operator, using the same equipment within short intervals of time.

Acceptance Criterion: The relative repeatability standard deviation should be below 25% over the whole dynamic range of the method.

Note: Estimates of repeatability submitted by the applicant should be obtained on a sufficient number of test results, at least 15, as indicated in ISO 5725-3 (1994).

Limit of Quantitation (LOQ)

Definition: The limit of quantitation is the lowest amount or concentration of analyte in a sample that can be reliably quantified with an acceptable level of precision and accuracy.

Acceptance Criterion: LOQ should be less than 1/10th of the value of the target concentration with an $\text{RSD}_r \leq 25\%$. Target concentration should be intended as the threshold relevant for legislative requirements. The acceptable level of accuracy and precision are described below.

Limit of Detection (LOD)

Definition: The limit of detection is the lowest amount or concentration of analyte in a sample, which can be reliably detected, but not necessarily quantified, as demonstrated by single laboratory validation.

Acceptance Criterion: LOD should be less than $1/20^{\text{th}}$ of the target concentration. Experimentally, quantitative methods should detect the presence of the analyte at least 95% of the time at the LOD, ensuring $\leq 5\%$ false negative results. Target concentration should be intended as the threshold relevant for legislative requirements.

Robustness

Definition: The robustness of a method is a measure of its capacity to remain unaffected by small, but deliberate deviations from the experimental conditions described in the procedure.

Acceptance Criterion: The response of an assay with respect to these small variations should not deviate more than $\pm 30\%$. Examples of factors that a robustness test could address are: use of different instrument type, operator, brand of reagents, concentration of reagents, and temperature of reaction.

Method Performance Requirements

Dynamic Range

Definition: In the collaborative trial the dynamic range is the range of concentrations over which the reproducibility and the trueness of the method are evaluated with respect to the requirements specified below.

Acceptance Criterion: The dynamic range of the method should include the $1/10$ and at least five times the target concentration. Target concentration should be intended as the threshold relevant for legislative requirements.

Reproducibility Standard Deviation (RSD_R)

Definition: The standard deviation of test results obtained under reproducibility conditions. Reproducibility conditions are conditions where test results are obtained with the same method, on identical test items, in different laboratories, with different operators, using different equipment. Reproducibility standard deviation describes the inter-laboratory variation.

Acceptance Criterion: The relative reproducibility standard deviation should be below 35% at the target concentration and over the entire dynamic range. An $RSD_R < 50\%$ is acceptable for concentrations below 0.2%.

Trueness

Definition: The closeness of agreement between the average value obtained from a large series of test results and an accepted reference value. The measure of trueness is usually expressed in terms of bias.

Acceptance Criterion: The trueness should be within $\pm 25\%$ of the accepted reference value over the whole dynamic range.