

JRC VALIDATED METHODS, REFERENCE METHODS AND MEASUREMENTS REPORT

Event-specific Method for the Quantification of Maize MON 87411 Using Real-time PCR

Validation report

European Union Reference Laboratory for
Genetically Modified Food and Feed

2016



This publication is a Validated Methods, Reference Methods and Measurements report by the Joint Research Centre (JRC), the European Commission's science and knowledge service. It aims to provide evidence-based scientific support to the European policy-making process. The scientific output expressed does not imply a policy position of the European Commission. Neither the European Commission nor any person acting on behalf of the Commission is responsible for the use which might be made of this publication.

Contact information

Molecular Biology and Genomics Unit

Address: Joint Research Centre, Via Enrico Fermi 2749, TP 201, 21027 Ispra (VA), Italy

E-mail: eurl-gmff@jrc.ec.europa.eu

Tel.: +39 0332 78 5165

Fax: +39 0332 78 9333

JRC Science Hub

<https://ec.europa.eu/jrc>

JRC102146

Ispra, Italy: European Commission, 2016

© European Union, 2016

Reproduction is authorised provided the source is acknowledged.

How to cite: European Union Reference Laboratory for GM Food and Feed; Event-specific Method for the Quantification of Maize MON 87411 Using Real-time PCR; JRC102146

All images © European Union 2016



EUROPEAN COMMISSION
JOINT RESEARCH CENTRE

Institute for Health and Consumer Protection
Molecular Biology and Genomics Unit



Event-specific Method for the Quantification of Maize MON 87411 Using Real-time PCR

Validation Report

23 June 2016

European Union Reference Laboratory for GM Food and Feed

Executive Summary

In line with its mandate^a the European Union Reference Laboratory for GM Food and Feed (EURL GMFF), in collaboration with the European Network of GMO Laboratories (ENGL), validated an event-specific real-time polymerase chain reaction (qPCR) method for detecting and quantifying maize event MON 87411 (unique identifier MON-87411-9). The validation study was conducted according to the EURL GMFF validation procedure [<http://gmo-crl.jrc.ec.europa.eu/guidancedocs.htm>] and the relevant internationally accepted guidelines⁽¹⁻⁵⁾.

In accordance with current EU legislation^b, Monsanto Company provided the detection method and the positive and negative control samples (genomic DNA from homogenised seeds of MON 87411 maize as positive control DNA, and genomic DNA from homogenised seeds of conventional maize as negative control DNA). The EURL GMFF verified the method performance data provided by the applicant, where necessary experimentally, prepared the validation samples (calibration samples and blind samples at different GM percentage [copies GM/total maize genome copies]), organised an international collaborative study and analysed the results.

The EURL GMFF in-house verification and the collaborative study confirmed that the method meets the method performance requirements as established by the EURL GMFF and the ENGL, in line with the provisions of Annex III-3.C.2 to Regulation (EU) No 503/2013, and it fulfils the analytical requirements of Regulation (EU) No 619/2011^c. This validation report is published at <http://gmo-crl.jrc.ec.europa.eu/StatusOfDossiers.aspx>.

^a Regulation (EC) No 1829/2003 of 22 September 2003 "on genetically modified food and feed".

^b Regulation (EC) No 503/2003 of 3 April 2013 "on applications for authorisation of genetically modified food and feed in accordance with Regulation (EC) No 1829/2003 of the European Parliament and of the Council and amending Commission Regulations (EC) No 641/2004 and (EC) No 1981/2006".

^c Regulation (EU) No 619/2011 of 24 June 2011 laying down the methods of sampling and analysis for the official control of feed as regards presence of genetically modified material for which an authorisation procedure is pending or the authorisation of which has expired.

Quality assurance

The EURL GMFF is ISO 17025:2005 accredited [certificate number: ACCREDIA 1172, (Flexible Scope for DNA extraction and qualitative/quantitative PCR)]

The original version of the document containing evidence of internal checks and authorisation for publication is archived within the EURL GMFF quality system.

The EURL GMFF is also ISO 17043:2010 accredited (proficiency test provider) and applies the corresponding procedures and processes for the management of ring trials during the method validation.

The EURL GMFF conducts its activities under the certification ISO 9001:2008 of the Institute for Health and Consumer Protection IHCP provided by SGS.

Address of contact laboratory:

European Commission, Joint Research Centre (JRC)
Institute for Health and Consumer Protection (IHCP)
Molecular Biology and Genomics Unit (MBG)
European Union Reference Laboratory for GM Food and Feed
Via E. Fermi 2749, 21027 Ispra (VA) – Italy
Functional mailbox: eurl-gmff@jrc.ec.europa.eu

Content

1. INTRODUCTION.....	4
2. STEP 1 (DOSSIER ACCEPTANCE) AND STEP 2 (SCIENTIFIC DOSSIER ASSESSMENT AND BIOINFORMATICS ANALYSIS).....	4
3. STEP 3 (EXPERIMENTAL TESTING OF THE SAMPLES AND METHOD).....	7
3.1 DNA EXTRACTION	7
3.2 METHOD PROTOCOL FOR THE PCR ANALYSIS.....	7
3.3 EURL GMFF EXPERIMENTAL TESTING.....	8
3.3.1 <i>Determination of the zygosity ratio in the positive control sample.....</i>	<i>8</i>
3.3.2 <i>In-house verification of the method performance against ENGL method acceptance criteria.....</i>	<i>9</i>
3.4 INTERNATIONAL COLLABORATIVE STUDY (STEP 4).....	10
3.4.1 <i>List of participating laboratories.....</i>	<i>10</i>
3.4.2 <i>Real-time PCR equipment used in the study.....</i>	<i>11</i>
3.4.3 <i>Materials used in the international collaborative study.....</i>	<i>12</i>
3.4.4 <i>Design of the collaborative study.....</i>	<i>13</i>
3.4.5 <i>Deviations reported from the protocol</i>	<i>13</i>
4. RESULTS	14
4.1 EURL GMFF EXPERIMENTAL TESTING.....	14
4.1.1 <i>Zygosity ratio in the positive control sample.....</i>	<i>14</i>
4.1.2 <i>In-house verification of method performance against ENGL method acceptance criteria.....</i>	<i>14</i>
4.2 RESULTS OF THE INTERNATIONAL COLLABORATIVE STUDY.....	16
4.2.1 <i>PCR efficiency and linearity</i>	<i>16</i>
4.2.2 <i>GMO quantification.....</i>	<i>18</i>
4.2.3 <i>Method performance requirements.....</i>	<i>19</i>
5. COMPLIANCE OF THE METHOD FOR DETECTION AND QUANTIFICATION OF EVENT MON 87411 WITH THE REQUIREMENTS OF REGULATION (EU) NO 619/2011.....	21
6. CONCLUSION.....	22
7. REFERENCES.....	22
ANNEX 1: EVENT-SPECIFIC METHOD FOR THE QUANTIFICATION OF MAIZE MON 87411 BY REAL-TIME PCR.....	24

1. Introduction

In line with Regulation (EC) No 1829/2003, Monsanto Company provided the EURL GMFF with an event-specific method for detection and quantification of maize event MON 87411 (unique identifier MON-87411-9) together with genomic DNA as positive and negative control samples.

The dossier was found complete (step 1 of the EURL-GMFF validation procedure) and the scientific dossier assessment (step 2) concluded that the reported method performance characteristics, assessed against the ENGL method acceptance criteria^d allowed moving the method forward to step 3 of the procedure (experimental testing), where the EURL GMFF verified the purity of the control samples provided and conducted an in-house testing of samples and method.

The positive and negative control DNA, submitted in accordance with Art 5(3)(j) and Article 17(3)(j) of Regulation (EC) No 1829/2003, were found of good quality.

Step 3 was completed with the conclusion that the method could be submitted to the collaborative study (step 4). This study confirmed that the method is well suited for quantifying genomic DNA of GM maize MON 87411, appropriately extracted from food or feed, down to a GM content level of 0.1% m/m.

The preparation of the report (step 5) was aligned with the timeline communicated by EFSA for its risk assessment.

2. Step 1 (dossier acceptance) and step 2 (scientific dossier assessment and bioinformatics analysis)

Documentation and data supplied by the applicant were evaluated by the EURL GMFF for completeness (step 1) and compliance with the ENGL method acceptance criteria (step 2).

The specificity of the event-specific assay was verified by the applicant and confirmed by the EURL GMFF by means of bioinformatics analysis, on the basis of the sequence data provided by the applicant.

Specificity assessment by the applicant

The specificity of the event-specific assay was assessed by the applicant in duplicate real-time PCR reactions, according to the method described (Tables 1, 2, 3 and 4), using genomic DNA extracted from 100% of: canola RT73, canola MON 88302; maize MON 87411, maize NK603, maize MON 810, maize MON 863, maize MON 88017, maize MON 89034, maize MON 87460, maize MON

^d EURL/ENGL guidance doc "Definition of Minimum Performance Requirements for Analytical Methods of GMO Testing" (<http://gmo-crl.jrc.ec.europa.eu/guidancedocs.htm>)

87427; cotton MON 531, cotton MON 15985, cotton MON 1445, cotton MON 88913, cotton MON 88701; soybean 40-3-2, soybean MON 89788, soybean MON 87769, soybean MON 87701, soybean MON 87705, soybean MON 87708; alfalfa J101, alfalfa J163, alfalfa KK179 and conventional canola, conventional maize, conventional cotton, conventional soybean, conventional wheat, conventional alfalfa, lentils, sunflower, peanut (shelled), quinoa, and millet.

According to the method developer the MON 87411 method did not react with any sample except the positive control.

The taxon-specific assay is based on the *hmg* and is also found to be fully specific to maize. In addition the HMG assay amplifies a single copy target sequence. Its specificity was assessed by the method developer in duplicate real-time PCR reactions, according to the method described (Tables 1, 2, 3 and 4), using genomic DNA extracted from 100% of: canola RT73, canola MON 88302; maize MON 87411, maize NK603, maize MON 810, maize MON 863, maize MON 88017, maize MON 89034, maize MON 87460, maize MON 87427; cotton MON 531, cotton MON 15985, cotton MON 1445, cotton MON 88913, cotton MON 88701; soybean 40-3-2, soybean MON 89788, soybean MON 87701, soybean MON 87705, soybean MON 87708; alfalfa J101, alfalfa J163, alfalfa KK179 and conventional canola, conventional maize, conventional cotton, conventional soybean, conventional wheat, conventional alfalfa, lentils, sunflower, peanut (shelled), quinoa, and millet.

According to the method developer the maize-specific reference system did not react with any sample except the positive control maize lines

Bioinformatics specificity assessment by the EURL GMFF

The detection method spans the 3' insert-to-plant junction in maize MON 87411. The forward primer "87411 QF" binds to the insert. The reverse primer "87411 QR" binding site was found in the 3' maize (*Zea mays*) genomic border adjacent to the insertion. The probe "87411 QP" binds to the junction between the insert and the 3' genomic region of *Zea mays*.

The amplicon size is expected to be 109 bp, consistent to what reported by the applicant. The sequence of the amplicon was analysed by BLAST (NCBI) against local copies of the "nt" and "patents" databases, and no significant similarity was found with any other published sequence. In addition, the primers were tested against the sequences of the other GMO events present in the Central Core Sequence Information System (CCSIS) of the JRC, as well as the whole genomes of more than 80 plants (including *Brassica rapa*, *Glycine max*, *Oryza sativa*, *Solanum lycopersicum* and *Zea mays*) using the e-PCR prediction tool (NCBI), and no potential amplicon was identified.

Verification of the ENGL acceptance parameter

The parameters of the calibration curves (slope, R^2 coefficient) were determined by the applicant by quantifying three test samples at different GM levels (see Table 1).

Table 1. Summary of the average slope and R² values obtained by the applicant

MON 87411		<i>hmg</i>	
Slope	R²	Slope	R²
-3.45	1.00	-3.41	1.00

According to the ENGL method acceptance criteria, the average value of the slope of the standard curve shall range from -3.1 to -3.6 and the R² coefficient shall be ≥ 0.98 .

Table 1 indicates that the slope and R² coefficient of the standard curves for the GM-system (MON 87411) and the maize-specific high mobility group (*hmg*) system, as established by the applicant, were within the ENGL acceptance criteria.

Also precision and trueness of the method were established by the applicant and 15 values for each of 3 GM levels (expressed as mass fraction of GM-material) were provided. Table 2 reports precision and trueness values for the three GM-levels as provided by the applicant. Both parameters were within the ENGL acceptance criteria (trueness $\pm 25\%$, RSD_r $\leq 25\%$ across the entire dynamic range).

Table 2. Mean %, precision and trueness values provided by the applicant estimated for single measurements

Expected GM %	Test results		
	0.085	1.0	10.0
Measured mean GM %	0.083	1.01	9.91
Precision (RSD _r %)	13.72	8.43	4.42
Trueness (bias %)	-1.85	1.02	-0.88

* Numbers are not rounded but are presented as reported by the applicant

3. Step 3 (experimental testing of the samples and method)

3.1 DNA extraction

Genomic DNA was isolated from ground maize seeds, using a "CTAB-based" protocol coupled with PEG purification previously submitted for detection of maize event MON 88017.

This protocol has already been validated in-house by the EURL GMFF. The protocol for DNA extraction and a report on testing were published in 2008 at http://gmo-crl.jrc.ec.europa.eu/summaries/MON88017_DNAExtr_report.pdf.

In agreement with the ENGL position, endorsing the modularity principle (see also Annex III to Reg. (EU) No 503/2013), and given the similarity in the matrix, the EURL GMFF considers the above mentioned DNA extraction protocol applicable in the context of the validation of the method for maize event MON 87411.

Whenever DNA is extracted from more complex and difficult matrices, a thorough control of the quality of the DNA is recommended in order to ensure that it has the required quality for the subsequent PCR analysis.

3.2 Method protocol for the PCR analysis

The PCR method provided by the applicant is an event-specific, quantitative, real-time TaqMan[®] PCR procedure for the determination of the relative content of GM event MON 87411 DNA to total maize DNA. The procedure is a simplex system, in which a maize specific assay targeting the endogenous gene high mobility group (*hmg*), and the GM target assay for MON 87411 are performed in separate wells. The validated method protocol is published by the EURL GMFF at <http://gmo-crl.jrc.ec.europa.eu/StatusOfDossiers.aspx>.

For the detection of GM event MON 87411, a 109 bp fragment of the region spanning the 3' insert-to-plant junction in maize MON 87411 is amplified using specific primers. PCR products are measured during each cycle (real-time) by means of a target-specific oligonucleotide probe labelled with FAM (6-carboxyfluorescein) as reporter dye at its 5' end and TAMRA (6-carboxytetramethylrhodamine) as quencher dye at its 3' end.

For the relative quantification of GM event MON 87411, a maize taxon-specific system amplifies a 79-bp fragment of a maize high mobility group (*hmg*) endogenous gene, using *hmg* gene-specific primers and a *hmg* gene-specific probe labelled with FAM as reporter dye at its 5' end and TAMRA as quencher dye at its 3' end.

Standard curves are generated for both the MON 87411 and the *hmg* systems by plotting the C_q values measured for the calibration points against the logarithm of the DNA copy numbers and by

fitting a regression line into these data. Thereafter, the standard curves are used to estimate the copy numbers in the test sample DNA by interpolation from the standard curves.

For relative quantification of event MON 87411 DNA in a test sample, the MON 87411 copy number is divided by the copy number of the maize reference gene (*hmg*) and multiplied by 100 to obtain the percentage value ($GM \% = MON\ 87411/hmg \times 100$).

The absolute copy numbers of the calibration curve samples are calculated by dividing the sample DNA mass (nanograms) by the published average 1C value for the maize genome (2.73 pg)⁽⁵⁾. The copy number values used in the quantification, the GMO contents of the calibration samples, and the total DNA quantity used in the PCR reactions are listed in Table 3.

Note: Numerical values presented in the following tables were rounded keeping two digits for values ≤ 1 , one digit for values between 1 and 10 and no digit for values ≥ 10 , unless otherwise stated. The calculations in the MS Excel files however were done over not rounded data. This approach might create small inconsistencies in the numerical values reported in the tables but it allows a higher precision in the final results.

Table 3. Copy number values of the standard curve samples

Sample code	S1	S2	S3	S4	S5
Total amount of maize DNA in the reaction (ng) *	200	50	13	3.1	0.78
Target taxon <i>hmg</i> copies	73260	18315	4579	1145	286
Target MON 87411 copies	7326	1832	458	114	29

3.3 EURL GMFF experimental testing

3.3.1 Determination of the zygosity ratio in the positive control sample

Annex II of Regulation (EU) No 619/2011 requires that "when results are primarily expressed as GM-DNA copy numbers in relation to target taxon-specific DNA copy numbers calculated in terms of haploid genomes, they shall be translated into mass fraction in accordance with the information provided in each validation report of the EURL GMFF." This requires knowledge of the zygosity of the event. In order to satisfy this requirement, the EURL GMFF assessed the zygosity (GM-target to reference target ratio) in the positive control sample submitted by the applicant.

The copy number of the MON 87411 and of the *hmg* targets in the positive control sample were determined by digital PCR (dPCR) performed on the BioMark HD System using the 12.765 digital arrays (Fluidigm).

Reaction mixes were prepared in order to test the zygosity in five replicates to a final volume of 9 μ L and contained 1X TaqMan® Universal PCR Master Mix (Applied Biosystems, Cat. number 4318157), 1X GE sample loading reagent (Fluidigm PN 85000746), primers and probes at concentrations indicated in the corresponding validated method (87411 QF and 87411 QR primers at 450 nM each, QP probe at 200 nM; *hmg* primer 1 and *hmg* primer 2 at 300 nM each, *hmg* probe at 160 nM), and 1 μ L of DNA at a concentration of 3.5 ng/ μ L; the DNA concentration was chosen in order to avoid panel saturation (optimal between 200<positive partitions<700).

Loading of the digital chip was performed according to the manufacturer's instructions by using the IFC controller (Fluidigm). Approximately 4.6 μ L of the reaction mixes were loaded into each well and distributed into the 765 partitions constituting one panel. The experiment was repeated three times for a total number of fifteen data sets for both targets. 'No template controls' were included. Amplification conditions were as reported in Annex 1 or in the Validated Method document at <http://gmo-crl.jrc.ec.europa.eu/StatusOfDossiers.aspx>. Data analysis and copy number calculations were performed using the BioMark digital PCR Analysis software. The range of Cq retention was from 20 to 45. Method precision (RSD, %) at 0.1% GM level in mass fraction was also estimated based on the fifteen replicates.

Calculations of means and variances were carried out according to the procedure outlined for random variables in the Annex 4 of the ENGL guidance document 'Verification of analytical methods for GMO testing when implementing inter-laboratory validated methods'^e.

3.3.2 In-house verification of the method performance against ENGL method acceptance criteria

The method performance characteristics were verified by quantifying on a copy number basis five blind test samples distributed over a range of GM levels (0.06%-5.0%). The experiments were performed on an ABI 7500, an ABI 7900 and a Roche LC480 real-time platform under repeatability conditions and followed the protocol provided by the applicant. Test samples with GM levels 0.55%, 0.9%, 2.0% and 5.0% were tested in two real-time PCR runs with two replicates for each GM-level on each plate (total of four replicates per GM-level). The test sample with GM-level 0.06% was tested in 15 replicates in an additional run for each platform. Average values of the slope and of the R² coefficient of the standard curves and method trueness and precision over the dynamic range were evaluated against the ENGL method acceptance criteria.

^e Verification of analytical methods for GMO testing when implementing inter-laboratory validated methods. European Network of GMO Laboratories (ENGL), 2011.
<http://gmo-crl.jrc.ec.europa.eu/doc/ENGL%20MV%20WG%20Report%20July%202011.pdf>

3.4 International collaborative study (step 4)

The international collaborative trial involved twelve randomly selected laboratories, all being "national reference laboratories, assisting the EURL GMFF for testing and validation of methods for detection", as listed in annex to Regulation (EC) No 1981/2006 who had expressed their interest in participation. The study was carried out in accordance with the following internationally accepted guidelines:

- The IUPAC "Protocol for the design, conduct and interpretation of method-performance studies" (Horwitz, 1995) ⁽⁶⁾
- ISO 5725 "Accuracy (trueness and precision) of measurement methods and results", Part 1 and Part 2 (ISO, 1994); ISO 5725-1:1994/Cor 1 (ISO 1998) and ISO 5725-2:1994/Cor 1 (ISO, 2002) ⁽⁷⁻¹⁰⁾

The objective of the international collaborative study was to verify in experienced laboratories the trueness and precision of the PCR analytical method provided by the applicant and verified in-house by the EURL GMFF.).

3.4.1 List of participating laboratories

The twelve laboratories in participating the MON 87411 international collaborative study were randomly selected from 26 national reference laboratories (NRL) that offered to participate.

Clear guidance was given to the selected laboratories for strictly following the validation protocol that was provided for the execution to them. The participating laboratories are listed in Table 4.

Table 4. Laboratories participating in the validation of the detection method for maize event MON 87411

Laboratory	Country
Scientific Institute of Public Health	BE
State Institute of Chemical and Veterinarian Analysis - Freiburg	DE
National Food Chain Safety Office, Food and Feed Safety Directorate, GMO Laboratory	HU
Crop Research Institute - Reference Laboratory for GMO Detection and DNA fingerprinting	CZ
Science and Advice for Scottish Agriculture	UK
Walloon Agricultural Research Centre - Department Valorization des productions (D4) - Unit 16 - Authentication and traceability	BE
National Food and Veterinary Risk Assessment Institute Molecular Biology and GMO Department	LT
Landeslabor Schleswig-Holstein - Food, Veterinary and Environmental Diagnostic Institute	DE
Plant Breeding and Acclimatization Institute – National Research Institute, GMO Controlling Laboratory	PL
Environment Agency Austria	AT
Center for Agricultural Technology Augustenberg	DE
Service commun des laboratoires du ministere de l'economie et des finances- Etablissement de Strasbourg	FR

3.4.2 Real-time PCR equipment used in the study

Laboratories involved in the collaborative study used the following real-time PCR equipment: five laboratories used ABI 7500, two used ABI 7900, one used ABI 7300, two used Stratagene Mx 3005P, one ABI Step One Plus and one Roche LC480.

This variability of equipment, with its known potential influence on PCR results, reflects the real-life situation in the control laboratories and provides additional assurance that the method is robust and usable under real conditions.

3.4.3 Materials used in the international collaborative study

For the validation of the quantitative event-specific method, test samples were provided by the EURL GMFF to the participating laboratories. They were derived from:

- i) genomic DNA extracted by the applicant from homogenized heterozygous maize seeds harbouring the MON 87411 event, and
- ii) genomic DNA extracted by the applicant from homogenized conventional maize seeds genetically similar to those harbouring the MON 87411 event.

The control samples were prepared by the EURL GMFF from the genomic DNA provided by the applicant in accordance to Regulation (EC) No 1829/2003, Art 2.11^f.

The control samples were used by the EURL GMFF to prepare standards (of known GMO content) and test samples (of undisclosed GM content = blind samples) by mixing MON 87411 maize DNA and non-GM maize DNA.

The calibration sample S1 was prepared by mixing the appropriate amount of MON 87411 DNA with control non-GM maize DNA to obtain a 10% (in copy number ratio) GM-sample. Calibration samples S2-S5 were prepared by 4.0-fold serial dilutions from the S1 sample.

The twelve NRLs participating in the validation study received the following materials:

- ✓ Five calibration samples with known concentrations of GM-event (140 µL of DNA solution each) labelled from S1 to S5 (Table 3).
- ✓ Twenty blinded test DNA samples (70 µL of DNA solution, each at 45 ng/µL) labelled from U1 to U20, representing five GM levels, each in four replicates (Table 5)

Table 5. MON 87411 blinded samples GM % contents

MON 87411 GM %
GM copy number/maize genome copy number x 100
5.0
2.0
0.90
0.55
0.06

^f Control sample defined as the GMO or its genetic material (positive sample) and the parental organism or its genetic material that has been used for the purpose of the genetic modification (negative sample).

- ✓ Reaction reagents:
 - TaqMan[®] Universal PCR Master Mix (2x), one vial: 12 mL
 - distilled sterile water, one vial: 7.1 mL
- ✓ Primers and probes (1 tube each) as follows:
 - hmg* taxon-specific assay
 - *hmg* primer 1 (10 µM): 240 µL
 - *hmg* primer 2 (10 µM): 240 µL
 - *hmg* probe (10 µM): 130 µL
 - MON 87411 assay
 - 87411 QF (10 µM): 720 µL
 - 87411 QR (10 µM): 720 µL
 - 87411 QP (10 µM): 300 µL

3.4.4 Design of the collaborative study

Participating laboratories received a detailed validation protocol that included the exact design of the PCR plates, ensuring that on each PCR plate the samples were analysed for the MON 87411 event-specific system and for the *hmg* taxon-specific system. In total, two plates were run by each participating laboratory.

The laboratories prepared the PCR master-mixes for the maize event MON 87411 and the *hmg* assay in accordance with the description provided in the validation protocol. Calibration and test samples were loaded on the PCR plates as per pre-determined plate layout.

The amplification reaction followed the cycling program specified in the protocol. Participants determined the GM % in the test samples according to the instructions and also reported the raw data to the EURL GMFF on an Excel sheet that was designed, validated and distributed by the EURL GMFF. All data are stored by the EURL GMFF on a dedicated and protected server.

The EURL GMFF analysed the data against the parameters and the limits set by the ENGL, i.e. trueness, precision, amplification efficiency and linearity.

3.4.5 Deviations reported from the protocol

Nine laboratories reported no deviations from the validation protocol.

One laboratory reported that reagents were not kept in ice after being thawed and that the master mix was mixed but not spun down.

Another laboratory reported pipetting errors in three wells.

One laboratory found at thawing of samples further to reception that the tube of the universal master mix had not properly closed.

One laboratory informed that the configuration of the real-time PCR machine (Agilent/Stratagene Mx3005P) used did not contain any filter corresponding to the quencher dye. Therefore, only filters for the reporter dye and the passive reference dye were used to perform both real-time PCR plates. However, this is not a deviation from the validation protocol but it is reported for information.

4. Results

4.1 EURL GMFF experimental testing

4.1.1 Zygosity ratio in the positive control sample

The results of the digital PCR analysis conducted by the EURL GMFF on the MON 87411 and *hmg* targets to determine the zygosity ratio in the positive control samples are shown in Table 6.

Table 6. Zygosity ratio of the MON 87411 and *hmg* targets in the positive control sample.

Mean ratio (MON 87411/ <i>hmg</i>)	0.62
Standard deviation	0.04
RSD _r (%)	6.3
Standard error of the mean	0.01
Upper 95% CI of the mean	0.64
Lower 95% CI of the mean	0.60

The mean ratio (MON 87411/*hmg*) equals 0.62%. The 95% confidence interval (CI) spans around 0.62% and therefore the mean ratio is not significantly different from the expected ratio for a maize heterozygous control sample whose GM parental contribution is of female origin and assuming a single-copy endogenous gene target, for an $\alpha = 0.05$.

Hence:

$$0.06 \text{ GM \% in DNA copy number ratio} = 0.1 \text{ GM \% in mass fraction}$$

4.1.2 In-house verification of method performance against ENGL method acceptance criteria

Test samples with GM levels from 0.55% to 5.0% were tested in three real-time PCR runs with two replicates for each GM-level on each plate (total of four replicates per GM-level). The sample at 0.06% GM-level was tested for its precision in quantification at 15 replicates in separate runs.

Tests were conducted on ABI 7500, on ABI 7900 and on Roche LC480 for robustness.

The standard curve parameters and the results of efficiency, linearity, trueness and precision obtained in the two real-time PCR runs with the test samples are shown in Tables 7A, 7B, 8, 9 and 10.

According to the ENGL method acceptance criteria, the average value of the slope of the standard curve shall range from -3.1 to -3.6 and the R^2 coefficient shall be ≥ 0.98 . Table 7A and 7B document that the slopes of the standard curves and the R^2 coefficients were within the limits established by the ENGL. The EURL GMFF in-house results confirm the data provided by the applicant.

Table 7A. Standard curve parameters of the real-time PCR tests, carried out on ABI 7500, ABI 7900, and Roche LC480 to quantify GM-levels in the range 0.55% to 5.0% in four replicates each. Slope and R^2 coefficient values were rounded to two digits.

	MON 87411 system			hmg system		
	Slope	PCR efficiency*	R^2	Slope	PCR efficiency*	R^2
Run A	-3.51	93	1.00	-3.40	97	1.00
Run B	-3.50	93	1.00	-3.40	97	1.00
Run C	-3.28	102	1.00	-3.32	100	1.00
Run D	-3.51	93	0.99	-3.41	97	1.00
Run E	-3.34	99	0.99	-3.30	101	1.00
Run F	-3.11	109	1.00	-3.30	101	1.00

* PCR efficiency (%) is calculated using the formula $\text{Efficiency} = (10^{(-1/\text{slope})} - 1) \times 100$

Runs A-B were carried out on ABI 7500; Runs C-D were carried out on ABI 7900; Runs E and F were carried out on Roche LC 480.

Table 7B. Standard curve parameters of the real-time PCR tests, carried out on ABI 7500, ABI 7900, and Roche LC480 to quantify the GM-level 0.06% in 15 replicates. Slope and R^2 coefficient values were rounded to two digits.

	MON 87411 system			hmg system		
	Slope	PCR efficiency*	R^2	Slope	PCR efficiency*	R^2
Run G	-3.50	93	1.00	-3.37	98	1.00
Run H	-3.30	101	1.00	-3.36	98	1.00
Run I	-3.24	103	1.00	-3.28	102	1.00

* PCR efficiency (%) is calculated using the formula $\text{Efficiency} = (10^{(-1/\text{slope})} - 1) \times 100$

Run G was carried out on ABI 7500; Run H was carried out on ABI 7900; Run I was carried out on Roche LC 480.

According to the ENGL method acceptance criteria the method trueness (measured as bias in % of the target GM level) should be within $\pm 25\%$ of the accepted reference value over the entire dynamic range and the precision, expressed as RSD_r % (relative standard deviation of repeatability), should be $\leq 25\%$, also over the entire dynamic range.

Tables 8, 9 and 10 show that trueness and precision of quantification were within the limits established by the ENGL for both PCR machines used.

Table 8. Values of trueness and precision as established by the EURL GMFF in its in-house verification using an ABI 7500. GM% in copy/copy.

Target GM-levels %	Measured GM-level %	Bias % of the target GM-level	Precision (RSD_r %)
5.0	4.4	-13	4.6
2.0	1.8	-12	2.7
0.90	0.85	-6.0	2.8
0.55	0.49	-11	7.6
0.06	0.06	2.5	17

Table 9. Values of trueness and precision as established by the EURL GMFF in its in-house verification using an ABI 7900. GM% in copy/copy.

Target GM-levels %	Measured GM-level %	Bias % of the target GM-level	Precision (RSD_r %)
5.0	5.3	4.9	11
2.0	1.9	-4.2	11
0.90	0.94	4.9	11
0.55	0.53	-2.8	12
0.06	0.05	-23	23

Table 10. Values of trueness and precision as established by the EURL GMFF in its in-house verification using a Roche LC 480. GM% in copy/copy.

Target GM-levels %	Measured GM-level %	Bias % of the target GM-level	Precision (RSD_r %)
5.0	4.9	-1.7	5.0
2.0	2.0	-1.6	5.1
0.9	0.91	0.63	8.0
0.55	0.51	-7.3	10
0.06	0.05	-19	15

4.2 Results of the international collaborative study

4.2.1 PCR efficiency and linearity

The PCR efficiency (%) and R² values (expressing the linearity of the regression) for the standard curve, reported by participating laboratories are displayed in Table 11. The PCR efficiency (%) was calculated from the standard curve slopes using the formula:

$$\text{Efficiency} = (10 * (-1/\text{slope})) - 1) \times 100$$

Table 11 indicates that the efficiency of amplification for the MON 87411 system ranges from 83% to 100% and the linearity from 0.99 to 1.00; the amplification efficiency for the maize-specific system ranges from 86% to 99% and the linearity from 0.99 to 1.00 as well. The mean PCR efficiency was 91% for MON 87411 assay and 92% for the *hmg* one. The average R^2 of the methods was 0.99 and 1.00 for the MON 87411 and *hmg* assays, respectively. Both PCR efficiency and linearity values were within the ENGL acceptance criteria.

Table 11. Values of slope, PCR efficiency and R^2 obtained during the international collaborative trial. Slope and R^2 coefficient values were rounded to two digits.

Lab	Plate	MON 87411			<i>hmg</i>		
		Slope	PCR Efficiency (%)	R^2	Slope	PCR Efficiency (%)	R^2
1	A	-3.72	86	0.99	-3.60	90	1.00
	B	-3.47	94	1.00	-3.55	91	1.00
2	A	-3.80	83	0.99	-3.58	90	1.00
	B	-3.82	83	0.98	-3.51	93	1.00
3	A	-3.45	95	0.99	-3.51	93	1.00
	B	-3.60	90	0.99	-3.49	93	0.99
4	A	-3.52	92	1.00	-3.56	91	1.00
	B	-3.31	100	1.00	-3.34	99	0.99
5	A	-3.49	93	0.99	-3.59	90	1.00
	B	-3.55	91	0.99	-3.54	91	1.00
6	A	-3.45	95	1.00	-3.44	95	0.99
	B	-3.41	96	0.99	-3.41	96	1.00
7	A	-3.41	96	1.00	-3.52	92	1.00
	B	-3.37	98	1.00	-3.47	94	1.00
8	A	-3.80	83	0.99	-3.50	93	1.00
	B	-3.65	88	0.98	-3.34	99	0.98
9	A	-3.82	83	0.99	-3.68	87	0.99
	B	-3.66	88	0.99	-3.71	86	1.00
10	A	-3.67	87	1.00	-3.63	89	1.00
	B	-3.64	88	1.00	-3.50	93	1.00
11	A	-3.57	91	1.00	-3.56	91	0.99
	B	-3.52	92	1.00	-3.50	93	1.00
12	A	-3.63	89	0.99	-3.56	91	1.00
	B	-3.57	91	1.00	-3.57	91	1.00
Mean		-3.58	91	0.99	-3.53	92	1.00

These results confirm the appropriate performance characteristics of the methods tested in terms of efficiency and linearity.

4.2.2 GMO quantification

Table 12 reports the values of quantification for the four replicates of each GM level as reported by each of the twelve participating laboratories.

Table 12. GM% values determined by laboratories for test samples


LAB	GMO content (%) *															
	0.06				0.55				0.90				2.0			
	REP 1	REP 2	REP 3	REP 4	REP 1	REP 2	REP 3	REP 4	REP 1	REP 2	REP 3	REP 4	REP 1	REP 2	REP 3	REP 4
1	0.06	0.07	0.10	0.08	0.60	0.61	0.58	0.66	0.98	0.95	1.10	1.19	2.06	2.14	2.06	2.07
2	0.08	0.09	0.09	0.09	0.80	0.73	0.73	0.76	1.11	1.08	1.20	1.14	2.27	2.28	2.26	2.19
3	0.06	0.09	0.07	0.06	0.68	0.68	0.66	0.65	1.07	1.12	1.04	1.04	2.39	2.16	2.27	2.07
4	0.04	2.63	0.07	0.08	0.76	0.54	0.54	0.70	0.83	0.82	0.93	1.02	1.88	1.86	1.55	2.09
5	0.08	0.06	0.08	0.07	0.70	0.60	0.62	0.77	1.01	0.99	1.13	0.93	2.32	2.20	2.31	2.31
6	0.05	0.06	0.07	0.07	0.60	0.54	0.65	0.60	0.93	0.93	0.79	0.94	1.73	1.94	1.87	2.16
7	0.08	0.06	0.06	0.07	0.60	0.56	0.62	0.55	0.93	0.97	0.94	1.01	1.96	1.99	2.12	2.26
8	0.09	0.08	0.04	0.05	0.68	0.56	0.47	0.41	1.48	0.33	0.89	0.95	1.06	0.97	1.69	2.17
9	0.07	0.07	0.10	0.09	0.67	0.59	0.69	0.67	1.07	1.02	1.18	1.08	2.28	2.45	2.41	2.24
10	0.05	0.06	0.05	0.06	0.60	0.61	0.48	0.49	0.77	0.92	0.97	1.00	1.97	2.00	1.81	1.42
11	0.07	0.08	0.07	0.07	0.65	0.74	0.74	0.57	1.27	1.01	1.03	1.20	2.33	2.35	2.23	2.14
12	0.07	0.07	0.08	0.06	0.64	0.68	0.62	0.59	0.95	0.95	0.98	1.03	2.23	2.01	2.10	2.17

* GMO% = (GMO copy number/maize genome copy number) x 100

n.a. not available

A graphical representation of the data reported in Table 12 is provided in Figure 1 that shows the relative deviation from the true value for each GM level tested for the participating laboratory. The coloured bars represent the deviation of the GM level measured in % of the true GM level; the green bar on the right represents the mean relative deviation over all data retained after eliminating outliers for each GM level.

Figure 1. Relative deviation (%) from the true value of GM level *



* Regarding PL04 and PL06 at level 0.9%, as well as PL09 and PL12 at level 5.0%, very small relative deviations from the true value were observed and therefore the corresponding histograms do not show up in Figure 1. PL = participating laboratory.

Overall, most laboratories' mean relative deviations from the true values were within a maximum of 25%. At GM-level 0.06% eight laboratories were within the limit, at GM-level 0.55% eleven and at GM-level 0.9% ten laboratories. At 2.0% and 5.0% eleven GM-level laboratories were within the limit. Three laboratories overestimated GM-level 0.06% by more than 25%, similarly, one laboratory overestimated GM-level 0.55% and 0.9% by more than 25%, with a trend for overestimation for all laboratories at all GM levels.

The mean bias generated by all laboratories ranged between +4.0% and +18.8%. The method is well within the accepted limits established by the ENGL ($\pm 25\%$ over the entire dynamic range).

4.2.3 Method performance requirements

Among the performance requirements established by ENGL and adopted by the EURL GMFF (<http://gmo-crl.jrc.ec.europa.eu/guidancedocs.htm>), repeatability and reproducibility are to be assessed through an international collaborative trial. Table 13 illustrates the estimation of repeatability and reproducibility at the various GM levels tested during the study (see Table 4 for a list of the participant laboratories).

According to the ENGL method performance requirements the relative reproducibility standard deviation (RSD_R), that describes the inter-laboratory variation, should be below 35% at the target concentration and over the majority of the dynamic range, while it should be below 50% at the lower end of the dynamic range.

As it can be observed in Table 13, the method satisfies this requirement at all GM levels tested. Indeed, the highest value of RSD_R (%) is 19% at the 0.06% GM level, thus within the acceptance criterion.

Table 13. Summary of validation results for the MON 87411 method, expressed as GM copy numbers in relation to target taxon copy numbers. Standard deviation values and absolute bias value are rounded to three digits.

	Test Sample Expected GMO %				
	0.06	0.55	0.9	2.0	5.0
Laboratories having returned valid results	12	12	12	12	12
Samples per laboratory	4	4	4	4	4
Number of outliers	1	0	1	1	1
Reason for exclusion	C*	-	C*	C*	C*
Mean value	0.07	0.63	1.0	2.1	5.2
Relative repeatability standard deviation, RSD_r (%)	16	10	7.9	6.8	9.1
Repeatability standard deviation	0.012	0.066	0.080	0.143	0.471
Relative reproducibility standard deviation, RSD_R (%)	19	13	11	11	11
Reproducibility standard deviation	0.013	0.084	0.112	0.222	0.596
Bias** (absolute value)	0.011	0.080	0.113	0.112	0.202
Bias (%)	19	14	13	5.6	4.0

* C= Cochran's test; identification and removal of outliers through Cochran and Grubbs tests, according to ISO 5725-2.

** Bias is estimated according to ISO 5725 data analysis protocol.

Table 13 also documents the relative repeatability standard deviation (RSD_r) estimated for each GM level. In order to accept methods for collaborative study, the EURL GMFF and ENGL require that the RSD_r value indicated by the applicant and confirmed by the EURL GMFF through in-house experiments, is below 25% (see ENGL document "Definition of Minimum Performance Requirements for Analytical Methods of GMO Testing" <http://gmo-crl.jrc.ec.europa.eu/guidancedocs.htm>). As it can be observed from the values reported, the repeatability standard deviation is below 25% at all GM levels, with the highest value of 16% at the 0.06% GM level.

The trueness of the method is estimated using the measures of the method bias for each GM level. According to ENGL method performance requirements, trueness should be $\pm 25\%$ across the entire dynamic range. The method satisfies this requirement across the dynamic range tested, with the highest value of bias (%) of 19% at the 0.06% GM level.

5. Compliance of the method for detection and quantification of event MON 87411 with the requirements of Regulation (EU) No 619/2011

To verify the compliance of the method under validation with the requirements of Regulation (EU) No 619/2011, the following steps were carried out and their outcome is summarised in Table 14:

- at step 2 of the validation process (scientific assessment of the dossier), the EURL GMFF acknowledged that the RSD_r value at the 0.085% level shown by the applicant's dossier (expressed as mass fraction of GM-material) was 9.71%, based on 15 replicates (Table 2), hence below the maximum value of 25% required by the ENGL. The EURL GMFF therefore concluded that it could accept the applicant's data on method performance;
- at step 3 of the validation process (experimental testing of samples and methods), the EURL GMFF determined the RSD_r % value at the level of 0.1% in mass fraction of GM-material (corresponding to 0.06% expressed in terms of copy number ratio). The experiments were carried out under repeatability conditions on fifteen replicates. The RSD_r resulted to range between 15% and 23% (Table 8, 9 and 10) depending on the qPCR platform applied, hence also below 25%;
- the collaborative study (step 4 of the validation process) established that over the twelve participating laboratories at the level of 0.1% related to mass fraction of GM-material the RSD_r of the method was 16 %, therefore also below 25% and well in line with the previous data.

The outcome of the different steps is summarised Table 14.

Table 14. Precision of the event-specific method for quantitative detection of MON 87411 at or around 0.1% level related to mass fractions of GM material.

Source	RSD_r %	GM %
Applicant's method optimisation	13.72%	0.085%
EURL GMFF tests	15-23%	0.1%
Collaborative study	16%	0.1%

Based on the results of the EURL GMFF in-house verification and of the international collaborative study, it is concluded that the method RSD_r % is lower than 25% at the level of 0.1% related to mass fraction of GM material, hence the method meets the requirement laid down in Regulation (EU) No 619/2011.

6. Conclusion

The method provided by the applicant has been validated in accordance to the EURL GMFF validation process, respecting all requirements of the relevant EU legislation and international standards for method validation.

This validation study confirmed that the method is applicable to the control samples provided by the applicant (see paragraph 3.4.3), in accordance with the requirements of Annex I-3.C.2 to Commission Regulation (EU) No 503/2013 and (EU) No 619/2011 and meets all method performance requirements established by the ENGL and the EURL GMFF. The method is therefore valid to be used for regulatory purposes, including the quantification of low level presence of 0.1% (m/m) of the GM event. It can be assumed that it is applicable to any appropriately extracted maize genomic DNA.

In any case the user of the method is advised to verify the quality of the extracted genomic DNA in order to ensure that it is suitable for the subsequent PCR analysis. This is particularly relevant for more complex matrices of samples from food and feed products.

The validated method is described in detail as "Validated Method" at <http://gmo-crl.jrc.ec.europa.eu/StatusOfDossiers.aspx> and in Annex 1.

7. References

1. Bauer T et al., 2003. The effect of processing parameters on DNA degradation in food. *Eur. Food Res. Technol.*, 217: 338-343.
2. Murray SR et al., 2007. Use of quantitative real-time PCR to estimate maize endogenous DNA degradation after cooking and extrusion or in food products. *J. Agric. Food Chem.*, 55: 2231-2239.
3. Demeke T and Jenkins GR, 2010. Influence of DNA extraction methods, PCR inhibitors and quantification methods on real-time PCR assay of biotechnology-derived traits. *Anal. and Bioanal. Chem.*, 396: 1977-1990.
4. Peano C et al., 2004. Qualitative and quantitative evaluation of the genomic DNA extracted from GMO and non-GMO foodstuffs with four different extraction methods. *J. Agric. Food Chem.*, 52: 6962, 6968.
5. Plant DNA C-values Database, <http://data.kew.org/cvalues/>
6. Horwitz W. Protocol for the design, conduct and interpretation of method- performance studies, *Pure and Appl. Chem.* 1995; 67: 331-343.

7. International Standard (ISO) 5725-1, 1994. Accuracy (trueness and precision) of measurement methods and results. Part 1: General principles and definitions. International Organization for Standardization, Genève, Switzerland.
8. ISO 5725-1:1994/Cor 1:1998.
9. International Standard (ISO) 5725-2, 1994. Accuracy (trueness and precision) of measurement methods and results. Part 2: Basic method for the determination of repeatability and reproducibility of a standard measurement method. International Organization for Standardization, Genève, Switzerland.
10. ISO 5725-2:1994/Cor 1:2002.

Annex 1: Event-specific Method for the Quantification of maize MON 87411 by Real-time PCR

Validated Method

Method development:

Monsanto Company

1. General information and summary of the methodology

This protocol describes an event-specific real-time quantitative TaqMan[®] PCR (polymerase chain reaction) procedure for the determination of the relative content of maize event MON 87411 DNA to total maize DNA in a sample.

Template DNA extracted by means of suitable methods should be tested for quality and quantity prior to use in PCR assays. Tests for the presence of PCR inhibitors (e.g. monitor run of diluted series, use of DNA spikes) are also recommended to ensure suitability of the extracted DNA.

For the detection of GM event MON 87411, a 109 bp fragment of the region spanning the 3' insert-to-plant junction in maize MON 87411 is amplified using specific primers. PCR products are measured during each cycle (real-time) by means of a target-specific oligonucleotide probe labelled with FAM (6-carboxyfluorescein) as reporter dye at its 5' end and TAMRA (6-carboxytetramethylrhodamine) as quencher dye at its 3' end.

For the relative quantification of GM event MON 87411, a maize taxon-specific system amplifies a 79-bp fragment of a maize high mobility group (*hmg*) endogenous gene (Accession number, GeneBank: AJ131373.1), using *hmg* gene-specific primers and a *hmg* gene-specific probe labelled with FAM as reporter dye at its 5' end and TAMRA as quencher dye at its 3' end.

The measured fluorescence signal passes a threshold value after a certain number of cycles. This threshold cycle is called the "Cq" value. For quantification of the amount of MON 87411 DNA in a test sample, Cq values for the MON 87411 and the *hmg* systems are determined for the sample. Standard curves are then used to estimate the relative amount of MON 87411 DNA to total maize DNA.

2. Validation and performance characteristics

2.1 General

The method was optimised for suitable DNA extracted from mixtures of genetically modified and conventional maize seeds. Precision and trueness of the method were tested through an international collaborative ring trial using DNA samples at different GM contents.

2.2 Collaborative trial

The method was validated in an international collaborative study by the European Union Reference Laboratory for GM Food and Feed (EURL GMFF). The study was undertaken with twelve participating laboratories in June 2015-July 2015.

A detailed validation report can be found at <http://gmo-crl.jrc.ec.europa.eu/StatusOfDossiers.aspx>.

2.3 Limit of detection (LOD)

According to the method developer, the relative LOD of the method is at least 0.040% (related to mass fraction of GM material) in 200 ng of total suitable maize DNA. The relative LOD was not assessed in the collaborative study.

2.4 Limit of quantification (LOQ)

According to the method developer, the relative LOQ of the method is at least 0.085% (related to mass fraction of GM material) in 200 ng of total suitable maize DNA. The lowest relative GM content of the target sequence included in the collaborative trial was 0.1% (mass fraction of GM-material).

2.5 Molecular specificity

The method exploits a unique DNA sequence in the region spanning the 3' insert-to-plant junction in maize MON 87411 and is therefore event-specific for the event MON 87411. This was confirmed in the validation study.

3. Procedure

3.1 General instructions and precautions

- The procedures require experience of working under sterile conditions.
- Laboratory organisation, e.g. "forward flow direction" during PCR-setup, should follow international guidelines, e.g. ISO 24276:2006.
- PCR reagents should be stored and handled in a separate room where no nucleic acids (with exception of PCR primers or probes) or DNA degrading or modifying enzymes have been handled previously. All handling of PCR reagents and controls requires dedicated equipment, especially pipettes.
- All the equipment should be sterilised prior to use and any residue of DNA should have been removed. All material used (e.g. vials, containers, pipette tips, etc.) must be suitable for PCR and molecular biology applications. They must be DNase-free, DNA-free, sterile and unable to adsorb protein or DNA.
- Filter pipette tips protected against aerosol should be used.
- Powder-free gloves should be used and changed regularly
- Laboratory benches and equipment should be cleaned periodically, with 10% sodium hypochlorite solution (bleach).
- Pipettes should be checked regularly for precision and calibrated, if necessary.

- All handling steps, unless specified otherwise, should be carried out at room temperature.
- In order to avoid repeated freeze/thaw cycles aliquots should be prepared.

3.2 Real-time PCR for quantitative analysis of maize event MON 87411

3.2.1 General

The real-time PCR set-up for the taxon (*hmg*) and the GMO (event MON 87411) target sequences are carried out in separate vials. Multiplex qPCR (using differential fluorescent labels for the probes) has not been tested or validated by the EURL GMFF.

The method is developed for a total volume of 50 µL per reaction mixture for the GM (event MON 87411) and of 25 µL per reaction mixture for the taxon (*hmg*) with the reagents as listed in Table 2 and Table 3.

3.2.2 Calibration

The calibration curves have to be established on at least five samples. The first point of the calibration curve (S1) should be established for a sample containing 10% maize MON 87411 DNA in a total of 200 ng of maize DNA (corresponding to 73260 maize genome copies with one haploid genome assumed to correspond to 2.73 pg of maize genomic DNA) ⁽¹⁾. Standards S2 to S5 are to be prepared by serial dilutions (dilution factor 4.0) according to Table 1 below.

Table 1. Copy number values of the standard curve samples

Sample code	S1	S2	S3	S4	S5
Total amount of maize DNA in reaction (ng) *	200	50	12.5	3.1	0.78
<i>hmg</i> copies	73260	18315	4579	1145	286
MON 87411 copies	7326	1832	458	114	29

* Total nanograms are rounded to the integral value

A calibration curve is to be produced by plotting the C_q values against the logarithm of the target copy number for the calibration points. This can be done by means of spreadsheet software, e.g. Microsoft Excel, or directly by options available with the software.

The copy number measured for each unknown sample DNA is obtained by interpolation from the standard curves.

3.2.3 Real-time PCR set-up

1. Thaw, mix and centrifuge the components needed for the run. Keep thawed reagents on ice.
2. In two tubes on ice, add the components in the order mentioned below (except DNA) to prepare the reaction mixes for the MON 87411 maize specific system (Table 2) and the *hmg* reference gene system (Table 3). Please note that additional volume is included in the total to cover pipetting variability due to the viscosity of the solution.

Table 2. Amplification reaction mixture in the final volume/concentration per reaction well for the MON 87411 assay.

Component	Final concentration	µL/reaction
TaqMan [®] Universal PCR Master Mix (2x)	1x	25
87411 QF (10 µM)	450 nM	2.25
87411 QR (10 µM)	450 nM	2.25
87411 QP* (10 µM)	200 nM	1.0
Nuclease free water	-	15.5
DNA (max 200 ng)	-	4.0
Total reaction volume:		50 µL

*TaqMan[®] probe is labeled with 6-FAM at its 5'-end and TAMRA at its 3'-end

Table 3. Amplification reaction mixture in the final volume/concentration per reaction well for the maize *hmg* assay.

Component	Final concentration	µL/reaction
TaqMan [®] Universal PCR Master Mix (2x)	1x	12.5
<i>hmg</i> primer 1 (10 µM)	300 nM	0.75
<i>hmg</i> primer 2 (10 µM)	300 nM	0.75
<i>hmg</i> probe* (10 µM)	160 nM	0.40
Nuclease free water	-	6.6
DNA (max 200 ng)	-	4.0
Total reaction volume:		25 µL

*TaqMan[®] probe is labeled with 6-FAM at its 5'-end and TAMRA at its 3'-end

3. Mix well and centrifuge briefly.
4. Prepare two 0.5 mL reaction tubes (one for the maize MON 87411 and one for the *hmg* system) for each DNA sample to be tested (standard curve samples, unknown samples and control samples).
5. Add into each reaction tube the amount of reaction mix for 3.5 PCR repetitions (73.5 µL for the *hmg* system and 161 µL for the MON 87411 maize system). Add to each tube the correct

amount of DNA for 3.5 PCR repetitions (14 µL DNA). The volume for the additional 0.5 repetition will ensure adequate volume when loading the samples. Vortex each tube for approx. 10 seconds. This step is mandatory to reduce to a minimum the variability among the repetitions of each sample.

6. Spin down the tubes. Aliquot 50 µL for MON 87411 system and 25 µL for the *hmg* reference system in each well.
7. Place an optical cover on the reaction plate and briefly centrifuge the plate.
8. Place the reaction plate in the real-time PCR apparatus (possibly apply a compression pad, depending on the model), according to the manufacturer's instructions and your Standard Operating Procedures and start the run.
9. Select FAM as reporter dye for the MON 87411 and FAM for the *hmg* reference system. Define TAMRA as quencher dye for MON 87411 specific system and TAMRA for *hmg* reference system. Select ROX as the passive reference dye. Enter the correct reaction volume (50 µL).
10. Run the PCR with the cycling program described in Table 4.

Table 4. Cycling program for MON 87411/*hmg* assays.

Step	Stage	T (°C)	Time (s)	Acquisition	Cycles	
1	UNG*	50	120	No	1	
2	Initial denaturation	95	600	No	1	
3	Amplification	Denaturation	95	15	No	40
		Annealing & Extension	60	60	Yes	

*UNG: Uracil-N-glycosylase

3.3 Data analysis

After the real-time PCR, analyse the run following the procedure below:

- a) Set the threshold: display the amplification curves of one assay (e.g. MON 87411) in logarithmic mode. Locate the threshold line in the area where the amplification profiles appear linear (exponential phase of PCR) and where there is no "fork effect" between repetitions of the same sample. Press the "update" button to ensure changes affect C_q values (only needed for some analysis software). Switch to the linear view mode by clicking on the Y axis of the amplification chart, and check that the threshold previously set falls within the exponential phase of the curves.
- b) Set the baseline: determine the cycle number at which the threshold line crosses the first amplification curve and set the baseline three cycles before that value (e.g. earliest C_q = 25, set the baseline crossing at C_q = 25 – 3 = 22).
- c) Save the settings.
- d) Repeat the procedure described in a), b) and c) on the amplification plots of the other system (e.g. *hmg*).
- e) Save the settings and export all the data for further calculations.

3.4 Calculation of results

After having defined a threshold value within the logarithmic phase of amplification as described above, the instrument's software calculates the C_q values for each reaction.

The standard curves are generated both for the *hmg* and the MON 87411 specific assays by plotting the C_q values measured for the calibration points against the logarithm of the DNA copy numbers and by fitting a linear regression line into these data.

Thereafter, the standard curves are used to estimate the DNA copy number in the unknown samples.

To obtain the percentage value of event MON 87411 DNA in the unknown sample, the MON 87411 copy number is divided by the copy number of the maize endogenous gene *hmg* and multiplied by 100 (GM% = MON 87411/*hmg* × 100).

4. Equipment and Materials

4.1 Equipment

- Real-time PCR instrument for plastic reaction vessels (glass capillaries are not recommended for the described buffer composition) and appropriate analysis software
- 96-well reaction plates
- Optical caps/adhesion covers
- Microcentrifuge
- Micropipettes
- Standard bench top centrifuge with rotor or standard microfuge fit for 0.5 mL reaction tubes, centrifuge for 96-Well reaction plates
- Vortex
- Racks for reaction tubes, also cooled
- 0.5, 1.5 mL and 5 or 15 mL DNase free reaction tubes

4.2 Reagents

- TaqMan® Universal PCR Master Mix. Applied Biosystems Part No 4318157.

4.3 Primers and Probes

Table 5. Primers and probes for the MON 87411 and *hmg* methods

	Name	DNA Sequence (5' to 3')	Length (nt)
MON 87411			
Forward primer	87411 QF	CTC TGT AAC AGA AAA CAC CAT CTA GAG	27
Reverse primer	87411 QR	ACA AAA GTG AAC TAG TTC TAG GGT AGA T	28
Probe	87411 QP	6-FAM- CCG CGT TTA AAC TAT CAG TGT TTA GAG AAT-TAMRA	30
<i>hmg</i>			
Forward primer	<i>hmg</i> primer 1	TTG GAC TAG AAA TCT CGT GCT GA	23
Reverse primer	<i>hmg</i> primer 2	GCT ACA TAG GGA GCC TTG TCC T	22
Probe	<i>hmg</i> probe	6-FAM- CAA TCC ACA CAA ACG CAC GCG TA-TAMRA	23

FAM: 6-carboxyfluorescein; TAMRA: carboxytetramethylrhodamine;

5. References

1. Plant DNA C-values Database. Royal Botanic Gardens, Kew, <http://data.kew.org/cvalues/>

Page intentionally left blank

Europe Direct is a service to help you find answers to your questions about the European Union

Free phone number (*): 00 800 6 7 8 9 10 11

(*) Certain mobile telephone operators do not allow access to 00 800 numbers or these calls may be billed.

A great deal of additional information on the European Union is available on the Internet.

It can be accessed through the Europa server <http://europa.eu>

How to obtain EU publications

Our publications are available from EU Bookshop (<http://bookshop.europa.eu>), where you can place an order with the sales agent of your choice.

The Publications Office has a worldwide network of sales agents.

You can obtain their contact details by sending a fax to (352) 29 29-42758.

JRC Mission

As the science and knowledge service of the European Commission, the Joint Research Centre's mission is to support EU policies with independent evidence throughout the whole policy cycle.



EU Science Hub

ec.europa.eu/jrc



@EU_ScienceHub



EU Science Hub - Joint Research Centre



Joint Research Centre



EU Science Hub