



EUROPEAN COMMISSION  
DIRECTORATE GENERAL JRC  
JOINT RESEARCH CENTRE  
INSTITUTE FOR HEALTH AND CONSUMER PROTECTION  
COMMUNITY REFERENCE LABORATORY FOR GM FOOD AND FEED



# Event-specific Method for the Quantification of Rice Line LLRICE62 Using Real-time PCR

## Protocol

### **Method development:**

Bayer CropScience

### **Method validation:**

Joint Research Centre – European Commission  
Biotechnology & GMOs Unit

## Contents

<b>1. GENERAL INFORMATION AND SUMMARY OF THE METHODOLOGY.....</b>	<b>4</b>
<b>2. VALIDATION STATUS AND PERFORMANCE CHARACTERISTICS.....</b>	<b>4</b>
2.1 GENERAL .....	4
2.2 COLLABORATIVE TRIAL .....	5
2.3 LIMIT OF DETECTION.....	5
2.4 LIMIT OF QUANTITATION .....	5
2.5 MOLECULAR SPECIFICITY .....	5
<b>3. PROCEDURE .....</b>	<b>6</b>
3.1 GENERAL INSTRUCTIONS AND PRECAUTIONS .....	6
3.2 REAL-TIME PCR FOR QUANTITATIVE ANALYSIS OF RICE LLRICE62.....	7
<b>3.2.1 General</b> .....	7
<b>3.2.2 Calibration</b> .....	7
<b>3.2.3 Real-time PCR set-up</b> .....	7
3.3 DATA ANALYSIS .....	9
3.4 CALCULATION OF RESULTS .....	10
<b>4. MATERIALS.....</b>	<b>10</b>
4.1 EQUIPMENT .....	10
4.2 REAGENTS.....	11
4.3 PRIMERS AND PROBES .....	11
<b>5. REFERENCES .....</b>	<b>11</b>

<b><u>Document Approval</u></b>		
<b>Name / Function</b>	<b>Date</b>	<b>Signature</b>
<b>Marco Mazzara</b> <i>Sector Head</i>	<b>09/06/2006</b>	<b>Signed</b>
<b>Stephane Cordeil</b> <i>Quality Manager</i>	<b>09/06/2006</b>	<b>Signed</b>
<b>Guy Van den Eede</b> <i>B&amp;GMOs Unit Head</i>	<b>09/06/2006</b>	<b>Signed</b>

**Address of contact laboratory:**

European Commission, Joint Research Centre  
Institute for Health and Consumer Protection (IHCP)  
Biotechnology and GMOs Unit – Community Reference Laboratory  
Via Fermi 1, 21020 Ispra (VA) - Italy

## 1. General information and summary of the methodology

This protocol describes an event-specific real-time quantitative TaqMan<sup>®</sup> PCR procedure for the determination of the relative content of event LLRICE62 DNA to total rice DNA in a sample.

The PCR assay was optimised for the use in real-time PCR instruments for plastic reaction vessels. Glass capillaries are not recommended for the buffer composition described in this method.

Template DNA extracted by means of suitable methods should be tested for quality and quantity prior to the use in PCR assay. Tests for the presence of PCR inhibitors (e.g. monitor run, use of DNA spikes) are recommended.

For specific detection of event LLRICE62 genomic DNA, an 88-bp fragment of the recombination area between the insert and the plant genome (located at the 3' flanking region) is amplified using two specific primers. PCR products are measured during each cycle (real-time) by means of a target-specific oligonucleotide probe labelled with two fluorescent dyes: FAM as a reporter dye at its 5' end and TAMRA as a quencher dye at its 3' end.

For relative quantitation of event LLRICE62 DNA, a rice-specific reference system amplifies a 64-bp fragment of the *phospholipase D* gene (PLD), a rice endogenous gene, using a pair of specific primers and a PLD gene-specific probe labelled with FAM and TAMRA as described above.

The measured fluorescence signal passes a threshold value after a certain number of cycles. This threshold cycle is called the "Ct" value. For quantitation of the amount of event LLRICE62 DNA in a test sample, the normalized  $\Delta$ Ct values of the calibration samples are used to calculate by linear regression a reference curve  $\Delta$ Ct-formula. The normalized  $\Delta$ Ct values of the unknown samples are measured and, by means of the reference  $\Delta$ Ct-formula, the relative amount of LLRICE62 event DNA is estimated.

## 2. Validation status and performance characteristics

### 2.1 General

The method has been optimised for suitable DNA extracted from rice leaves, grains or seeds, containing mixtures of genetically modified and conventional rice.

The reproducibility and trueness of the method was tested through collaborative trial using samples at different GMO contents.

## **2.2 Collaborative trial**

The method was validated in a collaborative trial by the Joint Research Centre (JRC) of the European Commission. The study was undertaken with 11 laboratories.

Each participant received twenty unknown samples containing LLRICE62 rice genomic DNA at five GM contents, between 0.15 % and 3.3 %.

Each test sample was analyzed by PCR in three repetitions. The study was designed as a blind quadruplicate collaborative trial; each laboratory received each level of GM LLRICE62 in four unknown samples. Two replicates of each GM level were analyzed on the same PCR plate.

A detailed validation report can be found under <http://gmo-crl.jrc.it/statusofdoss.htm>

## **2.3 Limit of detection**

According to the method developer, the relative LOD of the method is at least 0.045% in 200 ng of total rice DNA. The relative LOD was not assessed in a collaborative trial. The lowest relative GM content of the target sequence included in collaborative trial was 0.15%.

## **2.4 Limit of quantitation**

According to the method developer, the relative LOQ of the method is at least 0.09% in 200 ng of total rice DNA. The lowest relative GM content of the target sequence included in collaborative trial was 0.15%.

## **2.5 Molecular specificity**

The method utilizes a unique DNA sequence of the recombination region between the insert and the plant genome. The sequence is specific to LLRICE62 event and thus imparts event-specificity to the method.

The specificity was experimentally tested against DNA extracted from plant materials containing the specific targets of maize T25, MON 810, Bt11, Bt176, GA21, NK603, CBH351, Roundup Ready Soybean, OSR Ms1, Ms8, Rf1, Fr2, Rf3, Topas19/2, T45, soybean A2704/12 and LLCotton25.

None of the materials yielded detectable amplification in replicate experiments.

### **3. Procedure**

#### **3.1 General instructions and precautions**

- All handling of reagents and controls should occur in an ISO 9001:2000 or ISO 17025 environment or equivalent.
- The procedures require experience of working under sterile conditions.
- Laboratory organization, e.g. “flow direction” during PCR-setup, should follow the guidelines given by relevant authorities like e.g. ISO, CEN, Codex alimentarius commission.
- PCR-reagents shall be stored and handled in a separate room and freezer and in equipment where no nucleic acids (with exception of PCR primers or probes) or DNA degrading or modifying enzymes have been handled previously. All handling of PCR reagents and controls requires dedicated equipment – especially pipettes.
- All the equipment used must be sterilized prior to use and any residue of DNA has to be removed. All material used (e.g. vials, containers, pipette tips, etc.) must be suitable for PCR and molecular biology applications. They must be DNase-free, DNA-free, sterile and shall not adsorb protein or DNA.
- In order to avoid contamination, filter pipette tips protected against aerosol should be used.
- Use only powder-free gloves and change them frequently.
- Clean lab-benches and equipment periodically with 10% sodium hypochloride solution (bleach).
- Pipettes should be checked regularly for precision and calibrated, if necessary.
- All handling steps - unless specified otherwise - shall be carried out at 0 - 4°C.
- In order to avoid repeated freeze/thaw cycles aliquots should be prepared.

## 3.2 Real-time PCR for quantitative analysis of rice LLRICE62

### 3.2.1 General

The PCR set-up for the taxon specific target sequence (PLD) and for the GMO (LLRICE62) target sequence should be carried out in separate vials. Multiplex PCR (using differential fluorescent labels for the probes) has not been tested or validated.

The use of maximum 200 ng of template DNA per reaction well is recommended.

The method is developed for a total volume of 25 µl per reaction mixture with the reagents as listed in Table 1 and Table 2.

### 3.2.2 Calibration

The calibration curve consists of five samples containing fixed percentages of LLRICE62 DNA in a total amount of 200 ng rice DNA. The GM content of the standard samples ranges from 3.6% to 0.09%.

A calibration curve is produced by plotting the  $\Delta C_t$ -values of calibration samples against the logarithm of the respective GM % contents; the slope (a) and the intercept (b) of the standard curve ( $y = ax + b$ ) are then used to calculate the mean % GM content of the blind samples based on their normalised  $\Delta C_t$  values.

### 3.2.3 Real-time PCR set-up

1. Thaw, mix gently and centrifuge the required amount of components needed for the run. **Keep thawed reagents at 1-4°C on ice.**
2. In two reaction tubes (one for LLRICE62 system and one for the PLD system) on ice, add the following components (Tables 1 and 2) in the order mentioned below (except DNA) to prepare the master mixes.

**Table 1.** Amplification reaction mixture in the final volume/concentration per reaction well for the reference PLD specific system.

<b>Component</b>	<b>Final concentration</b>	<b>µl/reaction</b>
TaqMan <sup>®</sup> Universal PCR Master Mix (2x)	1x	12.5
KVM159 primer (10 µM)	200 nM	0.5
KVM160 primer (10 µM)	200 nM	0.5
TM013 probe (10 µM)	200 nM	0.5
Nuclease free water	#	6
Template DNA (max 200 ng, see 3.2.1 and 3.2.2)	#	5
<b>Total reaction volume:</b>		<b>25</b>

**Table 2.** Amplification reaction mixture in the final volume/concentration per reaction well for the LLRICE62 specific system.

<b>Component</b>	<b>Final concentration</b>	<b>µl/reaction</b>
TaqMan <sup>®</sup> Universal PCR Master Mix (2x)	1x	12.5 µl
MDB616 primer (10 µM)	400 nM	1
MDB694 primer (10 µM)	400 nM	1
TM019 Probe (10 µM)	200 nM	0.5
Nuclease free water	#	5
Template DNA (max 200 ng, see 3.2.1 and 3.2.2)	#	5
<b>Total reaction volume:</b>		<b>25 µl</b>

- Mix gently and centrifuge briefly.
- Prepare two reaction tubes (one for the LLRICE62 and one for the PLD master mix) for each DNA sample to be tested (reference curve samples, unknown samples and control samples).
- Add to each reaction tube the correct amount of master mix (e.g. 20 x 3 = 60 µl master mix for three PCR repetitions). Add to each tube the correct amount of DNA



(e.g.  $5 \times 3 = 15 \mu\text{l}$  DNA for three PCR repetitions). Low-speed vortex each tubes at least three times for approx 30 sec. This step is mandatory to reduce the variability among the repetitions of each sample to a minimum.

6. Spin down the tubes in a micro-centrifuge. Aliquot 25  $\mu\text{l}$  in each well. Seal the reaction plate with optical cover or optical caps. Centrifuge the plate at low speed (e.g. approximately  $250 \times g$  for 1 minute at 4 °C to room temperature) to spin down the reaction mixture.
7. Place the plate into the instrument.
8. Run the PCR with cycling conditions described in Table 3:

**Table 3.** Cycling program for maize LLRICE62/PLD systems

Step	Stage	T°C	Time (sec)	Acquisition	Cycles
1	UNG	50 °C	120	No	1
2	Initial denaturation	95 °C	600"	No	1
3	Denaturation	95 °C	15"	No	45
	Amplification Annealing & Extension	60 °C	60"	Measure	

### 3.3 Data analysis

Subsequent to the real-time PCR, analyse the run following the procedure below:

- a) Set the threshold: display the amplification curves of one system (e.g. LLRICE62) in logarithmic mode. Locate the threshold line in the area where the amplification profiles are parallel (exponential phase of PCR) and where there is no "fork effect" between repetitions of the same sample. Press the update button to ensure changes affect Ct values. Switch to the linear view mode by clicking on the Y axis of the amplification plot, and check that the threshold previously set falls within the geometric phase of the curves.
- b) Set the baseline: determine the cycle number at which the threshold line crosses the first amplification curve and set the baseline three cycles before that value (e.g. earliest Ct = 25, set the baseline crossing at Ct =  $25 - 3 = 22$ ).

- c) Save the settings
- d) Repeat the procedure described in a) and b) on the amplification plots of the other system (e.g. PLD system).
- e) Save the settings and export all the data into an Excel file for further calculations.

### **3.4 Calculation of results**

After having defined a threshold value within the logarithmic phase of amplification as described above, the instrument's software calculates the Ct-values for each reaction.

The Reference  $\Delta$ Ct-curve is generated by plotting the  $\Delta$ Ct-values measured for the calibration points against the logarithm of the GM% content, and by fitting a linear regression line into these data.

Thereafter, the reference  $\Delta$ Ct-curve formula is used to estimate the relative amount (%) of LLRICE62 event in the unknown samples of DNA.

## **4. Materials**

### **4.1 Equipment**

- Real-time PCR instrument for plastic reaction vessels (glass capillaries are not recommended for the described buffer composition)
- Plastic reaction vessels suitable for real-time PCR instrument (enabling undisturbed fluorescence detection)
- Software for run analysis (mostly integrated in the software of the real-time PCR instrument)
- Microcentrifuge
- Micropipettes
- Vortex
- Rack for reaction tubes
- 1.5/2.0 ml tubes

## 4.2 Reagents

- TaqMan® Universal PCR Master Mix (2X). Applied Biosystems Part No 4304437

## 4.3 Primers and Probes

Name	Oligonucleotide DNA Sequence (5' to 3')
<i>LLRICE62 target sequence</i>	
MDB616	5'-AGC TGG CGT AAT AGC GAA GAG G-3'
MDB694	5'-TGC TAA CGG GTG CAT CGT CTA-3'
TM019 (Probe)	FAM-5'-CGC ACC GAT TAT TTA TAC TTT TAG TCC ACC T-3'- TAMRA
<i>Reference gene PLD target sequence</i>	
KVM159	5'- TGG TGA GCG TTT TGC AGT CT-3'
KVM160	5'- CTG ATC CAC TAG CAG GAG GTC C-3'
TM013 (Probe)	FAM-5'- TGT TGT GCT GCC AAT GTG GCC TG - 3'-TAMRA

## 5. References

Arumuganathan, K., Earle, E.D. (1991). Nuclear content of some important plant species. Plant Mol Biol Reporter 9, 208-218.