

Validation of TaqMan screening qPCR method specific to the GMM protease1 in targeting the left transgene flanking region between the pUB110 shuttle vector and a *Bacillus velezensis* gene coding for a protease (GMM protease1 left border)

This dossier describes the validation of a TaqMan qPCR [construct-specific](#) method for the detection of an unauthorized GMM producing protease (GMM protease1). This method targets the left transgene flanking region between the pUB110 shuttle vector and a *Bacillus velezensis* gene coding for a protease.

The method is intended for use to detect the above mentioned GMM producing protease in food and feed microbial fermentation products.

Author: Marie-Alice Fraiture; Verified by: Nina Papazova and Nancy Roosens

1 qPCR method

The method described here uses the TaqMan® qPCR chemistry for amplification of a 94 bp fragment of the left transgene flanking region of the GMM producing protease (Fraiture et al., 2020).

1.1. *Oligonucleotides*

The primer and probe sequences are shown in table 1.

Table 1. Primer and probe sequences.

Oligonucleotide type	Oligonucleotide name	Sequence (5'-3')
Forward primer	GMM_protease1_left_border-F	CGAGAATGCAGCTGAAACAG
Reverse primer	GMM_protease1_left_border-R	CATATGCTCGGGGAATTTATCT
Probe	GMM_protease1_left_border-P	FAM-GGACGGACAGATCAAGAACTGTTATGG-TAMRA

1.2. *qPCR mix*

The qPCR mix is given in table 2.

Table 2. qPCR mix for the GMM protease1 left border marker.

Component	Stock concentration	Final concentration	µl per reaction
TaqMan mastermix	2X	1 X	12.50 µl
Primer GMM_protease1_left_border-F	20 µM	400 nM	0.50 µl
Primer GMM_protease1_left_border-R	20 µM	400 nM	0.50 µl
Probe GMM_protease1_left_border-P	10 µM	200 nM	0.50 µl
Nuclease-free water			6.00 µl
		Total volume	20 µl
		DNA 5 ng/µl	5 µl

1.3. qPCR programme

The qPCR programme is given in table 3.

Table 3. qPCR programme

Step	Stage	T (°C)	Time (s)	Acquisition	Cycles
1	Activation UNG	50	120	No	1x
2	Taq Activation	95	600	No	1x
3	Amplification Denaturing Annealing and Extension	95 60	15 60	No Mesure	45x

1.4. Control materials

No Certified Reference Materials (CRM) or other Reference Materials (RM) are commercially available for this method. As positive control material, gDNA extracted from an isolated strain of GM *B. velezensis* producing protease, related to RASFF2019.3332 (isolate indicated as 2019.3332) is used. In addition, a plasmid, artificially synthesized by geneCust, carrying one copy of the target sequence (Figure 1) can be used.

2 Scope of the validation

The TaqMan qPCR method is intended for use as a qualitative detection method. Therefore, the validation parameters tested are the specificity and the limit of detection (LOD_{95%}), as described in Fraiture et al., 2020.

3 Method development

The GMM protease left border marker was developed by Service TAG in Sciensano within the SPECENZYM project and is described in Fraiture et al. 2020. This TaqMan qPCR method is targeting the left transgene flanking region of an unauthorized GMM producing protease. More precisely, the target sequence corresponds to the junction between the pUB110 shuttle vector and a gene coding for a protease belonging to *Bacillus velezensis* species.

3.1. Experimental set-up for testing *in situ* specificity

3.1.1 Bacterial strains and isolates

The specificity tests were done using positive and negative materials:

As positive materials, three samples were used: 1/ gDNA from isolate 2019-3332 related to RASFF2019.3332 (GM *B. velezensis* RASFF2019.3332 strain); 2/ DNA extract from a food enzyme product (labelled as containing protease; Sample 1) contaminated by the GM *B. velezensis* RASFF2019.3332 strain; 3/ a plasmid, artificially synthesized by genecust (<https://www.genecust.com/en/>) carrying one copy of the targeted sequence (Figure 1) (PC plasmid).

Different types of negative materials were used: 1/ WT bacterial species reported as used for the production of food and feed fermentation products; 2/ WT fungal species reported as used for the production of food and feed fermentation products; 3/ GM *Bacillus subtilis* RASFF2014.1249 strain producing vitamin B2.

In addition, a no template control (NTC = water) was used.

The list of all tested materials is given in table 4.

Table 4: Materials used to test the specificity of the GMM protease1 left border qPCR method. All materials were used at 10 ng, excepted the PC plasmid, carrying one copy of the target sequence, that was used at 100 copies of the target (see Table 6).

Species	Origin	Reference
Positive samples		
<i>Bacillus velezensis</i> GM (2019-3332)	/	2019-3332
gDNA extracted from a food enzyme product, labelled as containing protease. This "real-life sample" was reported as being contaminated by the GMM producing protease1 (RASFF2019.3332).	/	sample 1
Positive control plasmid carrying one copy of the target sequence	/	PC plasmid carrying one copy of the target sequence
Negative samples		
Fungal strains		
<i>Aspergillus acidus</i>	BCCM	IHEM 26285
<i>Aspergillus aculeatus</i>	BCCM	IHEM 05796
<i>Aspergillus fijiensis</i>	BCCM	IHEM 22812
<i>Aspergillus melleus</i>	BCCM	IHEM 25956
<i>Aspergillus niger</i>	BCCM	IHEM 25485
<i>Aspergillus oryzae</i>	BCCM	IHEM 25836
<i>Candida cylindracea</i>	BCCM	MUCL 041387
<i>Candida rugose</i>	BCCM	IHEM 01894
<i>Chaetomium gracile</i>	BCCM	MUCL 053569
<i>Cryphonectria parasitica</i>	BCCM	MUCL 007956
<i>Disporotrichum dimorphosporum</i>	BCCM	MUCL 019341
<i>Fusarium venenatum</i>	BCCM	MUCL 055417
<i>Hansenula polymorpha</i>	BCCM	MUCL 027761
<i>Humicola insolens</i>	BCCM	MUCL 015010
<i>Kluyveromyces lactis</i>	BCCM	IHEM 02051
<i>Leptographium procerum</i>	BCCM	MUCL 008094

Species	Origin	Reference
<i>Mucor javanicus</i>	BCCM	IHEM 05212
<i>Penicilliumcamemberti</i>	BCCM	IHEM 06648
<i>Penicilliumchrysogenum</i>	BCCM	IHEM 03414
<i>Penicilliumcitrinium</i>	BCCM	IHEM 26159
<i>Penicilliumdecumbens</i>	BCCM	IHEM 05935
<i>Penicilliumfuniculosum</i>	BCCM	MUCL 014091
<i>Penicilliummulticolor</i>	CBS	CBS 501.73
<i>Penicilliumroqueforti</i>	BCCM	IHEM 20176
<i>Pichia pastoris</i>	BCCM	MUCL 027793
<i>Rhizomucor miehei</i>	BCCM	IHEM 26897
<i>Rhizopus niveus</i>	ATCC	ATCC 200757
<i>Rhizopus oryzae</i>	BCCM	IHEM 26078
<i>Saccharomyces cerevisiae</i>	BCCM	IHEM 25104
<i>Sporobolomyces singularis</i>	BCCM	MUCL 027849
<i>Talaromyces cellulolyticus/pinophilus</i>	BCCM	IHEM 16004
<i>Talaromyces emersonii</i>	DSMZ	DSMZ 2432
<i>Trametes hirsuta</i>	BCCM	MUCL 030869
<i>Trichoderma citrinoviride</i>	BCCM	IHEM 25858
<i>Trichoderma longibrachiatum</i>	BCCM	IHEM 00935
<i>Trichoderma reesei</i>	BCCM	IHEM 05651
<i>Trichoderma viride</i>	BCCM	IHEM 04146
<u>Bacterial strains</u>		
<i>Arthrobacter ramosus</i>	BCCM	LMG 17309
<i>Bacillus amyloliquefaciens</i>	BCCM	LMG 98140
<i>Bacillus brevis</i>	BCCM	LMG 7123
<i>Bacillus cereus</i>	ATCC	ATCC 14579
<i>Bacillus circulans</i>	BCCM	LMG 6926T
<i>Bacillus coagulans</i>	BCCM	LMG 6326
<i>Bacillus firmus</i>	BCCM	LMG 7125
<i>Bacillus flexus</i>	BCCM	LMG 11155
<i>Bacillus lentus</i>	Sciensano	TIAC 101
<i>Bacillus licheniformis</i>	BCCM	LMG 6933T
<i>Bacillus megaterium</i>	BCCM	LMG 7127
<i>Bacillus pumilus</i>	DSMZ	DSMZ 1794
<i>Bacillus smithii</i>	BCCM	LMG 6327
<i>Bacillus subtilis</i>	BCCM	LMG 7135 T
<i>Bacillus subtilis</i>	Sciensano	W04-510
<i>Bacillus subtilis</i>	Sciensano	E07-505
<i>Bacillus subtilis</i>	Sciensano	S10005
<i>Bacillus subtilis</i>	Sciensano	SUB033
<i>Bacillus subtilis</i>	Sciensano	BNB54

Species	Origin	Reference
<i>Bacillus velezensis</i>	BCCM	LMG 12384
<i>Bacillus velezensis</i>	BCCM	LMG 17599
<i>Bacillus velezensis</i>	BCCM	LMG 22478
<i>Bacillus velezensis</i>	BCCM	LMG 23203
<i>Bacillus velezensis</i>	BCCM	LMG 26770
<i>Bacillus velezensis</i>	BCCM	LMG 27586
<i>Cellulosimicrobiumcellulans</i>	BCCM	LMG 16121
<i>Corynebacteriumglutamicum</i>	BCCM	LMG 3652
<i>Enterococcus faecium</i>	BCCM	LMG 9430
<i>Escherichia coli</i>	BCCM	LMG2092T
<i>Geobacillus caldoproteolyticus</i>	DSMZ	DSMZ 15730
<i>Geobacillus pallidus</i>	BCCM	LMG 11159T
<i>Geobacillus stearothermophilus</i>	BCCM	LMG 6939T
<i>Klebsiella pneumonia</i>	BCCM	LMG 3113T
<i>Lactobacillus casei</i>	BCCM	LMG 6904
<i>Lactobacillus fermentum</i>	BCCM	LMG 6902
<i>Lactobacillus plantarum</i>	BCCM	LMG 9208
<i>Lactobacillus rhamnosus</i>	BCCM	LMG 18030
<i>Lactococcus lactis</i>	BCCM	LMG 6890T
<i>Leuconostoc citreum</i>	BCCM	LMG 9824
<i>Microbacteriumimperiale</i>	BCCM	LMG 20190
<i>Paenibacillus alginolyticus</i>	BCCM	LMG 18723
<i>Paenibacillus mæcerans</i>	BCCM	LMG 6324
<i>Protaminobacter pubrum</i>	CBS	CBS 574.77
<i>Pseudomonas amyloclavata</i>	ATCC	ATCC-21262
<i>Pseudomonas fluorescens</i>	BCCM	LMG1794T
<i>Pullulanibacillus naganoensis</i>	BCCM	LMG 12887
<i>Streptomyces aureofaciens</i>	BCCM	LMG 5968
<i>Streptomyces mobaraensis</i>	DSMZ	DSMZ 40847
<i>Streptomyces murinus</i>	BCCM	LMG 10475
<i>Streptomyces netropsis</i>	BCCM	LMG 5977
<i>Streptomyces rubiginosus</i>	BCCM	LMG20268
<i>Streptomyces violaceoruber</i>	BCCM	LMG 7183
<i>Streptoverticilliummobaraense</i>	CBS	CBS 199.75
GMM samples		
GM <i>Bacillus subtilis</i> producing vitamin B2 (RASFF2014.1249)	SCL	GMM producing vitamin B2
<i>gDNA extracted from a food enzyme product, labelled as containing protease. This "real-life sample" was reported as being contaminated by the GMM producing protease2 (RASFF2021.1641).</i>	Sciensano	GMM producing protease2
PCR controls		
No template control (NTC)		

BCCM: Belgian co-ordinated collections of micro-organisms (<http://bccm.belspo.be/>).
 CBS: Convention of Biological Diversity (<http://www.wi.knaw.nl/Collections/>).

DSMZ: Deutsche Sammlung von Mikroorganismen und Zellkulturen GmbH; (<https://www.dsmz.de/>).
 ATCC: American Type Culture Collection (<http://www.lgcstandards-atcc.org/>) ;
 Sciensano: Service "Foodborne pathogenes, J. Wytsmanstraat 14, 1050 Brussels.
 SCL : Service Commun des Laboratoires, France

3.1.2 Equipment:

qPCR: CFX96 Touch Real-Time PCR Detection System (BioRad, BE).

3.1.3 Reagents:

qPCR: Taq[®]Man 2x Master Mix (Diagenode BE, ref: GMO-MM2X-A300), Nuclease free water (Acros organics, BE), primers and probe (Eurogentec, BE).

3.1.4 Bacterial DNA isolation and measurement

Bacterial gDNA from GM *B. subtilis* and *velezensis* isolates was extracted as described in Fraiture et al. 2020. The wild-type microbial strains were provided as gDNA extracted by using "Quick-DNA™ Fungal/Bacterial Miniprep Kit" (Zymo, cat. No. D6005).

The DNA concentration was measured spectrophotometrically using Nanodrop® 2000 (ThermoFisher) and DNA purity was evaluated using the A260/A280 and A260/A230 ratios.

3.1.5 Taq[®]Man qPCR:

The qPCR conditions are described in part 1.1. – 1.3. The reactions were carried out in two replicates for each sample.

3.1.6 Analysis of the amplicon sequence

The PCR product from the GMM protease1 left border marker was purified using USB ExoSAP-IT PCR Product Cleanup (Affymetrix) to be sequenced on a Genetic Sequencer 3500 (ThermoFisher) using the Big Dye Terminator Kit v3.1 (Applied Biosystems). The amplicon sequence is shown on Figure 1.

**CGAGAATGCAGCTGAAACAGCACTTGGCGGAAAAAAGGAAGGACGGACAGATCAAGAACTGT
 TATGGCTACAAGATAAATTCCCCGAGCATATG**

Figure. 1 Amplicon sequence of GMM protease1 left border qPCR method. The primers and probe are indicated in bold.

3.1.7 Data processing

The following analyses of the data obtained by the qPCR specificity testing were done based on:

- Calculation of the false negatives rate. The false negatives rate was calculated according to the formula:

$$\text{False neg rate} = \frac{100 \times (\text{No. misclassified known pos samples})}{\text{total No. known pos samples}}$$

- Calculation of the false positives rate. The false positives rate was calculated according to the formula:

$$\text{False pos rate} = \frac{100 \times (\text{No. misclassified known neg samples})}{\text{total No. known neg samples}}$$

The GMM protease1 left border qPCR method is considered specific when the following conditions are obtained:

- Positive amplification signal in materials containing the targeted sequence from the left transgene flanking region of GM *B. velezensis* RASFF2019.3332 strain with a false negative rate of 0%.
- No positive amplification signal in materials that do not contain the targeted sequence from the left transgene flanking region of GM *B. velezensis* RASFF2019.3332 strain with a false positives rate of 0%.

3.2. Results

3.2.1 Specificity of the qPCR method

The positive materials (3 samples) gave a positive amplification signal. No misclassified positive samples were observed. The false negatives rate was 0%.

All negative materials (91 samples) and the negative PCR control (1 sample) gave no amplification signal. No misclassified negative samples were observed. The false positives rate was 0%.

The C_t values obtained for the positive materials is given in table 5. As no C_t was recorded for the negative materials, it is not included in the table 5. The full C_t dataset and amplification curves are shown in Annex 1.

Table 5 Results from the specificity test of the qPCR method for the positive materials. Information about the estimated target copy number is not available for the sample 1, as it is DNA extracted from a food enzyme product contaminated with an unknown amount of GMM protease 1. For the PC plasmid, the test was performed on 100 estimated target copies. For the 2019.3332 isolate, the test was performed on 2,200,000 estimated genome copies.

Bacterial species	C_t value
<i>Bacillus velezensis</i> GM (2019-3332)	15.5
sample 1	12.7
PC plasmid carrying one copy of the target sequence	33.2

3.3. Conclusions on specificity

The experimental results on specificity showed that:

- The false positives rate is 0%.
- The false negatives rate is 0%.

Based on this, it can be concluded that the GMM protease 1 left border qPCR method is specific for the detection of the unauthorized GM *Bacillus velezensis* producing protease (GMM protease1).

4 Limit of detection (LOD_{95%}) of the GMM protease 1 left border qPCR method

4.1. Experimental set-up

4.1.1 Material:

Plasmid DNA synthesized by genecust is carrying the target sequenced (Figure 1). The plasmid was used under the native/supercoiled form.

4.1.2 Preparation of the serial dilutions

The plasmid DNA was first diluted in nuclease free water (Acros organics, BE) in order to obtain 1,595 copies template in the final reaction.

The copy number calculation was done according to the formula hereunder taking into account that the

plasmid size is 2,903 bp:

$$N = \frac{m \times N_A}{MW \times L}$$

N is the number of plasmid molecules, m is the plasmid amount in g, N_A is Avogadro constant = $6,0221415 \times 10^{23} \text{ mol}^{-1}$. MW is average base pair weight = 649Da and L is the plasmid length in base pairs.

On this basis, dilution series starting from 100 copies to 0.1 copy as indicated in Table 6, were prepared in nuclease free water (Acros organics, BE). To determine the $\text{LOD}_{95\%}$, 12 replicates of each dilution indicated in Table 6 were tested.

Table 6: Preparation of dilution series.

DNA concentration (ng)	Copy number	Dilution factor	Volume from the previous dilution (µl)	Water (µl)	Total volume (µl)
0.000005	1595				
3.125E-07	100	16	100	1500	1600
6,25E-08	20	5	200	800	1000
3.125E-08	10	2	300	300	600
1.5625E-08	5	2	200	200	400
3.125E-09	1	5	100	400	500
3.125E-10	0.1	10	100	900	1000

4.1.3 Equipment

qPCR: CFX96 Touch Real-Time PCR Detection System (BioRad).

4.1.4 Reagents

qPCR: Taq[®]Man 2x Master Mix (Diagenode BE, ref: GMO-MM2X-A300), Nuclease free water (Acros organics, BE), primers and probe (Eurogentec, BE).

4.1.5 Taq[®]Man qPCR

The qPCR conditions are described part 1.1 – 1.3.

4.1.6 Data analysis

The $\text{LOD}_{95\%}$, defined as the number of copies of the target required to ensure a 95% probability of detection (POD), was determined by using Quodata web tool (Uhlig et al. 2015; Grohmann 2016; <https://quodata.de/content/validation-qualitative-pcr-methods-single-laboratory>) accessed in November 2020.

4.2. Results

The $\text{LOD}_{95\%}$ of the method was determined as described in part 4.1.6. The summarized data are shown in table 7. The full C_t data set and the amplification plots are shown in Annex 2. As indicated in Figure 2, the $\text{LOD}_{95\%}$ was set at 3 copies using DNA from a plasmid carrying one copy of the target sequence, complying with the acceptance criteria in the “Minimum Performance Requirements for Analytical Methods of GMO Testing” of the European Network of GMO Laboratories (ENGL, 2015). The plausibility check of the curve did not give any irregularities.

Table 7: Number of positive replicates and average C_t obtained for each dilution tested.

Copy number	Positive replicate number	Average C_t
100	12/12	33.2

20	12/12	35.7
10	12/12	36.8
5	12/12	38.2
1	8/12	40.2
0.1	0/12	/
0	0/12	/

Please enter your data

1st column: Number of copies of the target DNA sequence (= nominal copies)

2nd column: Number of positive test results

3rd column: Number of PCR replicates

[Click here to insert example data](#)

0.1	0	12
1	8	12
5	12	12
10	12	12
20	12	12
100	12	12

[Clear data](#)

Results

Plausibility Check

The plausibility check indicates no irregularities.

POD curve and LOD_{95%}

The LOD_{95%} is 3.037 with a 95 % confidence interval of [1.645, 5.594].

The figure below summarises the results. The blue diamonds characterise the laboratory-specific RODs. The blue curve denotes the mean POD curve along with the corresponding 95 % confidence range highlighted as the grey band. The POD curve under ideal conditions is displayed as the black dashed curve.




Figure 2 Calculation of LOD_{95%} based on the POD curve.

5 Final conclusions on the in-house validation

The results from the specificity and sensitivity testing showed that the tested parameters comply with the acceptance criteria described in the “Minimum Performance Requirements for Analytical Methods of GMO Testing” of the European Network of GMO Laboratories (ENGL, 2015).

As control material, gDNA extracted from the *B. velezensis* strain 2019-3332 as well as plasmid DNA harboring the targeted sequence (PC plasmid for GM protease1 left border marker) are recommended to be used. gDNA from the *B. velezensis* strain 2019-3332 is available in Sciensano. The plasmid was artificially synthesized by Genecust (<https://www.genecust.com/en/>). This plasmid consists in pUC57 where the sequence of interest (Figure 1) was introduced.

6 Method transferability

In order to evaluate whether the in-house validated GMM protease1 left border marker performs in the same manner when performed in another laboratory with different operators using different equipment and reagents, a transferability study was carried out. The experiments were performed in Unità Operativa Semplice a valenza Direzionale - Ricerca e controllo degli organismi geneticamente modificati at the Istituto Zooprofilattico Sperimentale del Lazio e della Toscana "M. Aleandri" (Roma, Italy).

6.1. *Experimental set-up*

The experimental setup was the same as described in section 4, for determination of the LOD_{95%}. Serial dilutions of the control plasmid DNA, from 100 to 0.1 estimated target copy number, used in-house was tested with the same real-time PCR protocol, were sent to the second laboratory. All reagents were ordered by the external laboratory, including oligonucleotides (Metabion International AG) and TaqMan™ Universal PCR Master Mix (Applied Biosystems®, 4304437). All runs were performed on a QuantStudio™ 7 Flex Real-Time PCR System (Life Technologies). The LOD_{95%} was determined based on the results of the second laboratory as described in the section 4.1.6.

6.2. *Results*

The summarized data are shown in table 8. The full C_t data set and the amplification plots are shown in Annex 3. [The number of the positive replicates per dilution was comparable to the results presented in Table 7.](#)

As indicated in Figure 3, the LOD_{95%} was set at 2 copies using DNA from a plasmid carrying one copy of the target sequence, complying with the acceptance criteria in the "Minimum Performance Requirements for Analytical Methods of GMO Testing" of the European Network of GMO Laboratories (ENGL, 2015).

Table 8: Number of positive replicates and average Ct obtained for each dilution tested.

Copy number	Positive replicate number	Average C _t
100	12/12	34.3
20	12/12	36.6
10	12/12	37.5
5	12/12	39.1
1	10/12	40.2
0.1	2/12	41.0
0	0/12	/

Please enter your data

1st column: Number of copies of the target DNA sequence (= nominal copies)

2nd column: Number of positive test results

3rd column: Number of PCR replicates

[Click here to insert example data](#)

0.1	2	12
1	10	12
5	12	12
10	12	12
20	12	12
100	12	12

[Clear data](#)

Start the calculation

Results

Plausibility Check

The calculated POD curve indicates sensitivity better than achievable according to the theoretical POD curve.

POD curve and LOD_{95%}

The LOD_{95%} is 1.664 with a 95 % confidence interval of [0.889, 3.108].

The figure below summarises the results. The blue diamonds characterise the laboratory-specific RODs. The blue curve denotes the mean POD curve along with the corresponding 95 % confidence range highlighted as the grey band. The POD curve under ideal conditions is displayed as the black dashed curve.




Figure 3 Calculation of LOD_{95%} based on the POD curve.

6.3. Conclusion on the method transferability

The results obtained by the external laboratory were comparable with the ones observed during the in-house validation. These results demonstrate that the method can be successfully performed in another laboratory using different operators, infrastructure, equipment and reagents.

7 References

Deckers M, Vanneste K, Winand R, De Keersmaecker SCJ, Denayer S, Heyndrickx M, Deforce D, Fraiture MA, Roosens NHC (2020a) Strategy for the identification of micro-organisms producing food and feed products: Bacteria producing food enzymes as study case. Food Chemistry 305:125431.

Deckers M, Vanneste K, Winand R, Hendrickx M, Becker P, De Keersmaecker SCJ, Deforce D, Fraiture MA, Roosens NHC (2020b) Screening strategy targeting the presence of food enzyme-producing fungi in food enzyme preparations. *Food Control* 117:107295.

ENGL (2015) Definition of Minimum Performance Requirements for Analytical Methods for GMO Testing. http://gmo-crl.jrc.ec.europa.eu/doc/MPR%20Report%20Application%2020_10_2015.pdf.

Fraiture et al. 2020. Next-generation sequencing: a key tool to identify unauthorized genetically modified microorganisms in food enzyme preparations. *Sci. Rep.* 10:7094.

Grohmann 2016; Guidelines for the validation of qualitative real-time PCR methods by means of a collaborative study.

<https://quodata.de/content/validation-qualitative-pcr-methods-single-laboratory>

International Standard: ISO5725-2 (1994). Accuracy (trueness and precision) of measurement methods and results - Part 2. International Organisation for Standardization, Genève, Switzerland.

Uhlig, S., Frost, K., Colson, B. et al. Validation of qualitative PCR methods on the basis of mathematical-statistical modelling of the probability of detection. *Accred Qual Assur* 20, 75–83 (2015). <https://doi.org/10.1007/s00769-015-1112-9>


CONFIDENTIAL

8 Annexes

Annex 1: Specificity: full C_t dataset and amplification plots. DNA extracts were previously assessed as amplifiable (Deckers et al., 2020a,b; Fraiture et al., 2020).

Samples		Ct values	
Positive samples			
<i>Bacillus velezensis</i> GM (2019-3332)		15.5	15.5
Sample 1		12.7	12.7
PC plasmid carrying one copy of the target sequence		33.2	33.4
Negative samples			
Fungal strains			
<i>Aspergillus acidus</i>	IHEM 26285	No signal	No signal
<i>Aspergillus aculeatus</i>	IHEM 05796	No signal	No signal
<i>Aspergillus fijiensis</i>	IHEM 22812	No signal	No signal
<i>Aspergillus melleus</i>	IHEM 25956	No signal	No signal
<i>Aspergillus niger</i>	IHEM 25485	No signal	No signal
<i>Aspergillus oryzae</i>	IHEM 25836	No signal	No signal
<i>Candida cylindracea</i>	MUCL 041387	No signal	No signal
<i>Candida rugosa</i>	IHEM 01894	No signal	No signal
<i>Chaetomium gracile</i>	MUCL 053569	No signal	No signal
<i>Cryphonectria parasitica</i>	MUCL 007956	No signal	No signal
<i>Disporotrichum dimorphosporum</i>	MUCL 019341	No signal	No signal
<i>Fusarium venenatum</i>	MUCL 055417	No signal	No signal
<i>Hansenula polymorpha</i>	MUCL 027761	No signal	No signal
<i>Humicola insolens</i>	MUCL 015010	No signal	No signal
<i>Kluyveromyces lactis</i>	IHEM 02051	No signal	No signal
<i>Leptographium procerum</i>	MUCL 008094	No signal	No signal
<i>Mucor javanicus</i>	IHEM 05212	No signal	No signal
<i>Penicillium camemberti</i>	IHEM 06648	No signal	No signal
<i>Penicillium chrysogenum</i>	IHEM 03414	No signal	No signal
<i>Penicillium citrinum</i>	IHEM 26159	No signal	No signal
<i>Penicillium decumbens</i>	IHEM 05935	No signal	No signal
<i>Penicillium funiculosum</i>	MUCL 014091	No signal	No signal
<i>Penicillium multicolor</i>	CBS 501.73	No signal	No signal
<i>Penicillium roqueforti</i>	IHEM 20176	No signal	No signal
<i>Pichia pastoris</i>	MUCL 027793	No signal	No signal
<i>Rhizomucor miehei</i>	IHEM 26897	No signal	No signal
<i>Rhizopus niveus</i>	ATCC 200757	No signal	No signal
<i>Rhizopus oryzae</i>	IHEM 26078	No signal	No signal
<i>Saccharomyces cerevisiae</i>	IHEM 25104	No signal	No signal
<i>Sporobolomyces singularis</i>	MUCL 027849	No signal	No signal
<i>Talaromyces cellulosilyticus/pinophilus</i>	IHEM 16004	No signal	No signal
<i>Talaromyces emersonii</i>	DSMZ 2432	No signal	No signal
<i>Trametes hirsuta</i>	MUCL 030869	No signal	No signal
<i>Trichoderma citrinoviride</i>	IHEM 25858	No signal	No signal
<i>Trichoderma longibrachiatum</i>	IHEM 00935	No signal	No signal
<i>Trichoderma reesei</i>	IHEM 05651	No signal	No signal
<i>Trichoderma viride</i>	IHEM 04146	No signal	No signal
Bacterial strains			
<i>Arthrobacter ramosus</i>	LMG 17309	No signal	No signal
<i>Bacillus amyloliquefaciens</i>	LMG 98140	No signal	No signal
<i>Bacillus brevis</i>	LMG 7123	No signal	No signal
<i>Bacillus cereus</i>	ATCC 14579	No signal	No signal
<i>Bacillus circulans</i>	LMG 6926T	No signal	No signal
<i>Bacillus coagulans</i>	LMG 6326	No signal	No signal
<i>Bacillus firmus</i>	LMG 7125	No signal	No signal
<i>Bacillus flexus</i>	LMG 11155	No signal	No signal
<i>Bacillus lentus</i>	TTAC 101	No signal	No signal
<i>Bacillus licheniformis</i>	LMG 6933T	No signal	No signal
<i>Bacillus megaterium</i>	LMG 7127	No signal	No signal
<i>Bacillus pumilus</i>	DSMZ 1794	No signal	No signal
<i>Bacillus smithii</i>	LMG 6327	No signal	No signal
<i>Bacillus subtilis</i>	W04-510	No signal	No signal
<i>Bacillus subtilis</i>	E07-505	No signal	No signal
<i>Bacillus subtilis</i>	LMG 7135 T	No signal	No signal
<i>Bacillus subtilis</i>	S10005	No signal	No signal
<i>Bacillus subtilis</i>	SUB033	No signal	No signal

<i>Bacillus subtilis</i>	BNB54	No signal	No signal
<i>Bacillus velezensis</i>	LMG 12384	No signal	No signal
<i>Bacillus velezensis</i>	LMG 17599	No signal	No signal
<i>Bacillus velezensis</i>	LMG 22478	No signal	No signal
<i>Bacillus velezensis</i>	LMG 23203	No signal	No signal
<i>Bacillus velezensis</i>	LMG 26770	No signal	No signal
<i>Bacillus velezensis</i>	LMG 27586	No signal	No signal
<i>Cellulosimicrobium cellulans</i>	LMG 16121	No signal	No signal
<i>Corynebacterium glutamicum</i>	LMG 3652	No signal	No signal
<i>Enterococcus faecium</i>	LMG 9430	No signal	No signal
<i>Escherichia coli</i>	LMG2092T	No signal	No signal
<i>Geobacillus caldoproteolyticus</i>	DSMZ 15730	No signal	No signal
<i>Geobacillus pallidus</i>	LMG 11159T	No signal	No signal
<i>Geobacillus stearothermophilus</i>	LMG 6939T	No signal	No signal
<i>Klebsiella pneumonia</i>	LMG 3113T	No signal	No signal
<i>Lactobacillus casei</i>	LMG 6904	No signal	No signal
<i>Lactobacillus fermentum</i>	LMG 6902	No signal	No signal
<i>Lactobacillus plantarum</i>	LMG 9208	No signal	No signal
<i>Lactobacillus rhamnosus</i>	LMG 18030	No signal	No signal
<i>Lactococcus lactis</i>	LMG 6890T	No signal	No signal
<i>Leuconostoc citreum</i>	LMG 9824	No signal	No signal
<i>Microbacterium imperiale</i>	LMG 20190	No signal	No signal
<i>Paenibacillus alginolyticus</i>	LMG 18723	No signal	No signal
<i>Paenibacillus macerans</i>	LMG 6324	No signal	No signal
<i>Protaminobacter pubrum</i>	CBS 574.77	No signal	No signal
<i>Pseudomonas amyloclavata</i>	ATCC-21262	No signal	No signal
<i>Pseudomonas fluorescens</i>	LMG1794T	No signal	No signal
<i>Pullulanibacillus naganensis</i>	LMG 12887	No signal	No signal
<i>Streptomyces aureofaciens</i>	LMG 5968	No signal	No signal
<i>Streptomyces mobaraensis</i>	DSMZ 40847	No signal	No signal
<i>Streptomyces murinus</i>	LMG 10475	No signal	No signal
<i>Streptomyces netropsis</i>	LMG 5977	No signal	No signal
<i>Streptomyces rubiginosus</i>	LMG20268	No signal	No signal
<i>Streptomyces violaceoruber</i>	LMG 7183	No signal	No signal
<i>Streptoverticillium mobaraense</i>	CBS 199.75	No signal	No signal
GMM samples			
<i>GM Bacillus subtilis</i> producing vitamin B2 (RASFF2014.1249)	/	No signal	No signal
GMM producing protease2*	/	No signal	No signal
PCR controls			
NTC	/	No signal	No signal






Validation of TaqMan screening qPCR method specific to the GMM protease1 in targeting the left transgene flanking region between the pUB110 shuttle vector and a *Bacillus velezensis* gene coding for a protease (GMM protease1 left border)

Annex 2: Sensitivity: full C_t dataset

Estimated full genome copy number	C _t values
100	33,41 32,94 33,01 32,7 33,38 33,15 33,02 33,2 32,99 33,57 33,53 33,77
20	35,23 35,51 35,1 35,5 35,79 35,61 36,19 36,06 36,53 35,62 36,24 35,3
10	35,51 37,15 36,66 37,03 36,88 36,61 37,34 36,45 36,67 36,8 37,37 37,61
5	36,79 39,21 37,45 38,22 38,47 39,45 37,76 38,35 37,66 37,62 37,6 39,74
1	41,05 No signal 39,36 No signal No signal 40,08 No signal 40,61 39,55 39,99 40,77 40,51
0.1	No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal
0	No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal



Annex 3: Transferability: full C_t dataset

Estimated full genome copy number	C_t values	
100	34,322	
	34,233	
	34,519	
	34,125	
	34,254	
	34,498	
	34,123	
	34,817	
	34,119	
	34,296	
	34,132	
	34,276	
	20	36,647
		36,035
36,468		
36,857		
36,275		
36,189		
36,718		
37,013		
36,692		
37,030		
36,073		
37,370		
10		37,590
		38,025
	37,342	
	37,356	
	37,416	
	37,763	
	36,882	
	37,596	
	37,567	
	37,701	
	37,772	
	37,315	
	5	38,506
		39,933
37,717		
40,245		
39,827		
39,063		
39,178		
38,911		
38,706		

		39,558
		39,501
		37,518
1		41,215
		40,263
		40,919
		38,596
		39,202
		41,056
		40,683
		No signal
		No signal
		40,045
		39,976
		39,927
0.1		40,945
		41,120
		No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal
0		No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal; No signal

