

Institute for Health and Consumer Protection

Comparative testing: 4th **EU-ASIA regional network meeting**

Diana Charels



www.jrc.ec.europa.eu

http://ihcpnet.jrc.it/

Serving society Stimulating innovation Supporting legislation



European Union Reference Laboratory for Genetically Modified Food and Feed (EU-RL GMFF)



Labelling of authorised Genetically Modified Organisms (GMOs) Prevent unauthorised GMOs on the European market

Two legal mandates:

- 1. Regulation (EC) No $1829/2003 \rightarrow validation of methods for detection and quantification of GM events$
- 2. Regulation (EC) No $882/2004 \rightarrow \text{official controls}$ applied to ensure the verification of compliance with feed and food law

23 May 2012 – 4th EU_ASIA Regional network meeting





Comparative testing (CT)

Article 32: Regulation (EC) No 882/2004

The EURLs for feed and food shall be responsible for:

Coordinating application by the NRLs of analytical methods, in particular by organising comparative testing and by ensuring appropriate follow-up of such comparative testing.

Proficiency testing (PT) = comparative testing





Aim of proficiency testing

For participating laboratories:

- Evaluation of laboratory performance
- Identification of problems in laboratories
- Education of participating laboratories

For customers, regulators and accreditation bodies:

- Ongoing confidence in laboratory performance
- ISO 17025 accredited laboratories need to show proof of participation in proficiency testing schemes to accreditation body





Launch of CT activities - 2010

- Two CT rounds per year
- On average 100 participants per round
- Each CT round: two test items containing different concentration

levels of one or more GM events

• Six weeks to conduct analyses





Preparation of test items








Preparation of test items (continued)





Final product

Bottling

23 May 2012 – 4th EU_ASIA Regional network meeting





Data analysis

- ISO 13528: Statistical methods for use in proficiency testing by interlaboratory comparisons
- Statistician
- Assigned value in-house versus calculation of robust mean
- z-scores





Data analysis (continued)

- Aim: obtain a simple and transparent result
- Submitted results are log-transformed
- Standard deviation for proficiency assessment σ
- Calculation of z-scores: reference value or consensus value from participants (i.e. robust mean $\hat{\mu}$)

$$z_{i} = \left(\log_{10} x_{i} - \hat{\mu}\right) / \hat{\sigma}$$

• Where $z_i = z$ -score of participant i $x_i = results$ of participant i





Participation of EU-ASIA in previous comparative testing rounds

- ILC-CRL-GMFF-CT-01/10: 5 participants (China, India, Malaysia, Republic of Korea, Singapore)
- ILC-CRL-GMFF-CT-02/10: 6 participants (India, Indonesia, Malaysia, Republic of Korea, Singapore, Vietnam)
- ILC-EURL-GMFF-CT-01/11: 5 participants (India, Malaysia, Republic of Korea, Singapore, Vietnam)





EURL-GMFF-CT-02/11

- Two maize test items containing different GM %ages
- List of 10 GM maize events: 3272, Bt11, Bt176, DAS59122, GA21, MIR604, MON 810, MON 863, NK603, TC 1507
- Qualitative PCR
- Detection of certain GM event \rightarrow quantify GM content
- Shipment: 24 October 2011
- Deadline submission results: 9 December 2011
- Test items: produced in-house by Marko Maras






Participants

- 158 laboratories invited
- 105 registered laboratories including also double registrations
- 15 laboratories of EU-ASIA network invited
- 8 participants of the EU-ASIA network: China, India, Indonesia, Malaysia, Singapore, Thailand, The Philippines, Vietnam





Distribution of participants



NRL - 882/2004

- **NRL 1981/2006**
- NRL both
- ENGL only
- Third country
- Official control lab only





Table 1: Overview of GM events present in
maize powder levels 1 and 2.

GM event	Maize powder level 1	Maize powder level 2
GA21	+	+
TC1507	+	+
MIR 604	+	+
NK603	Adv	Adv
59122	Adv	Adv
MON 810	Adv	Adv
MON 863	-	-
3272	-	-
Bt11	-	-
Bt176	-	-

+ indicates that the GM event was added to the test items, -indicates that the GM event was not added to the test items, -Adv: adventitious presence





Participants' performance Quantification of GA21, 1507 and MIR604

• Maize event GA21:

1. NRLs:	9					
2. ENGL only:	3					
3. Labs from 3 rd countries:	2					
4. Total:	14					
Maize event 1507:						
1. NRLs:	5					
2. ENGL only:	2					
3. Labs from 3 rd countries:	3					
4. Total:	10					
Maize event MIR604:						
1. NRLs:	8					
2. ENGL only:	3					
3 Labs from 3 rd countries:	1					

- 3. Labs from 3rd countries: 1
- 4. Total: 12





Participants' performance EU-ASIA network

- Six out of eight participants: qualitative analyses
- Two participants: maize events 1507 and MIR604
- Certificates of successful participation

Repetition of experimental work

- Shipment of test items to 18 laboratories of which 2 double sets
- Deadline for submission of results: 23 May 2012





Certificate of successful participation



Comparative Testing Certificate

The European Commission/ DG Joint Research Centre

certify that

Bureau of Plant Industry, Plant Quarantine Service

BPI Compound, Economic Garden 4030 Los Banos

has successfully participated in the comparative test : EURL-CT-02/11

on material(s): 1507 GA21 MIR604

launched in: June-2011

The laboratory has correctly detected the GM events GA21 and 1507

organised by the European Union Reference Laboratory for Genetically Modified Food and Feed

The EU-RL GMFF is accredited under ISO 17043 (Comparative Test Providers) by DAKKS (Accreditation D-EP-14322-01-00)





Reasons outlying z-scores

• Problems with calibration curve:

 $R^2 < 0.98$ and/or $-3.6 \le slope \le -3.1$

Calibrant		Ct-values	
3 % NK603		33.40	
3 % NK603		33.01	
5 % NK603		33.27	
5 % NK603		32.90	
Dilution factor	Theoretical ∆Ct	Actual ∆Ct	Theor. – actual ∆Ct
2	1	1.6	-0.6
10	3,322	2.322	1





Reasons outlying z-scores (continued)

• DNA quantity loaded on real-time plate = max. 200 ng


Overloading plate \rightarrow inhibition

- Overestimation of DNA concentration: (UV spectrophotometer) Not enough DNA loaded on plate → Ct values outside linear working range of calibration curve
- Swapping test items
- Use of one PCR replicate (i.e. one sample per well)





Calibration curve







Tips for CT participants

- Ensure labelling and traceability of test items \rightarrow avoid swapping of test items
- Calibration curve:
 - Dilution of standards
 - Quality criteria of calibration curve
 - Dilution of test items → within linear working range of calibration curve
- Good communication





Tips for CT participants (continued)

- Include controls:
 - Negative control \rightarrow contaminations
 - Positive control or quality control material \rightarrow check of quantification
- Reporting \rightarrow check:
 - Reported values
 - Measurement unit and conversion factors if applicable





Trick

- To avoid Ct values outside the linear working range of the calibration curve:
 - Concentration = 40 ng/ μ L for standards calibration curve
 - Concentration = 35 ng/µL for unknown samples and quality control materials





Final report

- Drafted after repetition of experimental work
- Results of repeated experiments are reported in a separate section → to see the effect on the z-score
- Draft report submitted to Advisory Board for comparative testing
- Meeting with Advisory Board to discuss report
- Finalisation of report
- Submission to PUBSY
- Each participant receives personal laboratory code along with report
- DG SANCO:
 - Only laboratory codes of NRLs are disclosed
 - Separate letter commenting on NRLs' performance





EU-RL GMFF Roberta Brustio **Diana Charels** Stéphane Cordeil Sara Jacchia Karolina Kolodziej Marko Maras Marco Mazzara Elena Nardini Dario Pardi Gregor Pinski Steven Price Maria-Grazia Sacco Cristian Savini Eleonora Scigliano Pierluigi Tenuta Rosemary Thompson Advisory Board Hermann Broll Bernard China Philippe Corbisier Hez Hird Lotte Hougs Martin Sandberg Manuela Schulze Isabel Taverniers IRMM FSQ Unit Fernando Cordeiro Raposo Beatriz De la Calle Inge Verbist





Commission

Thank you for your attention

